



APR 10, 2024

🌐 Assessment of implant accuracy using high-resolution postmortem MRI

Lucy Liang¹, Elvira Pirondini¹, Jonathan C Ho¹

¹University of Pittsburgh

ASAP Collaborative Research Network



Lucy Liang
University of Pittsburgh

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.14egn6x86l5d/v1

Protocol Citation: Lucy Liang, Elvira Pirondini, Jonathan C Ho 2024. Assessment of implant accuracy using high-resolution postmortem MRI. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn6x86l5d/v1>

MANUSCRIPT CITATION:

Liang, L., Zimmermann Rollin, I., Alikaya, A., Ho, J.C., Santini, T., Bostan, A.C., Schwerdt, H.N., Stauffer, W.R., Ibrahim, T.S., Pirondini, E., Schaeffer, D.J., 2024. An open-source MRI compatible frame for multimodal presurgical mapping in macaque and capuchin monkeys. *BioRxiv* <https://doi.org/10.1101/2024.02.17.580767>

ABSTRACT

Deep brain implant accuracy is important for successful experiments in non-human primates. In this protocol, we describe the steps to use postmortem imaging to assess the accuracy of an implant by visualizing implant location with a small thermal ablation and comparing its coordinates to image based pre-surgical planning.

This protocol is supplementary to the manuscript:

Liang, L., Zimmermann Rollin, I., Alikaya, A., Ho, J.C., Santini, T., Bostan, A.C., Schwerdt, H.N., Stauffer, W.R., Ibrahim, T.S., Pirondini, E., Schaeffer, D.J., 2024. An open-source MRI compatible frame for multimodal presurgical mapping in macaque and capuchin monkeys. *BioRxiv* <https://doi.org/10.1101/2024.02.17.580767>

PROTOCOL REFERENCES

Liang, L., Zimmermann Rollin, I., Alikaya, A., Ho, J.C., Santini, T., Bostan, A.C., Schwerdt, H.N., Stauffer, W.R., Ibrahim, T.S., Pirondini, E., Schaeffer, D.J., 2024. An open-source MRI compatible frame for multimodal presurgical mapping in macaque and capuchin monkeys. *BioRxiv* <https://doi.org/10.1101/2024.02.17.580767>

S.M. Smith, M. Jenkinson, M.W. Woolrich, C.F. Beckmann, T.E.J. Behrens, H. Johansen-Berg, P.R. Bannister, M. De Luca, I. Drobnjak, D.E. Flitney, R.K. Niazy, J. Saunders, J. Vickers, Y. Zhang, N. De Stefano, J.M. Brady, P.M. Matthews, Advances in functional and structural MR image analysis and implementation as FSL *Neuroimage*, 23 (2004), pp. S208-S219 <https://doi-org.pitt.idm.oclc.org/10.1016/j.neuroimage.2004.07.051>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 26, 2024

Last Modified: Apr 10, 2024

PROTOCOL integer ID: 97399

Keywords: ASAPCRN, postmortem, implant accuracy, non-human primate, MRI

Funders Acknowledgement:
Aligning Science Across Parkinson's
Grant ID: ASAP-020519

MATERIALS

Thermal Ablation

- Radiofrequency cannula (S-100 5 mm ActiveTip Straight cannula 22G, Abbott)
- Radiofrequency electrode (RF-SE-10 Reusable Stainless Steel Electrode, Abbott)
- Grounding pad (Abbott)
- Radiofrequency generator (NeuroTherm NT 1100, Abbott)

Imaging Equipment

- 9.4T/31 cm horizontal-bore Bruker AV3 HD animal scanner

Tissue Preparation

- Euthanasia drug (Fetal Plus or equivalent)
- Perfusion machine
- 1x phosphate buffered saline (PBS)
- 4% paraformaldehyde (PFA)
- Brain dissection tools

Image Processing Software

- FSL (Smith et al. 2004)

S.M. Smith, M. Jenkinson, M.W. Woolrich, C.F. Beckmann, T.E.J. Behrens, H. Johansen-Berg, P.R. Bannister, M. De Luca, I. Drobnjak, D.E. Flitney, R.K. Niazy, J. Saunders, J. Vickers, Y. Zhang, N. De Stefano, J.M. Brady, P.M. Matthews, Advances in functional and structural MR image analysis and implementation as FSL *Neuroimage*, 23 (2004), pp. S208-S219
<https://doi-org.pitt.idm.oclc.org/10.1016/j.neuroimage.2004.07.051>

BEFORE START INSTRUCTIONS

Follow this protocol after all terminal experiments have been completed, but before euthanizing the animal.

Note: Perfusion with 4% PFA may cause some degree of tissue shrinkage. We observed an average of 3.3% shrinkage with our brain tissue.

Thermal ablation & tissue preparation

- 1 Before euthanizing the animal according to protocol, perform a thermal lesion through the deep brain electrode.

We use a radiofrequency cannula (S-100 5 mm ActiveTip Straight cannula 22G, Abbott) and a radiofrequency electrode (RF-SE-10 Reusable Stainless Steel Electrode, Abbott), with a radiofrequency generator (NeuroTherm NT 1100).

Parameters: ~80°C, 1min

- 2 Remove electrode (not MR compatible) from the brain and euthanize the animal.
- 3 Perfuse transcardially with 1x phosphate buffered saline (PBS), followed by 4% paraformaldehyde (PFA).
- 4 Dissect out the brain with care. Carefully remove the dura and any blood clots from the surface of the brain.
- 5 Place the brain in 4% PFA for an additional 24 hours, then, move to 1x PBS.

Imaging and implant error calculation

- 6 Collect T2 or T2* weighted images of the brain in a 9.4 T/31 cm horizontal-bore Bruker AV3 HD animal scanner with the following parameters:
T2: 125 μ m isotropic, TR/TE = 1500/60ms, FOV = 52×80×56 mm
T2*: 80 μ m isotropic, TR/TE = 100/16 ms, FOV = 55×70×45 mm
- 7 Visualize the images in FSL, and align the anterior commissure to posterior commissure (ACPC) line with the x-axis of the mid sagittal plane with the FLIRT function.
- 8 Identify the tip of the electrode (at the center of the thermal lesion) within the internal capsule (or your target of interest), and record the coordinates.

- 9 Calculate the X_i , Y_i , and Z_i distance (in pixels) from the posterior commissure (PC) point. Convert to metric distance using pixel resolution.

(coordinate subscript "i" indicates implant coordinates, "t" indicates target coordinates)
- 10 Calculate target implant error in each axis by subtracting implant coordinates and target coordinates (e.g. $X_e = X_i - X_t$)
- 11 Finally, calculate the overall target implant error (IE) using Euclidean distance $IE = \sqrt{X_e^2 + Y_e^2 + Z_e^2}$