



NOV 16, 2022

DNA Extraction of Placenta Tissue

This protocol is published without a DOI.

Cayla Mason¹

¹University of California, San Diego



Cayla Mason

University of California, San Diego

COMMENTS 0

WORKS FOR ME

1

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

This protocol describes extracting DNA from placenta tissue, starting with dissociation.

PROTOCOL CITATION

Cayla Mason 2022. DNA Extraction of Placenta Tissue . **protocols.io**
<https://protocols.io/view/dna-extraction-of-placenta-tissue-cjejujcn>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 16, 2022

LAST MODIFIED

Nov 16, 2022

PROTOCOL INTEGER ID

72875

MATERIALS TEXT

1.0 mm diameter zirconia/silica beads, BioSpec Products, catalog 11079110z

DNeasy Blood and Tissue Kit, Qiagen, catalog 69504

100% ethanol

Nuclease-free water

2 mL screw-cap tubes with O-ring gaskets

Nuclease-free low-retention 1.5 mL microcentrifuge tubes

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

2h

Tissue Dissociation and Lysis

- 1 Place placenta tissue samples in screw cap tubes with O-ring gaskets on ice to thaw.
- 2 Label a 1.5 mL microcentrifuge tube for each placenta sample. Aliquot 1 mL silica beads per tube.
- 3 Add 540 uL Buffer ATL and 20 uL Proteinase K to each tissue. Pour in prepared beads.
- 4 Tightly close tubes. Place tubes in Bead Beater, and tightly close screw-on parts. Homogenize for 2 minutes.
- 5 Incubate lysates with beads in 56C water bath or heat block for 2 minutes.

- 6 Transfer lysate to pre-labeled microcentrifuge tubes. Do not discard beads.
- 7 Add 600 uL Buffer AL to screw-cap tubes with beads. Vortex and quick spin. Transfer supernatant to microcentrifuge tube with lysate. Discard screw-cap tube with beads.
- 8 Vortex lysate until homogenized.
- 9 Incubate lysates 56C for 10 minutes. Quick spin.
- 10 Freeze at -80C or proceed with DNA extraction.

1h 30m

DNA Extraction

- 11 Thaw lysates at room temperature.
- 12 Prepare columns and buffers. Add 25 mL 100% ethanol to Buffer AW1 and 30 mL 100% ethanol to Buffer AW2.
- 13 Add one-half volume (approximately 500 uL) 100% ethanol to lysate. Mix thoroughly by vortexing.
- 14 Pipet maximum 700 uL of the mixture to column. Centrifuge at $\geq 6000g$ for 1 minute. Discard flow through. Repeat this step until all lysate has passed through column.

- 15 Add 500 uL Buffer AW1. Centrifuge for 1 minute at $\geq 6000g$. Discard flow-through and collection tube.
- 16 Place the spin column in a new 2 mL collection tube. Add 500 uL Buffer AW2. Centrifuge for 3 minutes at 20000g. Discard flow-through and collection tube.
 - 16.1 While spinning, label 1.5 mL microcentrifuge tubes for elution.
- 17 Transfer column to pre-labeled microcentrifuge tubes.
- 18 Elute DNA by adding 100 uL Buffer AE or nuclease-free water to the center of the membrane. Incubate for 1 minute at room temperature. Centrifuge for 1 minute at $\geq 6000g$.
- 19 Repeat step 18.
- 20 Store eluted DNA at -20C.