

APR 02, 2023

OPEN BACCESS

DOI:

dx.doi.org/10.17504/protocol s.io.5jyl8jky8g2w/v1

Protocol Citation: Michael J Hurley 2023. STC-1 cell culture. **protocols.io** https://dx.doi.org/10.17504/p rotocols.io.5jyl8jky8g2w/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Feb 27, 2023

Last Modified: Apr 02, 2023

PROTOCOL integer ID:

77704

Keywords: ASAPCRN

STC-1 cell culture

Michael J Hurley^{1,2}

¹Department of Clinical and Movement Neurosciences, UCL Queen Square Institute of Neurology, London, UK;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Maryland, USA.



Michael J Hurley

ABSTRACT

Cell culture of STC-1 mouse enteroendocrine cells.

STC-1 cell culture

- 1 Obtain Mouse (*Mus musculus*) neuroendocrine duodenal adenoma STC-1 cells (RRID:CVCL_J405) from the American Type Culture Collection (ATCC-CRL-3254™).
- Grow STC-1 cells in DMEM:F12 Ham with GlutaMax™ media supplemented with 10 % charcoal absorbed foetal bovine serum (LabTech), 1 mM sodium pyruvate, non-essential amino acids (100 mM) and penicillin (100 U/ml) streptomycin (100 mg/ml) on 12 well plates at 37°C under saturating humidity in a 5% CO₂/95% air mixture.
- 3 Subculture cells when 80-90 % confluent using Accutase™ and replate at 1:10 dilution for maintenance or at an appropriate density for experiments.