



Mar 13,  
2020

# Live/Dead qPCR of *B. pertussis* IS481

Iain MacArthur<sup>1</sup>, Stacy Ramkissoon<sup>1</sup>, Andrew Preston<sup>1</sup>

<sup>1</sup>[Milner Centre for Evolution, Department of Biology and Biochemistry, University of Bath, Bath, UK]

1 Works for me dx.doi.org/10.17504/protocols.io.bc5niy5e

Stacy Ramkissoon

## ABSTRACT




This protocol describes the use of qPCR for the quantitative determination of live *B. pertussis* bacteria in a sample. The method uses propidium monoazide, a compound that can cross the cell wall of dead bacteria and binds DNA following photo activation to inhibit amplification of DNA from dead bacteria. This method uses a TaqMan probe specific to IS481 to ensure maximum sensitivity due to the high copy number of the target. Although IS481 is not specific to *B. pertussis* for experimental purposes it provides acceptable specificity.

## ATTACHMENTS

[SOP\\_001.3\\_live\\_dead\\_qPCR.docx](#)

## MATERIALS

NAME	CATALOG #	VENDOR
Water		
Centrifuge		
Microcentrifuge		
Microcentrifuge tubes	C2170	Denville Scientific Inc.
StepOnePlus™ Real-Time PCR System	4376600	Thermo Fisher Scientific
MicroAmp™ Fast Optical 96-Well Reaction Plate, 0.1 mL	4346907	Thermo Fisher
TaqMan™ Gene Expression Master Mix	4369016	Thermo Fisher
Adhesive PCR Plate Seals	AB0558	Thermo Fisher
PMA-Lite™ LED Photolysis Device	E90002	Biotium
1.5 mL Crystal Clear Microcentrifuge Tube	E1415-1500	StarLab
QIAamp DNA Mini Kit	51304	Qiagen
RNase A	R4642	Merck Millipore Sigma
PMA Dye 20 mM in H <sub>2</sub> O	40019	Biotium
1000 µl Filter Tip (Sterile) Racked	S1126-7810	StarLab
200 µl Graduated Filter Tip (Sterile) Racked	S1120-8810	StarLab
20 µl Bevelled Filter Tip (Sterile) Racked	S1120-1810	StarLab

NAME 	CATALOG # 	VENDOR 
10 µl Graduated Filter Tip (Sterile) Racked	S1121-3810	StarLab
Primers (IS481F/R)		
Probe (IS481)		
Control DNA for standard curve (Purified DNA from B1917 or strain being tested)		

#### MATERIALS TEXT

	Sequence (5'-3')
IS481 Forward Primer	ATCAAGCACCGCTTTACCC
IS481 Reverse Primer	TTGGGAGTTCTGGTAGGTGTG
IS481 Probe	FAM-AATGGCAAGGCCGAACGCTTCA-BHQ1

Primers (IS481F/R) and Probe (IS481)

#### SAFETY WARNINGS

Good Laboratory Practice must be followed. Wear laboratory coats and gloves.

#### BEFORE STARTING

##### Principle

To use propidium monoazide (PMA) to inhibit the PCR amplification signal of DNA from dead bacteria. Purified DNA will be used as a template for the quantification of *IS481* using a TaqMan probe. A standard curve consisting of DNA from pure B1917 will be used to determine absolute quantity in unknowns.

##### PMA treatment

- 1 Pellet fresh samples by centrifuging at 2,000 x g for 10 minutes and resuspend in 1.2 ml of PBS. 10m
- 2 Transfer 200µl of resuspended sample into two clear microfuge tubes (continue with 1 sample to Step #3 and 1 sample skip to Step #4).
- 3 Add 0.5µl of PMA to one of the 200µl of samples from Step #2.
- 4 Incubate microfuge tubes in the dark for 10 min at room temperature. Cover samples with aluminium foil and incubate on a rocker. 10m
- 5 Expose samples to light using the PMA-Lite™ LED Photolysis Device for 5 min. 5m

##### DNA purification

- 6 Add 20µl of QIAGEN Protease and 4 µl of RNase A.
- 7 Add 200µl of Buffer AL and vortex for 15 s. 15s
- 8 Incubate at 56°C for 10 min. 10m
- 9 Briefly spin tube to recover condensation.

- 10 Add 200µl of ethanol (96-100%), vortex and spin briefly.
- 11 Add to spin column and centrifuge for 1 min at 6000xg. 1m
- 12 Transfer to a new collection tube.
- 13 Add 500µl AW1 buffer and spin for 1 min at 6000xg. 1m
- 14 Transfer to a new collection tube.
- 15 Add 500µl AW2 buffer and spin for 3 min at 17000xg 3m
- 16 Transfer to a 2ml microfuge tube.
- 17 Spin at max speed for 1 min. 1m
- 18 Transfer to a 1.5ml microfuge tube.
- 19 Add 200µl of Buffer AE and incubate at room temperature for 1 min. 1m
- 20 Elute by centrifuging for 1 min at 6000xg. 1m

#### qPCR

- 21 Dilute stock primers (50µM) 9 in 50µl, to give a reaction concentration of 900nM.
- 22 Dilute stock probe (50µM) 3 in 100µl, to give a reaction concentration of 150nM.
- 23

Prepare Mastermix as follows:

Taqman MM	10µl	1x
Forward primer	2µl	900nM
Reverse primer	2µl	900nM
Probe	2µl	150nM

Mastermix per reaction

- 24 Add 16µl of master mix per well.
- 25 Serial dilute 9 times, 5:50 positive control DNA (25ng/ul).
- 26 Use dilutions 3-9 for 1000, 100, 10, 1, 0.1, 0.01, 0.001pg.
- 27 Load in triplicate 4µl of DNA for each control dilution, unknown and water samples.

#### PCR and analysis

- 28 Cover plate with film and spin for 1 min at 1000xg 1m

29 Set up qPCR program as follows:

	50°C	2min
	95°C	10min
40 cycles:	95°C	15sec
	60°C	1min

qPCR Cycle

## Results

30 Run analysis on data using StepOnePlus™ Software v2.3.

31 Outliers of triplicate samples should be ignored (omit).


32 Confirm standard curve has an  $r^2 > 0.95$ .

33 Convert pg of DNA to copy number:  $B. pertussis/\text{sample} = (50 \times (\text{mass of template in pg}) \times 6.022 \times 10^{23}) / (\text{length of genome in bp} \times 1 \times 10^{12} \times 650)$ .

34 Untreated sample (no PMA), represents total *B. pertussis* in the sample.

35 PMA treated sample gives number of live *B. pertussis* in the sample.

36 The number of dead can be determined from the difference between untreated and treated.

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited