

Sep 25, 2024 Version 2

C Lentivirus production & concentration V.2

This protocol is a draft, published without a DOI.

Allan JW Lui¹

¹Cancer Research UK, Cambridge Institute

Basic techniques



Allan JW Lui

Cancer Research UK Cambridge Institute

OPEN ACCESS



Protocol Citation: Allan JW Lui 2024. Lentivirus production & concentration. protocols.io https://protocols.io/view/lentivirusproduction-amp-concentration-dmzm4746 Version created by Allan JW Lui

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's

working

Created: September 30, 2021

Last Modified: September 25, 2024

Protocol Integer ID: 108301

Keywords: Lentivirus production, lentivirus concentration, lentivirus

Abstract

Lentivirus production protocol based on official protocol for Lipofectamine 3000 and TransIT-Lenti



Attachments



Lipofectamine3000-Le...

1MB

Materials

Protocol materials

GlutaMAX™ Supplement Thermo Fisher Catalog #35050038

Sodium Pyruvate (100 mM) **Thermo Fisher Catalog #**11360070 Step 7

Fetal Bovine Serum, qualified, One Shot™ format, Brazil Thermo Fisher Catalog #A3160801 Step 7

Opti-MEM™ I Reduced Serum Medium Thermo Fisher Catalog #31985047

Lipofectamine 3000 Thermo Fisher Scientific Catalog #L3000015

Before start

Prepare culture media for HEK293T:

MEM, low glucose, GlutaMAX™ Supplement, pyruvate Thermo Fisher Catalog #21885025 with 10% Fetal Bovine Serum, qualified, One Shot™ format, Brazil Thermo Fisher Catalog #A3160801



HEK293T seeding density titration

1d

Detach, count & seed HEK293T cells into T75 flasks at:

1. 4 x 10E6 cells

1h

2. 6 x 10E6 cells

- 3.8 x 10E6 cells
- 4. 10 x 10E6 cells

Note

>28 x 10E6 cells required in total; recommend detaching cells from 1x T175 at 70-80% confluence.

2 Incubate plates at 37 °C incubator (5% CO2, humidified) for up to 24:00:00 Observe confluence of cells under a microscope at 18 - 24 hours after plating.





The optimal plating density for transfection yield 80-90% confluence 18 - 24 hours after plating.

Lentivirus production (Lipofectamine 3000)

3 Detach, count & seed HEK293T cells at a density that yields 80-90% confluence 18 - 24 hours after plating

1d

Incubate at 37 °C incubator (5% CO2, humidified) for up to 24:00:00

Note

Lentivirus packaging medium (step 7 below) can be used when plating cells & during transfection the next day, but DMEM+10% FBS works well too

4 When HEK293T cells are at 80-90% confluence, transfect cells with lentiviral transfer and packaging plasmids using lipofectamine 3000



| | | Growth | Lipofectamine mix | | DNA mix | | | |
|------------------|-------------------|---------------|-------------------|----------------------------|------------------|----------|------------|--|
| Culture plate | Scaling factor | media (ml) | Opti-MEM (ml) | Lipofectamine 3000 (ul) | Opti-MEM (ml) | DNA (ug) | P3000 (ul) | |
| 6-well | 1 | 2 | 0.25 | 7 | 0.25 | 3 | 6 | |
| 100mm | 6 | 10 | 1.5 | 42 | 1.5 | 18 | 36 | |
| T75 | 8 | 15 | 2.0 | 56 | 2.0 | 24 | 48 | |
| 150mm | 16 | 30 | 4.0 | 112 | 4.0 | 48 | 96 | |
| T175 | 18 | 35 | 4.5 | 126 | 4.5 | 54 | 108 | |

Table 1.

If transfecting multple vessels, make a master mix of 1.1x the total volume required

Source document: Lipofectamine3000-LentiVirus-AppNo...

Calculate the amount of each plasmid required 4.1

Example using lentiCas9-Blast, psPAX2 and pMD2.G:

| A | В | С | D | E | F | G | Н |
|---------------------|---------------|----------------------|--------|----------------------------------|-------|-------|-------|
| | | | Mass c | Mass of plasmids per vessel (ug) | | | |
| Plasmid | Size (b p) | Equimolar ratio (ug) | 6-well | 10cm | T75 | 150mm | T175 |
| lentiCas9-Blas t | 12859 | 2.209 | 1.313 | 7.88 | 10.50 | 21.00 | 23.63 |
| psPAX2 | 10709 | 1.839 | 1.093 | 6.56 | 8.75 | 17.49 | 19.68 |
| pMD2.G | 5822 | 1.000 | 0.594 | 3.57 | 4.75 | 9.51 | 10.70 |
| | | | | | | | |
| Total | | | 3 | 18 | 24 | 48 | 54 |



Optimisation of plasmid ratios may be required for each transfer plasmid; an equimolar ratio is a good starting point. Zhang Lab uses a 4:3:2 transfer: psPAX2: pMD2.G ratio that approximates an equimolar ratio for lentiCas9-blast, lentiGuide-Puro & lentciCRISPR v2 plasmids. With lentiCas9-blast, 4:2:1, 9:8:1, 7.5:1.5:4 ratios all result in similar titres (+/-20% of 4:3:2), around 1-2e6 TU/ml.

4.2 Make up DNA mix by diluting transfer and packaging plasmids into Opti-MEM, then adding P3000 reagent (2ul / 1ug DNA)

3s

Mix well by vortexing (5) 00:00:03

4.3 Make up lipofectamine mix by diluting Lipofectamine 3000 (7ul / 3ug DNA) into Opti-MEM (volume specified in Table 1)

3s

Mix well by vortexing (5) 00:00:03

Note

Lipofectamine 3000 reagent diluted in Opti-MEM medium should be used within 15 minutes of dilution. Longer times can result in a loss of transfection efficiency

5 Combine both lipofectamine and DNA mixes and incubate for 00:20:00

20m

Room temperature to allow DNA-lipid complexes to form

Remove 50% volume of media from HEK293T culture vessels intended for transfection

| A | В | С | D | E | | |
|---|------------|-----|------------|------|--|--|
| Volume of media to remove per vessel (ml) | | | | | | |
| 6-well | 100mm dish | T75 | 150mm dish | T175 | | |
| 1 | 5 | 7.5 | 15 | 17.5 | | |

6 Gently add DNA-lipid complexes to cells

6h

| A | В | С | D | E | | |
|--|---|-----|------------|------|--|--|
| Volume of DNA-lipid complex to add per vessel (ml) | | | | | | |
| 6-well 10cm dish | | T75 | 150mm dish | T175 | | |
| 0.5 | 3 | 4 | 8 | 9 | | |



Gently rock culture vessel back-and-forth and from side-to-side to evenly distribute Incubate at 37 °C incubator (5% CO2, humidified) for 606:00:00

7 Prepare Lentivirus packaging medium (Opti-MEM + 1x GlutaMAX + 1mM sodium pyruvate + 5% FBS)

Supplement

500 mL

Opti-MEM™ I Reduced Serum Medium Thermo Fisher Catalog #31985047 with:

I 5 mL

GlutaMAX™ Supplement Thermo Fisher Catalog #35050038

I 5 mL

Sodium Pyruvate (100 mM) Thermo Fisher Catalog #11360070

I 25 mL

Fetal Bovine Serum, qualified, One Shot™ format, Brazil Thermo Fisher Catalog #A3160801

8 5-6 hours after transfection, replace media with pre-warmed packaging medium

Incubate plates at 37 °C incubator (5% CO2, humidified)

Consider supplementation with [M] 1 millimolar (mM) sodium butyrate or [M] 2 millimolar (mM) caffeine to improve titer.

Supplements will be present in viral supernatant if not concentrated

24 hours post-transfection, collect supernatant from culture vessels and store at 4 °C

Replace media with pre-warmed packaging medium

Incubate plates at 37 °C incubator (5% CO2, humidified) for up to 28:00:00

48-52 hours post-transfection, collect supernatant from culture vessels and combine with supernatant harvested the day before.

Centrifuge at 300 x g, Room temperature, 00:05:00 to separate any detached cells, and collect the supernatant.

11 Filter the supernatant through a 0.45µm PES filter to remove any remaining cellular debris

5m



Smaller filters (0.22µm) improve purity but can lower the viral titer.

12 Prepare single-use aliquots in cryovials and snap freeze virus in a dry ice-ethanol bath, then store at -80°C.

Alternatively, proceed with virus concentration or infection of cells.

Note

It is recommended to avoid subjecting lentiviral preparations to multiple freeze-thaw cycles, since each cycle can result in a 10%-20% loss in functional titers.

Lentivirus production (TransIT-Lenti)

13 Detach, count & seed HEK293T cells at a density that yields 80-90% confluence 18 - 24 hours after plating

1d

- Incubate at \$\mathbb{\
- 14 Change media for HEK293T cells
- 15 When HEK293T cells are at 80-90% confluence, transfect cells with lentiviral transfer and packaging plasmids using TransIT-Lenti

| Culture vessel | Scaling factor | Growth media (ml) | Opti-MEM (ml) | DNA (ug) | TransIT- Lenti (ul) |
|-------------------|-------------------|----------------------|------------------|----------|------------------------|
| 6-well | 1.0 | 2 | 0.2 | 2 | 6 |
| 10cm | 6.0 | 10 | 1.2 | 12 | 36 |
| T75 | 7.5 | 15 | 1.5 | 15 | 45 |
| 150mm | 15.0 | 30 | 3.0 | 30 | 90 |
| T175 | 17.5 | 35 | 3.5 | 35 | 105 |



If transfecting multple vessels, make a master mix of 1.1x the total volume required

15.1 Calculate the amount of each plasmid required

Example using lentiCas9-Blast, psPAX2 and pMD2.G:

| A | В | С | D | E | F | G | Н |
|---------------------|---------------|----------------------|--------|----------------------------------|------|-------|-------|
| | | | Mass c | Mass of plasmids per vessel (ug) | | | |
| Plasmid | Size (b p) | Equimolar ratio (ug) | 6-well | 10cm | T75 | 150mm | T175 |
| lentiCas9-Blas t | 12859 | 2.209 | 0.875 | 5.25 | 6.56 | 13.13 | 15.31 |
| psPAX2 | 10709 | 1.839 | 0.729 | 4.37 | 5.47 | 10.93 | 12.75 |
| pMD2.G | 5822 | 1.000 | 0.396 | 2.38 | 2.97 | 5.94 | 6.93 |
| | | | | | | | |
| Total | | | 2 | 12 | 15 | 30 | 35 |

Note

Optimisation of plasmid ratios is required for each transfer plasmid but an equimolar ratio is a good starting point, Zhang Lab uses a 4:3:2 ratio that approximates an equimolar ratio for lentiCas9-blast, lentiGuide-Puro & lentciCRISPR v2 plasmids.

15.2 Warm TransIT-Lenti reagent to room temperature and vortex 00:00:03

3s

15.3 Combine plasmids in quantities determined in step 14.1 Add Opti-MEM (Volume specified in Table 2), mix by pipetting Add TransIT reagent (3ul / 1ug DNA), mix by pipetting

10m

Incubate at Room temperature 00:10:00 to allow transfection complexes to form



Precipitation may be observed when excess DNA is used during complex formation. This

negatively impact transfection efficiency.

15.4 Add transfection copmlexes drop-wise to different areas of vessel Gently agitate culture vessel back-and-forth and from side-to-side to evenly distribute

2d

Incubate at 37 °C incubator (5% CO2, humidified) for 348:00:00

16 Harvest virus 48 hours post-transfection by following steps 10 - 12.

Lentivirus concentration

1h

1h

17

X LentiFuge Cellecta Catalog #LFVC1 Add ∡ 1 µL

∆ 1 mL lentiviral supernatant

Incubate at 👢 4 °C 01:00:00

18 Centrifuge at 12000 rpm, 4°C, 01:00:00

1h

Note

Pellets generated from small supernatant volumes may be difficult to see. It is advised to mark the centrifuge tube or bottle with a marker at the site where you expect the virus pellet to be.

- 19 Resuspend the pelleted lentiviral preparation in 1/100 of original volume using sterile phosphate buffered saline (PBS).
- 20 Prepare single-use aliquots and snap freeze virus in dry ice, then store at -80°C.