



DEC 27, 2022

Protocol 5: Agrobacterium-mediated transformation of *Spizellomyces punctatus* (Sp)

In 1 collection

Sarah M Prostack¹, Edgar M Medina¹, Erik Kalinka¹, Lillian Fritz-Laylin¹

¹University of Massachusetts at Amherst



fritzl原因linlab.umass

ABSTRACT

Once you have completed protocols 1-4, you are ready for transformation. This protocol is lengthy, requires accurate timing, and takes a minimum of 4 days to see results. Zoospores will need to be harvested from the plates prepared in Protocol 2; the process of harvesting is described below. Once Agro is at the appropriate OD, transformation is carried out by co-culturing Agro and Sp on induction media ("IM") to induce the expression of Agro's virulence genes and transformation of Sp cells. Co-culturing must be done on IM plates with premade depressions in the agar (see Protocol 4: Creating depressions in induction media plates).

ATTACHMENTS

[Spizellomyces transformation steps.pdf](#)

GUIDELINES

This section outlines the transformation process from growing Agro to the end of co-culturing. For harvesting the co-culture and selecting transformants, see *Protocol 6: Selecting for Spizellomyces punctatus transformants* and *Protocol 7: Picking colonies of transformed Spizellomyces punctatus* (Sp). All steps of this protocol, except necessary centrifugation, must be carried out in a sterile environment, either in a laminar flow hood or in the sterile area around an open flame.

MATERIALS TEXT

Materials

Overnight liquid cultures of desired *Agrobacterium tumefaciens* strains (see Protocol 3: Growing liquid cultures of *Agrobacterium* prior to transformation day for transforming Agro)

Active plates of wild-type *Spizellomyces punctatus*, grown without antibiotics (see Protocol 2: Culturing *Spizellomyces punctatus* (Sp) prior to transformation day)

IM liquid (see recipe)

IM plates with depressions (see Protocol 4: Creating depressions in induction media plates)

DS solution (see recipe)

Culture tubes, sterile (such as

VWR® Culture Tubes Plastic with Dual-Position Caps VWR Avantor Catalog #60818-703

)

50 mL conical tubes, sterile (such as CellTreat

Centrifuge Tubes-Bag 50mL Centrifuge Tube - Bag Sterile Cell Treat Scientific Products Catalog #229421

)

1.5 mL centrifuge tubes, sterile (such as

Fisherbrand™ Low-Retention Microcentrifuge Tubes Fisher Scientific Catalog #02-681-331

)

OPEN ACCESS

DOI:

dx.doi.org/10.17504/protocols.io.n92ldpp5x15b/v1

Protocol Citation: Sarah M Prostack, Edgar M Medina, Erik Kalinka, Lillian Fritz-Laylin 2022. Protocol 5: Agrobacterium-mediated transformation of *Spizellomyces punctatus* (Sp). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ldpp5x15b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Nov 03, 2022

Last Modified: Dec 27, 2022


PROTOCOL integer ID: 72235

Keywords: Agrobacterium, *Spizellomyces punctatus*, transformation


Pipette aid

Serological pipettes

100-1,000 µL micropipette such as


 Eppendorf Research Plus Single Channel pipette 100-1000 uL blue operating button for use with 100 **pipette.com Catalog #3123000063 ES-1000**

20-200 µL micropipette such as


 Eppendorf Research Plus single channel pipette 2-20 uL yellow operating button for use with 20 **pipette.com Catalog #3123000039**

Filter tips for the micropipettes, sterile such as

 TIPONE® FILTER TIPS **USA Scientific Catalog #1122-1830**

M4 repeat pipettor, such as  Repeater® M4 Starter Kit **Eppendorf Catalog #4982000322**


Combitips for the repeat pipettor, sterile such as

 Combitips® advanced - Positive Displacement Pipette Tips **Eppendorf Catalog #0030089804**


or similar 40 µm mesh filters (

 40 µm cell strainer **Fisher Scientific Catalog #22-363-547**)

25 mm syringe filter holders (

 Cole-Parmer Advantec 43303010 Polypropylene Filter Holder for 25 mm Membranes **Fisher Scientific Catalog #NC9972954**

), preloaded with Grade 1 Whatman paper (

 Cytiva Whatman™ Qualitative Filter Paper: Grade 1 Circles **Fisher Scientific Catalog #09-927-223**

), sterile

Plastic storage containers, or equivalent (e.x. plastic shoe box)

Centrifuge capable of holding 50 mL conicals

Equipment

Eppendorf™ 5810R Centrifuge

NAME

Centrifuge

TYPE

Eppendorf

BRAND

02-262-8187

SKU

<https://www.fishersci.com/shop/products/eppendorf-5810r-centrifuge-rotor-packages-16/022628187>

LINK

Laminar flow hood and/or open flame, for maintaining sterility.

70% (v/v) ethanol for maintaining sterility (if using laminar flow hood)

Shaking incubator at  28 °C

BEFORE START INSTRUCTIONS




On transformation day, Agro must be diluted to an OD₆₆₀ of 0.15 and then grown to an OD₆₆₀ of 0.6. Under our laboratory conditions, this takes about 4-6 hours, but this time should be empirically tested for each lab before starting to ensure proper and repeatable results. Be sure to take this growth time into account when planning all steps on and prior to transformation day.



[Spizellomyces transformation steps.pdf](#)

Growing *Agro* to the appropriate OD₆₆₀

- 1 Dilute each overnight *Agro* culture to an OD₆₆₀ = 0.15 using IM.

Note

In our experience, this works out to be  500 µL of overnight culture that has been spun for  4500 rcf, 00:05:00 and then resuspended into  4 mL of IM. Dilution conditions should be empirically tested for each laboratory.

- 2 Grow the *Agro* at  28 °C and shaking at  225 rpm until OD₆₆₀ = 0.6.

Note

Under our laboratory conditions, this takes about 4-6 hours. Growth time to reach optimal OD₆₆₀ should be empirically tested for each laboratory.


Harvesting *Sp* zoospores for transformation

- 3 Flood all wild type *Sp* plates with  1 mL of DS or IM for  01:00:00 .

1h

Note

Do this about one hour before *Agro* has reached an OD₆₆₀ of 0.6.
The solution used should have little to no effect on transformation success.

- 4 Harvest *Sp* zoospores by running  1 mL of fresh DS or IM along the agar of one plate, holding the plate at a 45 degree angle.


Note

Do this  00:20:00 before Agro has reached an OD₆₆₀ of 0.6.

5 Take the liquid from the first plate and use it to harvest the zoospores on a second plate in the same manner.

6 Do this for all plates of *Sp* that were flooded in step 3.

Note

If the volume becomes too much to hold in the serological pipette, empty the contents into a 50 mL conical tube (or other appropriately sized tube), then use  1 mL of this zoospore suspension to continue harvesting the remaining plates.

7 Pool all zoospores into a 50 mL conical (or other appropriately sized tube).

8 Pass the zoospore suspension through a 40 µm mesh filter into a new 50 mL conical tube.


Note

This removes clumps of cells that could clog a syringe filter.

9 Use a 18 gauge needle and an appropriately sized syringe to take up all of the 40 µm-filtered zoospore suspension.

Note

Be very gentle, the needle may damage the cells.

Alternatively: Directly load a 25 mL syringe with the 40 µm-filtered zoospore suspension by using a p1000 set to  1000 µL and gently pushing the cells into the syringe through the tip.

- 10 Use a sterile 25 mm syringe filter preloaded with grade 1 Whatman paper to further filter the spores into a new 50 mL conical tube (or other appropriately sized tube).



Note

PRESS SLOWLY, or risk damaging spores. As a rule of thumb, do not filter more than 5 mL of zoospores for one syringe filter.

This suspension is now to be used in co-culturing and should have a concentration of zoospores between $1-10 \times 10^6$ cells/mL.


4d 1h


Co-culturing Agro and Sp in varying ratios

- 11 Per each control and plasmid to be transformed, prepare and label four microcentrifuge tubes according to Figure 3.
- 12 Pipette the appropriate amounts of IM and the final solutions from protocols section "Growing Agro to the appropriate OD₆₆₀" and "Harvesting Sp zoospores for transformation" into the appropriate microcentrifuge tube in the order they appear in Figure 3.

Note

This will prevent cross contamination between strains.

Pipetting time can be sped up by using an M4 repeat pipettor (Eppendorf #4982000322) to dispense  50 µL volumes at a time, ensuring the proper amount is dispensed and that the liquid does not splash due to the force of dispensing. This is best done by dispensing along the tube walls just above the water line.

- 13 Gently pipette to mix the contents of the microfuge tube.
- 14 Transfer all  200 µL of an *Agro-Sp* co-culture to one of the premade depressions on a room temperature IM plate.

Note

All 4 co-culture tubes for one plasmid can be placed onto one IM plate. One tube per quadrant. Start on the left side and then go to the right side to lower the risk of contamination (or the opposite for left handed people).

- 15 GENTLY slide the IM plates to a place where they will not be disturbed.



Note

Avoid picking up the plates, as this may result in the pools of co-culture to combine, disturbing the ratios established and potentially lowering transformation efficiency or success.

- 16 Leave the plates to dry at 🌡 Room temperature for 12-24 hours.

- 17 Seal plates with parafilm, invert, and grow in a closed chamber (such as a plastic storage container) at 🌡 Room temperature for 4 days.

Note

Evidence of growth appears as an opaque, roughly circular area on the agar. If not all of the quadrants have obvious growth, that is okay. As long as there is at least one quadrant that shows growth, *Protocol 6: Selecting for Spizellomyces punctatus transformants* can be done, though success at gaining selectable transformants may be lower.