



Dec 03, 2020

qPercerus

Forked from qPercerus

Robert P Loewe¹, Benedikt Kirchner¹¹GeneSurge GmbH**1** Works for me dx.doi.org/10.17504/protocols.io.bp42mqye

Benedikt Kirchner

ABSTRACT

Protocol for qPercerus SARS-CoV-2 Test (attached document is giving more details).

Final product will be a commercial kit that is consisting of four components (lysis buffer, primer/probe mix identification, primer/probe mix confirmation, qPCR master mix). The initial test for SARS-CoV-2 is the identification assay (FAM assay for S gene with HEX assay as internal control) for patient screening. The confirmation assay (FAM assay for Orf3a with HEX assay as internal control) is used to re-evaluate the results.

ATTACHMENTS

[qPercerus SARS-CoV-2.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bp42mqye

PROTOCOL CITATION

Robert P Loewe, Benedikt Kirchner 2020. qPercerus. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp42mqye>

FORK NOTE

FORK FROM

[Forked from qPercerus, Robert Loewe](#)

KEYWORDS

SARS-CoV-2 Detection, RT-qPCR

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 26, 2020

LAST MODIFIED

Dec 03, 2020

PROTOCOL INTEGER ID

44922

ATTACHMENTS



[qPercerus SARS-CoV-2.pdf](#)

qPercerus Lysis/Extraction Buffer
qPercerus Primer-Probes Identification Assay
qPercerus Primer-Probes Confirmation Assay
qPercerus qPCR Mastermix

ABSTRACT



Protocol for qPercerus SARS-CoV-2 Test (attached document is giving more details).




Final product will be a commercial kit that is consisting of four components (lysis buffer, primer/probe mix identification, primer/probe mix confirmation, qPCR master mix). The initial test for SARS-CoV-2 is the identification assay (FAM assay for S gene with HEX assay as internal control) for patient screening. The confirmation assay (FAM assay for Orf3a with HEX assay as internal control) is used to re-evaluate the results.

- 1 Add  **200 µl Lysis Buffer** to reaction tube
- 2 Add **Nasopharyngeal** or **Throat Swab** to **Lysis Buffer** and **swivel** for  **00:00:10** to lysate and inactivate^{10s} virus. Remove and **discard** the **Swab** afterwards.

If **Nasopharyngeal** and **Throat Swab** are available from the same patient, the sensitivity can be increased by using both in one qPercerus reaction.

- 3 Briefly vortex and incubate at  **Room temperature** for  **00:10:00** 10m

Sample is stable in **Lysis Buffer** for up to  **48:00:00** at  **Room temperature** allowing for storage and transport.

- 4 Mix RT-qPCR Master Mix for first **Identification Assay** with 10% surplus and distribute on suitable qPCR microtiter plate:
 **10 µl ID Primer/Probe Mix** per reaction
 **5 µl Master Mix** per reaction
- 5 Add virus lysate to each well on qPCR microtiter plate:
 **10 µl Virus Lysate**

6 Briefly spin down qPCR microtiter plate
🌀 **1000 x g, Room temperature , 00:00:30**

7 Run **One-Step RT-qPCR** with the following temperature protocol:
🔥 **50 °C** ⌚ **00:10:00 reverse transcription**
🔥 **95 °C** ⌚ **00:03:00 RT Inactivation/ Polymerase activation**
For **45 cycles**:
🔥 **95 °C** ⌚ **00:00:10**
🔥 **60 °C** ⌚ **00:00:30 acquiring on FAM and HEX**

FAM channel detects **SARS-CoV-2** related fluorescence
HEX channel detects **human Internal Control** related fluorescence

8 Evaluate **Cqs** of **HEX** channel for swab validity:
Cq of **HEX** <30 sample valid for analysis
Cq of **HEX** >30 inconclusive -> perform **Confirmation Assay** for validation

Evaluate **Cqs** of **FAM** channel for SARS-CoV-2 presence:
Cq of FAM <35 positive for SARS-CoV-2
40 > Cq of FAM >35 inconclusive -> perform **Confirmation Assay** for validation
Cq of FAM >40 or NA negative for SARS-CoV-2

If national testing regulations demand validation of positive SARS-CoV-2 presence on an independent target sequence, perform **Confirmation Assay** for Cq of FAM <35 .

9 Mix RT-qPCR Master Mix for first **Confirmation Assay** with 10% surplus and distribute on suitable qPCR microtiter plate:
📦 **10 µl CONF Primer/Probe Mix** per reaction
📦 **5 µl Master Mix** per reaction

10 Add virus lysate to each well on qPCR microtiter plate:
📦 **10 µl Virus Lysate**

11 Briefly spin down qPCR microtiter plate
🌀 **1000 x g, Room temperature , 00:00:30**

12 Run **One-Step RT-qPCR** with the following temperature protocol:
🔥 **50 °C** ⌚ **00:10:00 reverse transcription**
🔥 **95 °C** ⌚ **00:03:00 RT Inactivation/ Polymerase activation**
For **45 cycles**:
🔥 **95 °C** ⌚ **00:00:10**
🔥 **60 °C** ⌚ **00:00:30 acquiring on FAM and HEX**

FAM channel detects **SARS-CoV-2** related fluorescence
Cy5 channel detects **human Internal Control** related fluorescence

- 13 Evaluate **Cqs** of **HEX** channel for swab validity:
Cq of **HEX** <30 sample valid for analysis
Cq of **HEX** >30 sample invalid, repeat sampling of

Evaluate **Cqs** of **HEX** channel for SARS-CoV-2 presence:
Cq of **FAM** <40 positive for SARS-CoV-2
Cq of **FAM** >40 or NA negative for SARS-CoV-2