



Ultracentrifuge clean-up of phage lysate for TEM V.1

Kathy N Lam¹

¹University of California, San Francisco

Version 1

Sep 01, 2020

1

Works for me

This protocol is published without a DOI.

Kathy Lam

University of California, San Francisco

PROTOCOL CITATION

Kathy N Lam 2020. Ultracentrifuge clean-up of phage lysate for TEM. **protocols.io**
<https://protocols.io/view/ultracentrifuge-clean-up-of-phage-lysate-for-tem-bkmkku4w>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 01, 2020

LAST MODIFIED

Sep 01, 2020

PROTOCOL INTEGER ID

41356

- 1 Take 1.1 ml of existing plate lysate
- 2 Take 1 ml to clean up; save remainder to titer pre- versus post-clean up
- 3 Transfer 1 ml to ultracentrifuge tube
- 4 Centrifuge 25,000 rpm (~25,000g) 60 min 4C
- 5 Pour off supernatant away from pellet; pipette off residual
- 6 Wash with 0.1 M ammonium acetate
- 7 Centrifuge 25,000 rpm (~25,000g) 60 min 4C

- 8 Pour off supernatant away from pellet; pipette off residual
- 9 Wash with 0.1 M ammonium acetate
- 10 Centrifuge 25,000 rpm (~25,000g) 60 min 4C
- 11 Pour off supernatant away from pellet; pipette off residual
- 12 Resuspend in 0.1 M ammonium acetate
- 13 Titer pre- and post-clean up samples
- 14 Store at 4C for up to 2 weeks