



Apr 15, 2020

Detection of five bacterial pathogens of rice by multiplex PCR

Martine Bangratz¹, Charlotte Tollenaere¹

¹IRD, Montpellier, France

1 Works for me dx.doi.org/10.17504/protocols.io.bcpaivie



ABSTRACT

This multiplex PCR protocol is used for the detection of five bacterial pathogens of rice : *Pseudomonas fuscovaginae*, *Burkholderia glumae* and *gladioli*, *Pantoea* spp, *Xanthomonas oryzae*, *Sphingomonas* spp.

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

M. Bangratz, I. Wonni, K. Kini, M. Sondo, C. Brugidou, G. Béna, F. Gnacko, M. Barro, R. Koebnik, D. Silué, C. Tollenaere (2020) Design of a new multiplex PCR assay for rice pathogenic bacteria detection and its application to infer disease incidence and detect co-infection in rice fields in Burkina Faso. PLoS One.

MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Agarose		Sigma
5X Hot Firepol Multiplex Mix Ready to load		
(NH ₄) ₂ SO ₄ 160 mM		
TBE buffer		
ethidium bromide		

MATERIALS TEXT

- Applied Biosystems Veriti 96-Well Thermal Cycler

- Oligonucleotides :

Pseudomonas fuscovaginae (amplicon size: 710 pb)

Pfs207-F : CAGTTCGATGGTCTGGGAAT

Pfs207-R : GGGACTGGTAAAGCACGGTA

Burkholderia glumae and *B. gladioli* (amplicon size: 508 pb)

tox_B_F : GCATTTGAAACCGAGATGGT

tox_B_Rd : TCGCATGCAGATAACCRAAG

Sphingomonas spp. (amplicon size: 435 pb)

Sphingo_KK_F1 : CGGCTGCTAATACCGGATGAT

Sphingo_KK_R1 : AGGCAGTTCTGGAGTTGAGC

Xanthomonas oryzae (amplicon size: 331 pb)

Xo3756F : CATCGTTAGGACTGCCAGAAG

Xo3756R : GTGAGAACCACCGCCATCT

Pantoea spp (amplicon size: 263 pb)
PAN_KK263F : GCGAGCCAATCGACATTA
PAN_KK263R : CGAGTAACCTGAGTGTTTCAG

BEFORE STARTING

- Wear clean gloves
- Clean and disinfect the PCR cabinet
- Defreeze DNA samples and reagents
- Gently mix the DNA samples and the 5X Hot Firepol master mix
- Vortex the (NH₄)₂SO₄ and the oligonucleotides
- Spin down the DNA sample and all reagents and keep them on ice
- Mark the 1.5 ml "Mix" tube and the 0.2 ml tubes for the PCR . Keep them on ice


- 1 The amounts described in this protocol are for one sample. Calculate the quantity you need based on your number of sample.
Always prepare 10 % more mix.

2

Prepare the PCR reaction mixture following the specifications below

Component	Volume (µl)
Template	2
H2O	14,55
Pfs207-F (5 µM)	0,2
Pfs207-R (5 µM)	0,2
ToxB_F (100 µM)	0,2
ToxB_Rd (100 µM)	0,2
Sphingo_KK_F1 (10 µM)	0,1
Sphingo_KK_R1 (10 µM)	0,1
Xo3756F (10 µM)	0,3
Xo3756R (10 µM)	0,3
PAN_KK263F (100 µM)	0,3
PAN_KK263R (100 µM)	0,3
(NH ₄) ₂ SO ₄	1,25
Hot Firepol Multiplex Mix RTL	5

- 3 Perform the amplification in a thermocycler in the following conditions

Hot start	95°C	12 min	
Amplification cycles			
Denaturation	94°C	30 sec	
Annealing	58°C	30 sec	
Extension	72°C	45 sec	
Final extension	72°C	7 min	
Hold	4°C		

30 cycles

- 4 Prepare a 2% agarose gel. Load 10µl PCR product and perform the electrophoresis at 100 Volts for 90min
- 5 Observe the gel under UV light after a ethidium bromide bath



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited