

Aug 09, 2024

## Absorbance assay from an M2P flower plate

 Forked from [Absorbance assay from an m2p flower plate](#)

DOI

**[dx.doi.org/10.17504/protocols.io.x54v9pr51g3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v9pr51g3e/v1)**

Tijana Radivojevic<sup>1,2,3</sup>, Matthew Incha<sup>1,2,3</sup>, Apostolos Zournas<sup>1,2,3</sup>, vblayroger<sup>1,2</sup>, Stephen Tan<sup>1,2,3</sup>, Hector Garcia Martin<sup>1,2,3,4</sup>

<sup>1</sup>Lawrence Berkeley National Laboratory; <sup>2</sup>Joint BioEnergy Institute; <sup>3</sup>Agile BioFoundry;

<sup>4</sup>Basque Center for Applied Mathematics

JBEI



**Matt Incha**

Agile BioFoundry

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.x54v9pr51g3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v9pr51g3e/v1)**

**Protocol Citation:** Tijana Radivojevic, Matthew Incha, Apostolos Zournas, vblayroger, Stephen Tan, Hector Garcia Martin 2024. Absorbance assay from an M2P flower plate. **protocols.io** **<https://dx.doi.org/10.17504/protocols.io.x54v9pr51g3e/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 18, 2023

**Last Modified:** August 09, 2024

**Protocol Integer ID:** 87946



#### **Funders Acknowledgement:**

**Hector Garcia Martin**

Grant ID: Department of

Energy DE-AC02-05CH11231

**Hector Garcia Martin**

Grant ID: Basque Center for


Applied Mathematics

CEX2021-001142-S

## Abstract

This protocol explains the steps for preparing samples and taking OD600 and OD340 measurements from a BioLector plate.

## Safety warnings

 Wear proper PPE.

## Before start


Make sure you booked the equipment in the calendar:

- 978-4-BIOMEK-NX-S8\_ SPECTRAMAX (4148) EQ (1) (or any other "plate-reader" capable of measuring OD340 and OD600)
- 978-4-BIOMEK-NXP (4148) EQ (1)

Make sure you received proper training before operating on the Biomeks.

## Required equipment/labware

### 1 Destination plate:



Flat-black clear-bottom Tecan plate x2

### Water plate:



deep reservoir (available at Robotics lab)

### Supernatant plate:



96-deep-well plate x2

### 2 Liquid handler:



## Equipment

### new equipment

NAME

Beckman Coulter

BRAND

Biomek NXp

SKU

Pipette tips needed (available at Robotics lab):

- tips s200 (green box)
- tips p1000 (yellow box)

Spectrophotometer/plate-reader:

- Molecular Devices SpectraMax M2
- or any other spectrophotometer capable of measuring OD340 and OD600

### 3 Additional components:

- H2O
- m2p plate with cultures

Centrifuge options:

- Avanti-15R-centrifuge
- Allegra-25R
- Eppendorf 5810R
- (any swinging-bucket centrifuge capable of spinning a 96-deepwell plate)
- NOT the Eppendorf 5430

## Biomek setup

### 4

**Note**

If you wish to set up a method in advance of a run, you may access the Biomek remotely.

- Follow the instruction provided in [this file](#).
- Request a password for a Windows Active Directory (AD) account from Arthur Panganiban (ahpanganiban@lbl.gov)
- Find your Biomek, use password **Robotp@ss978**

We will use two Biomek methods for this assay.


The **first** method will prepare:

- a Tecan plate for OD600 measurements, with 10x dilution of the samples
- two 96-deepwell plates as an intermediate step for preparing a Tecan plate for OD340 measurements - one with 1 mL of samples (supernatant plate) and the other with 1 mL of water for balance in the centrifuge


The **second** method will prepare:

- a Tecan plate for OD340 measurements (aliquoting cell supernatants from the 96-deepwell plate)


- 4.1 Open the Biomek Software application on the Biomek. At the top toolbar, select Project > Open Project to open the "Automation Dev" project.



4.2 At the top toolbar of the Biomek Software, select File > Open to open the first method.



To open the second choose:





## Biomek run



5

### Note

Do not start a Biomek run remotely if it is not in the Simulation mode.

Hit the play button (green arrow) to run the method.





- 5.1
- The first method takes around  00:05:00 before the pause and  00:05:00 after the pause
  - When the method pauses, remove the two 96-deepwell plates and continue the Biomek run

10m

5.2 Centrifuge run:

5m

- Place the two 96-deepwell plates sealed in the centrifuge and spin at max speed for  00:05:00

- 5.3
- Carefully remove the supernatant plate from the centrifuge, remove the seal and place it back in the Biomek
  - Start the second Biomek method to prepare a Tecan plate for OD340 measurements (takes around  00:05:00 )


5m

## SpectraMax run

- 6
- Start the Spectramax software "SoftMax"
  - Click "Open a data file" in the top left of the SoftMax window



- Open the method for measuring OD600 named **OD600.sda** in the "My\_Documents/Keasling\_Lab/inch" subfolder on the Spectramax (NX-S8) computer



- Open tray and place plate
- Measure OD600 from the first prepared 96-well Tecan plate
- Save the data (see step 7)
- Click "Open a data file"
- Open a method for measuring OD340 named **OD340.sda** in "My\_Documents/Keasling\_Lab/inch" subfolder on the Spectramax (NX-S8) computer
- Measure OD340 from the second prepared 96-well Tecan plate
- Save the data (see step 7)

## Saving data

- 7
  - Copy the values for OD600 from the SoftMax window into an empty Excel file and rename the spreadsheet to "600"
  - Copy the values for OD340 to another spreadsheet of the same file and rename to "340"
  - Save the Excel file as "OD.xlsx"

AutoSave OFF

Home Insert Draw Page Layout Formulas Data Tell me Share Comments

Clipboard Font Alignment Number Conditional Formatting Format as Table Cell Styles Cells Editing Analyze Data Sensitivity

A1

	A	B	C	D	E	F	G	H	I	J	K	L	M	N
1		1	2	3	4	5	6	7	8	9	10	11	12	
2	A	0.3323	0.3433	0.3387	0.3707	0.3337	0.3306	0.3287	0.3313	0.1255	0.136	0.1326	0.1312	
3	B	0.3354	0.3338	0.3427	0.4422	0.3399	0.3362	0.3371	0.3399	0.1392	0.1386	0.1345	0.1318	
4	C	0.3365	0.3143	0.3308	0.1652	0.3318	0.3335	0.3396	0.3415	0.1393	0.1381	0.1363	0.132	
5	D	0.3369	0.3379	0.3472	0.1636	0.3332	0.3284	0.3311	0.3378	0.1387	0.1388	0.1365	0.1315	
6	E	0.3381	0.3664	0.3386	0.1432	0.3423	0.3346	0.3341	0.3348	0.1389	0.139	0.1367	0.1318	
7	F	0.3419	0.3791	0.1615	0.1745	0.3456	0.3395	0.3387	0.3365	0.1392	0.1388	0.136	0.132	
8	G	0.1324	0.1368	0.1385	0.1399	0.1391	0.1379	0.1382	0.1397	0.1387	0.1374	0.1349	0.1314	
9	H	0.1308	0.1336	0.1352	0.1347	0.1351	0.1339	0.1342	0.1348	0.1346	0.1352	0.1333	0.1307	
10														

600 340 +

Ready 100%