



FEB 26, 2024

🌐 SOP for LBP ELISA after DSS-induced injury

Malu G Tansey¹

¹College of Medicine |University of Florida



Senthilkumar Karuppagounder

ABSTRACT

SOP for LBP ELISA after DSS-induced injury

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.q26g7pk73gwz/v1

Protocol Citation: Malu G Tansey 2024. SOP for LBP ELISA after DSS-induced injury.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.q26g7pk73gwz/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: Jan 29, 2024

Last Modified: Feb 26, 2024

PROTOCOL integer ID: 94307

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across
Parkinson's (ASAP)
Grant ID: ASAP-020527

Protocol

- 1 Concentrations of plasma LPS-binding protein (LBP) were examined via Mouse LBP SimpleStep ELISA according to the manufacturer's protocol (Abcam, #ab269542).
- 2 Plasma samples were prepared at a 1:16,000 dilution.
- 3 Standards, antibody cocktail containing detector and capture antibodies, and 1x wash buffer were prepared according to the manufacturer's protocol.
- 4 In the plate, 50µL of standard or sample was added to appropriate wells in duplicate.
- 5 50µL of antibody cocktail was added to all wells.
- 6 The plate was then incubated at RT for 1 hour on a plate shaker at 400 rpm.

- 7 Wells were decanted and washed three times with 150µL 1x wash buffer.
- 8 100µL of TMB Development Solution was added to each well and allowed to incubate in the dark for 12 minutes with shaking at 400rpm
- 9 100µL of Stop Solution was added to each well and placed on the shaker for 1 minute at 400 rpm to mix.
- 10 Absorbance of 450nm light was recorded with a FLUOstar Omega microplate reader (BMG Labtech).
- 11 The standard curve was created by plotting the blank-subtracted average LBP standard absorbance against the standard LBP concentrations.
- 12 Sample LBP concentrations were calculated by interpolating the blank-subtracted average absorbance values against the standard curve, multiplied by the dilution factor.