

JAN 30, 2024

# OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.k qdg3xo71g25/v1

Protocol Citation: mariangela.m assarocenere 2024. Optical densitometry of tyrosine hydroxylase fibers. protocols.io https://dx.doi.org/10.17504/protoc ols.io.kqdg3xo71g25/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Jan 26, 2024

Last Modified: Jan 30, 2024

### Optical densitometry of tyrosine hydroxylase fibers



In 1 collection

## mariangela.massarocenere<sup>1,2,3</sup>

<sup>1</sup>Department of Experimental Neuroscience, Santa Lucia Foundation IRCCS, Rome, Italy;

<sup>2</sup>Department of Systems Medicine, University of Rome Tor Vergata, Rome, Italy;

<sup>3</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, United States



Nicola Biagio Mercuri

#### **ABSTRACT**

Protocol for quantifying TH+ fibers in rat striatum and substantia nigra pars reticulata (SNpr)

### PROTOCOL integer ID: 94192

- 1 Acquire images with a light microscope of immunostained serial coronal sections covering the caudo-rostral of the brain regions (3sections/animal for striatum and SNpr)
- 2 Open ImageJ (RRID:SCR\_003070,https://imagej.net/)
- 3 Calibrate grey mean values for standard optical density values
- Open 8-bit images in ImageJ, trace the region of interest (striatum, SNpr), and measure the optical density.

  Then, in the same image, trace the region of interest in a blank area (cortex for striatum, surrounding neuropil for SNpr) and measure the optical density.
- Calculate the optical density by subtracting the background optical density from the striatal optical density (if traced region of interest areas are the same) or by the formula: SNpr Integrated density-[(SNpr area) x (background optical density)]