



Version 3

Jul 21, 2021

An X-HTDC method for estimating particulate phosphorus from microalgae V.3

Yingyu YY Hu¹, Zoe V. Finkel¹¹Dalhousie University

1 Works for me



Share

This protocol is published without a DOI.

Marine Microbial Macroecology Lab
Tech. support email: ruby.hu@dal.ca



Yingyu Hu
Dalhousie University

ABSTRACT

Here we describe a protocol to estimate particulate phosphorus associated with microalgae using an extra-high temperature dry combustion method (X-HTDC). An 800°C-combustion can completely decompose intracellular macromolecules, improving estimates of particulate phosphorus relative to methods that use lower combustion temperatures. Based on tests on *Prorocentrum triestinum* in exponential growth in f/2 medium, estimates of intracellular particulate phosphorus are increased by 10% compared with the traditional 500 °C combustion method.

We recommend microalgae be collected on polycarbonate filters and then combusted at 800 °C before digestion by 0.2 M HCl for 30 minutes at 90°C. The resulting orthophosphate is detected by mixing the digested sample with a mixture of molybdate and ascorbic acid to produce molybdenum blue (Chen 1956).

The working range of this assay is 0.074 to 200 µM orthophosphate. Minimum sampling volume (mL)=130/Chl-a_{ug}/L.

In order to assess the intracellular phosphorus in microalgae, we recommend an oxalate reagent (Tovar-Sanchez 2003) to wash the microalgae collected on the filter to remove surface adsorbed phosphorus.

P.S. Chen, T.Y. Toribara and Huber Warner. Microdetermination of Phosphorus. Anal. Chem..

<https://doi.org/10.1021/ac60119a033>

AntonioTovar-Sanchez, Sergio A Sañudo-Wilhelmy, Manuel Garcia-Vargas, Richard S Weaver, Linda C Popels, David A Hutchins. A trace metal clean reagent to remove surface-bound iron from marine phytoplankton. Marine Chemistry.

[https://doi.org/10.1016/S0304-4203\(03\)00054-9](https://doi.org/10.1016/S0304-4203(03)00054-9)

PROTOCOL CITATION

Yingyu YY Hu, Zoe V. Finkel 2021. An X-HTDC method for estimating particulate phosphorus from microalgae. **protocols.io**

<https://protocols.io/view/an-x-htdc-method-for-estimating-particulate-phosph-bwskpecw>

Version created by Yingyu Hu

KEYWORDS

particulate phosphorus, intracellular phosphorus, phosphomolybdenum-ascorbic reduction, orthophosphate, oxalate reagent, adsorbed phosphorus, X-HTDC, High temperature dry combustion

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 21, 2021

LAST MODIFIED

Jul 21, 2021

PROTOCOL INTEGER ID

51756

SAFETY WARNINGS

Polycarbonate filter can release toxic gas and smoke during combustion. An exhaust system is required for muffle furnace while using the X-HTDC method.

BEFORE STARTING

Crucibles are cleaned with phosphate free detergent and then rinsed with reverse osmosis water. After the crucibles were dried, 0.2 M HCl was added into the crucibles. The crucibles were covered and heated in the oven at **90 °C** for **00:30:00**. After acid was removed from the crucibles, rinse crucibles with reverse osmosis water and dry completely. Crucibles are then combusted at **500 °C** for **06:00:00** at ramp-rate **150 °C /h**. Some crucibles may shatter, but the remainders can intact after being combusted with samples at **800 °C**. Crucibles are numbered with pencil at the bottom of crucible cups.

Sampling

1 Sampling microalgae for total particulate phosphorus (i.e. intracellular phosphorus and adsorbed phosphorus)

- 1.1 Filter microalgae in liquid media onto polycarbonate filters, using gentle vacuum pressure (5 inches Hg).

Filter forceps
blunt end, stainless steel
Millipore XX6200006P

- 1.2 Rinse samples with filtered seawater

- 1.3 Place sample filters in 2 mL Cryogenic Vials.

Cryogenic Vials with Closures
Polypropylene, 2 mL
Corning® 66021-974



1.4 Filter blank media (without cells) through polycarbonate filter as blank.

1.5 Flash freeze filters and store at -20°C .

2 Sampling microalgae for intracellular particulate phosphorus

2.1 Filter microalgae in liquid media onto polycarbonate filters, using gentle vacuum pressure (5 inches Hg).

Filter forceps
blunt end, stainless steel
Millipore XX6200006P

2.2 Add  5 mL oxalate reagent onto the filter, and let oxalate reagent sit in the filter funnel for 5m
 00:05:00




Preparation of oxalate reagent
by Yingyu Hu,
Dalhousie University

PREVIEW

RUN

2.2.1 Add  50 mL MilliQ water in a 250 mL beaker.

2.2.2 Weigh  40 g NaOH and slowly pour into the beaker.

2.2.3 Use squeeze bottle to rinse the weighing boat and transfer rinse water into the same beaker.

2.2.4 Use glass rod to gently stir and fully dissolve NaOH.

The solution is very hot and corrosive. It can cause skin burns and eye damage.

- 2.2.5 Carefully transfer NaOH solution into 100 mL volumetric flask by using glass rod.
- 2.2.6 Rinse beaker with small amount of MilliQ water three times, transfer rinse water into the flask.
- 2.2.7 Mix the solution by gently shaking the capped volumetric flask and top to 100 mL with MilliQ water.
- 2.2.8 Transfer the prepared reagent into a 250 mL PP bottle.
- 2.2.9 Label the bottle with SDS pictogram.



- 2.10 In a 1000 mL beaker with stir bar, add **600 mL** MilliQ water.
- 2.11 Add **18.6 g** EDTA, **14.7 g** sodium citrate, **0.74 g** KCl and **5 g** NaCl into the beaker, stir until all ingredients are dissolved. **pH5.7**
- 2.12 **[M]10 Molarity (M)** NaOH is added dropwise to bring pH in between 6 to 7 by using a transfer pipet
- 2.13 Add **12.6 g** oxalic acid to the solution, stir the mixture while heating.
- 2.14 After oxalic acid is completely dissolved, stop heating and let it cool to room temperature. A water bath filled with tap water can be used to speed up cooling. **pH3.3**
- 2.15 Add **[M]10 Molarity (M)** NaOH dropwise to bring pH to **pH8**

2.16 Top to 1 L in volumetric flask with MilliQ water.

2.17 Filter oxalate reagent by rapid flow to a 1 L PP bottle.

Sterile Disposable Filter Units with PES
Membrane
Thermo Scientific™ Nalgene™ Rapid- 5964
Flow™ 520

2.18 Label the bottle and keep it at 🧊 **Room temperature** .

2.3 Drain and then rinse the sample with filtered seawater once

2.4 Place sample filters in 2 mL Cryogenic Vials.

Cryogenic Vials with Closures
Polypropylene, 2 mL
Corning® 66021-974

2.5 Filter blank media (without cells) through polycarbonate filter as blank.

2.6 Flash freeze filters and store at 🧊 **-20 °C** .

X-HTDC-ing

- 3 Mark number at the bottom of each crucible with pencil, log the following information:
- (1) The number of crucible
 - (2) The code of sample in the crucible

Porcelain crucibles
40 mL
VWR 89037-996

Crucible cover
VWR 71000-146


4 Transfer sample to crucible with clean filter forceps and lay filter at the bottom.

5 [M]0.17 M MgSO_4 reagent:

Dissolve  1.023 g MgSO_4 in 50 mL MilliQ water

 Magnesium sulfate anhydrous Fisher

Scientific Catalog #M65500

6 Add  200 μL [M]0.17 M MgSO_4 directly onto each sample and blank filter.

Sing-use pipet tip to avoid cross-contamination.

7 Cover the crucibles and place in the oven at  90 °C until samples are completely dry.

Forced air oven
VWR 89511-410

Remove samples out of the oven as soon as they are dried. If muffle furnace is not available, remove crucibles to vacuum desiccator.

8 Combust dried samples at  800 °C for  09:00:00

9h

Muffle furnace
F30428C
Thermo 10-505-13

Map the location of crucibles in the oven, in case pencil mark disappears under 800°C.

Ramp rate should be controlled at < **200 °C** /hour or follow the instruction provided by manufacture, otherwise the crucibles might shatter.
SP.RAT: 150/PAMPU: hour
Or
SP.RAT: 2.5/PAMPU: minute

- 9 Allow samples to gradually cool down in the muffle furnace.
- 10 Pencil mark on crucibles should be still visible, however, it can be easily removed by water. Therefore, when removing samples out of the furnace, label the lid and crucible with sharpie immediately.

Digesting

- 11 **0.2 M** HCl reagent:

In a reagent bottle, dissolve one part of **12 N** HCl in 59 parts of MilliQ water

12 N Hydrochloric acid Contributed by users

Volume of HCl_0.2M_mL = (5_mL) X (#Sample + #Blank)

- 12 Preheat oven to **90 °C**
- 13 Add 5 mL 0.2 M HCl to each crucible.
- 14 Gently swirl the crucible.
- 15 Cover the crucibles and place crucibles in the oven for **00:30:00** 30m
- 16 Cool samples down to **Room temperature**
- 17 Gently swirl the crucible and then transfer 500 ul solution to 2 mL microtube. Duplicate each sample and blank.

Maxymum Recovery® Snaplock
Microcentrifuge Tube
2.0 mL, Polypropylene, Clear, Nonsterile,
Axygen® MCT-200-L-C

18 Standard working solutions and reagents can be prepared during sample digestion.

19 KH_2PO_4 primary standard stock solution ($\approx 1 \text{ mM}$)

[Potassium dihydrogen orthophosphate ACP](#)

Chemicals Catalog #P-4550

19.1 Transfer about 1 g KH_2PO_4 into a beaker, cover the beaker with foil

19.2 Place the beaker into an oven, dry KH_2PO_4 at 110°C for at least 02:00:00

2h

19.3 Move KH_2PO_4 into a vacuum desiccator, allow KH_2PO_4 to cool to room temperature

19.4 Dissolve around 0.136 g dried KH_2PO_4 in 1 L milliQ water.

- Use 1 L volumetric flask
- Take notes of the actual weight of KH_2PO_4 for final concentration of standard stock solution

19.5 Transfer standard stock solution into a 1 L bottle and store in the fridge.

This stock solution lasts quite a long time, unless there is evidence for growth of algae or other extraneous biotic material.

20 Standard working solution

KH_2PO_4	Primary (ul)	MilliQ (ul)
S1	0	1000
S2	5	995
S3	10	990
S4	20	980
S5	50	950
S6	100	900
S7	150	850
S8	200	800

21 Transfer 500 ul of each standard working solution to 2 mL microtube.

22

All reagents are freshly prepared before colorimetric measurement.

23 [M]6 N (3 M) sulfuric acid reagent:

Carefully add 1 part [M]18 M concentrated sulfuric acid into 5 part MilliQ water

[\[M\]18M sulfuric acid](#) Contributed by users

24 [M]2.5 % ammonium molybdate reagent:

Weigh [\[M\]0.25 g](#) ammonium molybdate in a Falcon tube and top to [\[M\]10 g](#) with MilliQ water.

Cap and shake until totally dissolved.

[\[M\]Ammonium molybdate](#) Sigma

Aldrich Catalog #09878-100G

25 [M]10 % ascorbic acid reagent:

Weigh [\[M\]1 g](#) ascorbic acid in a Falcon tube and top to [\[M\]10 g](#) with MilliQ water;

Cap and shake until all dissolved.

[\[M\]Ascorbic acid](#) Sigma

Aldrich Catalog #A5960-100G

Wrap the tube with foil if the reagent is not used right after prepared.

26 Calculate the volume of molybdate-ascorbic reagent:


Total volume of reagent_mL = (0.5 mL) X (#standard working solution + #samples + #blanks)

27 Mix the reagents into Falcon tube:

A	B
Reagent	Part(s) as in volume
MilliQ	2
6N sulphuric acid	1
2.5% ammonium molybdate	1
10% ascorbic acid	1

28 Preheat incubator/shaker to [\[M\]37 °C](#)

SHAKING INCUBATOR
71L
Corning® LSE™ 6753

- 29 Add  500 µl reagent to each standard, sample and blank, starting from blanks, including blank for standards and blank for samples.

Finntip Stepper Tips
5 mL
Thermo Scientific 9404200

Before dispensing the reagent, wipe or dab the liquid drop on the outside of the tip, avoid wiping the open tip.

- 30 Vortex each tube.

- 31 Incubate at  37 °C for  03:00:00 while shaking at 200 rpm

3h

- 32 Load microplate with 250 µl reactant from each tube, duplicate.

	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>7</u>	<u>8</u>	<u>9</u>	<u>10</u>	<u>11</u>	<u>12</u>
<u>A</u>	S1	S1	Samples and sample blanks: 40 with duplicate									
<u>B</u>	S2	S2										
<u>C</u>	S3	S3										
<u>D</u>	S4	S4										
<u>E</u>	S5	S5										
<u>F</u>	S6	S6										
<u>G</u>	S7	S7										
<u>H</u>	S8	S8										

Example of loading the microplate

96-Well Microplates, Polystyrene, Clear,
Greiner Bio-One 655101

- 33 Read plate in microplate reader

A	B
Shake duration	00:00:05
Shaking type	Continuous
Shaking force	High
Shaking speed [rpm]	600
Wavelength [nm]	820
Use transmittance	No
Pathlength correction	No
Measurement Time [ms]	100

Varioskan LUX Multimode Microplate
Reader
Thermo Fisher VL0L00D0

Calculating

- 34 Subtract the average absorbance at 820 nm of the blank standard replicates from the absorbance at 820 nm of all other standard working solutions.
- 35 Subtract the average absorbance at 820 nm of the blank sample (i.e. blank filter) replicates from the absorbance at 820 nm of all other individual samples.
- 36 Prepare a standard curve by plotting the average blank-corrected 820 nm absorbance for each standard working solution versus its concentration in μM .
- 37 Use the standard curve to determine the orthophosphate concentration of each unknown sample by using its blank-corrected 820 nm absorbance.
- 38 $(\text{P per sample})_{\mu\text{g}} = (\text{orthophosphate})_{\mu\text{M}} \times (\text{V}_{\text{HCl}})_{\text{mL}} \times (0.001) \times (30.97)$