



JUL 28, 2023

immunofluorescent staining with anti-GFP and anti-CD63 antibodies

rosanne.wouters¹, Peter Vangheluwe¹

¹KU Leuven



rosanne.wouters

ABSTRACT

This protocol was used for immunofluorescent staining in fixed HeLa cells with anti-GFP and anti-CD63 antibodies, followed by confocal imaging.

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.14egn2pypg5d/v1

Protocol Citation: rosanne.wouters, Peter Vangheluwe 2023. immunofluorescent staining with anti-GFP and anti-CD63 antibodies.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.14egn2pypg5d/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jun 01, 2023

Last Modified: Jul 28, 2023

PROTOCOL integer ID:
82751

Keywords: ASAPCRN

- 1 plate cells on glass coverslips and grown until they reach 60% confluency

- 2 fix cells with 4% paraformaldehyde for 20 min at room temperature

- 3 permeabilize cells with 0.1% Triton X-100 in PBS for 5 min


- 4 block for 1 h with blocking buffer (PBS with 0.5% Tween20, 0.1% BSA, 0.2% FBS)

- 5 incubate coverslips with primary antibodies for 2 h at room temperature
(anti-CD63, exbio, 11-343-C100, mouse; anti-GFP, abcam, ab13970, chicken)

- 6 wash coverslips 3 times with PBS-T

- 7 incubated coverslips 30 min with secondary antibodies
(goat-anti-mouse-AlexaFluor647, goat-anti-chicken-AlexaFluor488)

- 8 wash coverslips 3 times with PBS-T

- 9 incubate coverslips with DAPI
- 10 wash coverslips 3 times with PBS-T
- 11 mount coverslips using
 FluorSave™ Reagent Merck MilliporeSigma (Sigma-Aldrich) Catalog #345789
- 12 images were acquired using an LSM780 confocal microscope (Zeiss) with a 10x or 40x objective
- 13 colocalization analysis was performed with Fiji plugin JacoP