




Feb 11, 2022

🌐 Measuring urea concentrations in water samples

Jacob Waldbauer¹, Amy Amy Zimmerman²¹University of Chicago; ²Pacific Northwest National Lab

1

dx.doi.org/10.17504/protocols.io.b4wmqxc6 Jacob Waldbauer

Colorimetric assay for direct (vs. enzymatic) measurement of urea to a detection limit of 0.4 μ M concentration. The reaction of urea with diacetylmonoxime (DAM) to form a colored product is enhanced by addition of thiosemicarbazide (TSC). The original direct method was adapted for use with a single mixed reagent and incubation at room temperature (vs. the “high temperature direct method” that incubates at 85°C for 30 minutes). *NOTE: This protocol is written for measurement in 24-well plates.*

DOI

dx.doi.org/10.17504/protocols.io.b4wmqxc6

Jacob Waldbauer, Amy Amy Zimmerman 2022. Measuring urea concentrations in water samples. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b4wmqxc6>

 protocol ,

Feb 10, 2022

Feb 11, 2022

58029

Spectrophotometer (or plate reader),
200 and 1000 μL pipettes,
Tube racks,
Vortexer,
200 and 1000 μL filter tips,
5 mL polypropylene tubes,
24-well microplate with lid, clear,
fresh Reagent A,
fresh COLDER reagent.

1 Making Standards.

- 1.1 Prepare 200 μM stock solution. Dilute 1:500 from 0.1 M solution \Rightarrow 20 μL + 9.980 mL nanopure water.
- 1.2 Dilute the stock solution to the following concentrations in nanopure water: 0, 0.5, 1.25, 2.5, 5, 7.5, 10 μM .

2 Making Reagents and Solutions

- 2.1 1. *Diacetylmonoxime (DAM) solution*: Dissolve 3.4 g in 100 mL nanopure or LC-MS water (34 g L^{-1} or 0.3363 M stock). Store solution at 4°C in dark. Stable for at least 1 month.
 - a. Also known as 2,3-butanedione monoxime.
 - b. Dissolve using rotisserie (hybridization oven) set to room temp (prop open door).
- 2.2 *Thiosemicarbazide (TSC) solution*: Dissolve 0.19 g in 20 mL nanopure water (9.5 g L^{-1} or 0.104235 M stock). Store solution at 4°C in dark. Stable for at least 1 month. a. Dissolve using rotisserie (hybridization oven) set to room temp (prop open door).
- 2.3 *Ferric chloride (hexahydrate) solution*: Dissolve 0.15 g ferric chloride ($\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$) in 10 mL nanopure water (15 g L^{-1} or 0.0554877 M stock). Store solution at 4°C in dark.
- 2.4 *Reagent A solution*: Mix 25 parts DAM with 1 part TSC. *Make fresh prior to*

each analysis. (32.7 g L⁻¹ or 0.3234 M DAM and 0.365 g L⁻¹ or 0.004 M TSC)

- 2.5 Reagent B solution: Add 300 mL concentrated sulfuric acid (~98% or 18.4 M) to 535 mL nanopure water (final concentration H₂SO₄ = 35.2% or 6.6 M). Add 0.5 mL ferric chloride solution to diluted acid (8.977 mg L⁻¹ or 0.03321 mM stock). Store solution at 4°C in dark. Stable for at least 1 month.
- 2.6 Color developing reagent (COLDER): Mix 1 part of Reagent A with 3.2 parts of Reagent B. **Use within 15 minutes**. a. Turbidity blank solution: Substitute nanopure water for Reagent A above for determination of optical turbidity blank.

3 Assay set-up.

- 3.1 Label four polypropylene tubes for each sample and all standards (7) (includes triplicate reactions and single turbidity blanks). NOTE: Assay requires >8 mL of each sample and standard (includes replication).
- 3.2 Calculate total volume of COLDER reagent needed to run triplicate reactions for each standard and sample: (7 standards + # samples + 1 extra) x 3 x 0.6 mL = total vol (mL). a. *If quantifying urea from 2 samples: 7+2+1 = 10 x 3 x 0.6 = 18 mL reagent.*
- 3.3 Calculate total volume of turbidity blank solution needed to run duplicates for each standard and sample: (7 standards + # samples + 1 extra) x 2 x 0.6 mL = total vol (mL). a. *If quantifying urea from 2 samples: 7+2+1 = 10 x 2 x 0.6 = 12 mL blank solution.*
- 3.4 Aliquot 2 mL of each sample and all standards into corresponding reaction tubes.
- 3.5 Prepare the volumes of COLDER reagent and turbidity blank solution needed by mixing 1 part of Reagent A (or nanopure water) with 3.2 parts of Reagent B as described below. **Use reagent within 15 minutes**.
a. If need 18 mL, mix 4.5 mL of Reagent A with 14.4 mL of Reagent B (18.9 mL).
b. If need 12 mL, mix 3 mL of nanopure water with 6 mL of Reagent B (12.6 mL).

Add 0.6 mL of COLDER reagent (or turbidity blank) to each reaction tube.

- 3.6 a. DAM: 1.7963 g L^{-1} or 0.01777 M final concentration
b. TSC: 0.02008 g L^{-1} or 0.0002203 M final concentration
c. FeCl_3 : $0.001578 \text{ g L}^{-1}$ or 0.00000584 M ($5.84 \text{ }\mu\text{M}$)
d. H_2SO_4 : 6.2% or 1.16 M final concentration
- 3.7 Mix (vortex) and incubate in the dark at room temperature ($\sim 22^\circ\text{C}$) for 3 days.
- 3.8 After 72 hours, transfer/pour entire 2.6 mL volume of each reaction into corresponding wells of 24-well plates and measure absorbance on plate reader.

4 Reading plates.

- 4.1 Turn on Tecan Infinite 200 PRO plate reader 20-30 minutes prior to use.
- 4.2 Once warmed up, open the iControl software on MLCLab-PC.
- 4.3 Open file "Revilla_urea_24well".
- 4.4 Load the plate—check whether the "plate with cover" box is checked (since using clear plates for this, can be read with lid on).
- 4.5 Read absorbance at $520 \pm 9 \text{ nm}$ (25 flashes). Program automatically opens an Excel file that documents read parameters and data. 'Save as' before measuring 2nd plate.

5 Analyzing data.

- 5.1 Subtract the average absorbance of sample turbidity blanks from the absorbance of the samples treated with COLDER reagent (= corrected sample absorbance).
- 5.2 Subtract the average absorbance of the nanopure water tubes (i.e., 0 μ M urea) treated with turbidity blank solution from the absorbances of all the standards.
- 5.3 Plot corrected absorbance (y) vs. concentration (x) for all standards to establish a standard curve with linear regression.
- 5.4 Use the equation of the standard curve to calculate sample concentration from absorbance.