



Jun 27, 2024

🌐 Ultra-high field, multi-echo MRI acquisition for resting-state functional connectivity

DOI

dx.doi.org/10.17504/protocols.io.dm6gpbbj1lzp/v1

Elizabeth Rizer¹, Dr. Scott Grafton²

¹UC Santa Barbara; ²University of California, Santa Barbara



Elizabeth Rizer

UC Santa Barbara

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.dm6gpbbj1lzp/v1

External link: <https://socalkinesia.org/high-resolution-connectomics/>

Protocol Citation: Elizabeth Rizer, Dr. Scott Grafton 2024. Ultra-high field, multi-echo MRI acquisition for resting-state functional connectivity. protocols.io <https://dx.doi.org/10.17504/protocols.io.dm6gpbbj1lzp/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 03, 2022

Last Modified: June 27, 2024

Protocol Integer ID: 57737

Keywords: Functional Magnetic Resonance Imaging (f..., resting state functional connectivity (r..., Parkinson's Disease, motor circuits, basal ganglia, ASAPCRN, ultra-high field, multi-echo

Funders Acknowledgement:

Aligning Science Across

Parkinson's

Grant ID: ASAP-020-519



Abstract

This is an MRI acquisition sequence used to examine resting-state functional connectivity of open and closed loop motor circuits in healthy controls and those with Parkinson's Disease. Scans include a Localizer, anatomical T1 MPRAGE, CMRR spin-echo AP and PA echo-planar imaging (EPI) images for field map generation, and a modified CMRR multi-echo EPI sequence for resting-state imaging. All participants complete a 30-minute MRI protocol at the University of Southern California's Center for Image Acquisition, where they are scanned with a Siemens 7T Terra scanner and a research-only 8-channel transmit 32-channel receive Nova 8Tx/32Rx head coil (Nova Medical Inc.)

Attachments



highres_fc_acquisiti...

73KB

Materials

Equipment

- Siemens 7T MAGNETOM Terra
- Nova Medical 8-channel pTx system and 8Tx/32Rx head coil



Protocol references

Puckett, A.M., Bollmann, S., Poser, B. A., Palmer, J., Barth, M., & Cunnington, R. (2018). Using multi-echo simultaneous multi-slice (SMS) EPI to improve functional MRI of the subcortical nuclei of the basal ganglia at ultra-high field (7T). *NeuroImage*, 172, 886–895. <https://doi.org/10.1016/j.neuroimage.2017.12.005>

Shao, X., Yan, L., Ma, S. J., Wang, K., & Wang, D. J. (2021). High-Resolution Neurovascular Imaging at 7T: Arterial Spin Labeling Perfusion, 4-Dimensional MR Angiography, and Black Blood MR Imaging. *Magnetic Resonance Imaging Clinics*, 29(1), 53-65.

Wu, X., Auerbach, E. J., Vu, A. T., Moeller, S., Van de Moortele, P. F., Yacoub, E., & Uğurbil, K. (2019). Human Connectome Project-style resting-state functional MRI at 7 Tesla using radiofrequency parallel transmission. *Neuroimage*, 184, 396-408.