



JUL 25, 2023

T cell depletion

Connor Monahan¹

¹Columbia University



Connor Monahan
Columbia University

ABSTRACT

This protocol details T cell depletion by antibodies targeting CD4 and CD8.

ATTACHMENTS

[786-2014.pdf](#)

MATERIALS

Reagents

Ultra-LEAF™ Purified anti-mouse CD4 Antibody BioLegend Catalog #100457

Ultra-LEAF™ Purified Rat IgG2b κ Isotype Ctrl Antibody BioLegend Catalog #400644

Ultra-LEAF™ Purified anti-mouse CD8a Antibody BioLegend Catalog #100746

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.rm7vzx372gx1/v1

Protocol Citation: Connor Monahan 2023. T cell depletion. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzx372gx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working





Created: Jul 13, 2023

Last Modified: Jul 25, 2023

PROTOCOL integer ID:
84950

Keywords: antibody, T cell, CD4 T cell, CD8 T cell, T cell depletion, flow cytometry

Procedure

- 1 Inject mouse intraperitoneally with  250 µg of Ultra-LEAF Purified anti-mouse CD4 or Rat IgG2b isotype control antibodies (Biolegend, Cat #100457 and Cat #400644; San Diego, CA) 3 and 1 days prior to immunization.
- 2 After immunization, inject mice with  250 µg of antibody weekly.
- 3 Assess cell depletion by flow cytometric analysis of mouse blood serum.
- 4 For depletion of CD8+ T cells, inject the mouse intraperitoneally with  250 µg of Ultra-LEAF Purified CD8+ or Rat IgG2b isotype control antibodies (Biolegend Cat # 100746 and Cat #400644; San Diego, CA) 2 days prior to immunization.
- 5 After immunization, inject mice with  250 µg of each antibody 9- and 19-days post immunization.
- 6 Assess cell depletion by flow cytometric analysis of mouse blood serum.