



Apr 13, 2021

## Colorimetric determination of urea

noah.langenfeld<sup>1</sup>, Laurenpayne<sup>1</sup>, Bruce Bugbee<sup>1</sup><sup>1</sup>Crop Physiology Laboratory, Utah State University**1** Works for me [dx.doi.org/10.17504/protocols.io.bsz3nf8n](https://dx.doi.org/10.17504/protocols.io.bsz3nf8n)

USU Crop Physiology Laboratory

noah.langenfeld

### ABSTRACT

This protocol measures the absorbance of urea in solution in complexation with diacetyl monoxime at 520 nm and is linearly proportional to concentration up to 4.6 mM urea.

### DOI

[dx.doi.org/10.17504/protocols.io.bsz3nf8n](https://dx.doi.org/10.17504/protocols.io.bsz3nf8n)

### PROTOCOL CITATION

noah.langenfeld , Laurenpayne , Bruce Bugbee 2021 . Colorimetric determination of urea. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bsz3nf8n>

### KEYWORDS

urea, urea assay, colorimetric urea, urea test, thiosemicarbazide, diacetyl monoxime

### LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Mar 04, 2021

### LAST MODIFIED

Apr 13, 2021

### PROTOCOL INTEGER ID

47899

### MATERIALS TEXT

thiosemicarbazide, diacetyl monoxime, water, sulfuric acid, phosphoric acid, ferric chloride

### SAFETY WARNINGS

This protocol utilizes strong acids. Ensure proper precautions and safety equipment are used when dealing with these chemicals.

#### Mixed Acid Reagent Preparation

1

Dissolve **2.5 mg** ferric chloride in **45 mL** DI H<sub>2</sub>O in a 250 mL volumetric flask.

2 Add **80  $\mu$ l** phosphoric acid.

3



Caution: Preparation of **18 Molarity (M)**  $\text{H}_2\text{SO}_4$  is highly exothermic. Slowly add acid to **100 mL** water and stir frequently to avoid spattering.

Prepare **18 Molarity (M)**  $\text{H}_2\text{SO}_4$  by diluting **65.25 mL** concentrated  $\text{H}_2\text{SO}_4$  up to **250 mL** with DI  $\text{H}_2\text{O}$ .

4 Dilute ferric chloride, DI water, and phosphoric acid mixture to **250 mL** with **18 Molarity (M)**  $\text{H}_2\text{SO}_4$ .

5 Mix until dissolved.

#### Mixed Color Reagent Preparation

6 Add **20.9 mg** diacetyl monoxime and **52.4 mg** thiosemicarbazide to a 250 mL volumetric flask.

7 Dilute to volume ( **250 mL** ) with DI  $\text{H}_2\text{O}$ .

8 Mix until dissolved.




#### Urea Assay

9 Fill a 600 mL beaker with **200 mL** DI  $\text{H}_2\text{O}$ .

10 Place the 600 mL beaker on hot plate and bring to a boil.

11 To prepare sample, aliquot **1 mL** into a 20 mL glass test tube.

12 Add **2 mL** Mixed Acid Reagent.

- 13 Add  2 mL Mixed Color Reagent.
- 14 Mix by inverting test tube at least 3 times.
- 15 To prepare blank, repeat steps 11-14 using  1 mL of DI H<sub>2</sub>O instead of the sample.
- 16 Label test tubes and place in boiling water bath for  00:20:00
- 17 Remove tubes from water bath and let cool to room temperature.
- 18 Fill plastic cuvette to line with cooled solution.
- 19 Place sample cuvette in spectrophotometer and record absorbance at 520 nm. Blank against the solution with both reagents, but no urea.