



Sep 06, 2022

Standard operating procedure for the isolation of genetically engineered hPSCs clones in a high-throughput way

Hanqin Li¹, Yogendra Verma¹, Dirk Hockemeyer¹, Frank Soldner²¹University of California, Berkeley; ²Albert Einstein College of Medicine

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.b4mmqu46

Devin E Snyder

ABSTRACT

This collection describes a standard procedure for isolating single human pluripotent stem cell (hPSC) clones in a high-throughput way. This collection follows nucleofection of hPSCs as described in detail in the collection "Nucleofection (Amaxa) and electroporation (Biorad) of hPSCs;" dx.doi.org/10.17504/protocols.io.b4qnvve

Collection overview

Seeding nucleofected hPSCs in 96-well plates using limited dilution

Duplicating 96-well plate-cultured hPSCs clones

Subcloning of genotype-confirmed hPSCs clones

General notes

1. Throughout these protocols, the term hPSC is used to collectively refer to both hiPSCs and hESCs. All described procedures have been tested and work equally well for hiPSCs and hESCs.

DOI

dx.doi.org/10.17504/protocols.io.b4mmqu46

COLLECTION CITATION

Hanqin Li, Yogendra Verma, Dirk Hockemeyer, Frank Soldner 2022. Standard operating procedure for the isolation of genetically engineered hPSCs clones in a high-throughput way. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b4mmqu46>



FUNDERS ACKNOWLEDGEMENT

Aligning Science Across Parkinson's
Grant ID: ASAP-000486

KEYWORDS

ASAPCRN

LICENSE

————— This is an open access collection distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 03, 2022

LAST MODIFIED

Sep 06, 2022

COLLECTION INTEGER ID

57741

MATERIALS TEXT

A	B	C
Item	Vendor	Catalog #
DMEM/F12	Thermo Fisher	11320082
DPBS w/o Calcium and magnesium (DPBS)	Corning	MT21031CV
Fetal Bovine Serum (FBS)	Corning	35-011-CV
Knockout Serum Replacement	Thermo Fisher	10828-028
L-Glutamine	Sigma	G8540
Penicillin & Streptomycin (100X)	Thermo Fisher	15140163
MEM Non-Essential Amino Acids (100X)	Thermo Fisher	11140050
Heat Stable Recombinant Human FGF2 *	Thermo Fisher	PHG0360
Y-27632	Chemdea	CD0141
2-Mercaptoethanol	Sigma	M3148
0.25% Trypsin with EDTA (Trypsin)	Thermo Fisher	25200114
DMSO	Fisher Scientific	BP231-100
Proteinase K	Sigma	P6556

*While we prefer Heat Stable Recombinant Human FGF2, we also have used regular FGF2

ABSTRACT

This collection describes a standard procedure for isolating single human pluripotent stem cell (hPSC) clones in a high-throughput way. This collection follows nucleofection of hPSCs as described in detail in the collection "Nucleofection (Amaxa) and electroporation (Biorad) of hPSCs;" dx.doi.org/10.17504/protocols.io.b4qnqvve



Collection overview



Seeding nucleofected hPSCs in 96-well plates using limited dilution
 Duplicating 96-well plate-cultured hPSCs clones
 Subcloning of genotype-confirmed hPSCs clones



General notes

1. Throughout these protocols, the term hPSC is used to collectively refer to both hiPSCs and hESCs. All described procedures have been tested and work equally well for hiPSCs and hESCs.

FILES

- 

Seeding nucleofected hPSCs in 96-well plates using limited dilution
Version 1
by Devin E Snyder
- 

Duplicating 96-well plate-cultured hPSCs clones
Version 1
by Devin E Snyder
- 

Subcloning of genotype-confirmed hPSCs clones
Version 1
by Devin E Snyder