



JUL 21, 2023

## Baiting *Pythium myriotylum* from Infested Soil

Nimalka M  
Weerasuriya<sup>1</sup>

<sup>1</sup>Oklahoma State University



Nimalka M Weerasuriya  
Oklahoma State University, USDA Agricultural Research Servic...

### DISCLAIMER

### OPEN ACCESS

**Protocol Citation:** Nimalka M Weerasuriya 2023. Baiting *Pythium myriotylum* from Infested Soil . **protocols.io** <https://protocols.io/view/baiting-pythium-myriotylum-from-infested-soil-cvxkw7kw>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development  
We are still developing and optimizing this protocol

**Created:** Jun 19, 2023

**Last Modified:** Jul 21, 2023

**PROTOCOL integer ID:**  
83660

### DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

### ABSTRACT

Baiting *Pythium* from Ft. Cobb or other infested soil.

## MATERIALS

P<sub>5</sub>ARP plates

CMA/PDA+amp plates

Soil with *Pythium* infestation signs

Seeds of susceptible peanut cultivar; we used PI 378012 from the USDA Germplasm Collection

Fume hood and tools for sterile plating technique

Sieve and beakers for surface sterilization

RO water

Forceps and scalpel/blade

95% ethanol

Flame

For seedling baiting:

4" pots with catch tray

greenhouse or adequate outdoor conditions

Mature banker plants (Alyssum and Peppers) for IPM

## Preparation

- 1 Prep P<sub>5</sub>ARP Plates 1-2 days before plating.  
Prepare working culture plates (CMA, PDA + amp) up to 1 week before plating.

## Seed Pod or Seedling Baiting

8m

- 2 Select whether you want to bait from diseased hulls in soil or from seedlings grown in diseased soil.

Step 2 includes a Step case.

**Peanut Seedlings**

**Diseased Hulls**

step case

### Peanut Seedlings




To bait from soil infested with *Pythium* using susceptible germinated seedlings.

- 3 Germinate susceptible seedlings using of  Sample by wrapping in RO-dampened paper towels and incubating for 2-4 days at  28 °C in the dark.



Check daily to re-wet paper towels as necessary.

- 4 Plant germinated seeds in infested soil, in triplicate, in 4" diameter pots with a drip tray. Maintain

greenhouse at  70 % humidity with temperatures between  24 °C and  30 °C .

**4.1** Water to maintain a level of 1-2" in drip tray for the first 2 days, and then to maintain soil moisture without excessive water until seedlings reach 4-6" in height.

**5** At 4 DAP, trim to 1 seed/pot, or carefully repot in separate pots.

**6** Grow for 15-20 days and check for rot symptoms.

**Note**

Brown rot will start at 4-6 DAP

Wilt at 6-18 DAP

Those exposed to higher temperatures were more susceptible to rot.

## Harvest and Surface-Sterilization

**7** Harvest seedlings. Rinse plants thoroughly in RO water to remove soil particles.

**8** Select appropriate samples with black rotted segments.



Rotted segments of peanut seedling suitable for plating

**8.1** Cut up segments into 1/4 inch (1-2 cm) pieces.

**9** Try with or without [M] 95 % (v/v) ethanol soak for ⌚ 00:04:00 or rinse in sterile RO water. Use a sieve in a small beaker.

4m

## CITATION

Garren, K. H. (1996). Peanut (Groundnut) Microfloras and Pathogenesis in Peanut Pod Rot. *Journal of Phytopathology*. 55(4), 359–367.

LINK

<https://doi.org/10.1111/j.1439-0434.1966.tb02238.x>

**Do not use bleach** due to *Pythium* sensitivity.

## CITATION

Stanghellini, M. E., & Kronland, W. C. (1985). Detrimental effect of surface sterilization on isolation of *Pythium* spp. from feeder roots (Abstr.). *Phytopathology*.

LINK

[https://www.apsnet.org/publications/phytopathology/backissues/Documents/1987Articles/Phyto77n08\\_1192.pdf](https://www.apsnet.org/publications/phytopathology/backissues/Documents/1987Articles/Phyto77n08_1192.pdf)

**9.1** Blot dry on sterile paper towels.

**10** Plate on P<sub>5</sub>ARP or other media with bacterial inhibition (CMA or PDA + amp).

## CITATION

Jeffers, S. N., & Martin, S. B. . Comparison of Two Media Selective for *Phytophthora* and *Pythium* Species. *Plant Disease*. 70(11), 1038–1043.

LINK

<https://doi.org/10.1094/PD-70-1038>

## Transfer Culture

**11** Check plates after 24-48 h.



- 12 Hyphal tip transfer using sterile loop to clean plates of CMA (V8 or PDA).

## Oospore Check

- 13 Check for "gold coin" oospores at plate edges to indicate *Pythium myriotylum*.



- Use CMA for oospore production
- Use full strength PDA to increase hyphal growth for DNA extraction
- Use V8 to increase hyphal growth for DNA extraction (untested)