

Oct 05, 2024

Cell culture, transfection, immunocytochemistry, and imaging

DOI

dx.doi.org/10.17504/protocols.io.eq2lyp55mlx9/v1

Cell culture,
transfection,
immunocytochemistry,
and imaging

Michael G G Hanna^{1,2}, Pietro De Camilli^{1,2}

¹Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815



Michael G G Hanna

Yale University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.eq2lyp55mlx9/v1

Protocol Citation: Michael G G Hanna, Pietro De Camilli 2024. Cell culture, transfection, immunocytochemistry, and imaging. protocols.io <https://dx.doi.org/10.17504/protocols.io.eq2lyp55mlx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: July 19, 2021

Last Modified: October 05, 2024

Protocol Integer ID: 51695

Keywords: Cell culture, transfection, immunocytochemistry, imaging, ASAPCRN

Abstract

This protocol is to help with the maintenance, transfection, immunocytochemistry, and imaging of adherent mammalian cells associated with the publication DOI above.

Attachments



Cell culture, transf...

67KB



Materials

☒ DMEM, high glucose **Thermo Fisher Catalog #11965118**

☒ Penicillin-Streptomycin (10,000 U/mL) **Gibco - Thermo Fisher Catalog #15140122**

☒ L-Glutamine (200 mM) **Gibco - Thermo Fischer Catalog #25030081**

☒ DMEM/F-12, HEPES **Thermo Fisher Catalog #11330057**

☒ ProLong[®]; Gold Antifade Mountant with DAPI **Thermo Fisher Catalog #P36935**

☒ Live Cell Imaging Solution **Thermo Fisher Catalog #A14291DJ**

Required buffers

Fixation buffer

400mM Sucrose

5mM KCl

1mM NaH₂PO₄

2mM MgCl₂

5mM Glucose

1mM EGTA, 5mM PIPES

Stored at 4°C

Dilute PFA to 4% (Electron Microscopy Sciences, #15710) in fixation buffer just before use

☒ 16% Paraformaldehyde Aqueous Solution EM Grade **Fisher Scientific Catalog #50-980-487**

PBS

Diluted to 1x from 10x stock; Sigma # D1408-500mL

☒ Dulbecco's Phosphate Buffered Saline **Sigma Aldrich Catalog #D1408-500ML**

Stored at room temperature

Antibody dilution (AbDil) buffer

1xPBS

4% BSA

0.1% Triton X-100

0.1% NaN₃

Stored at 4°C

Solubilization buffer

1xPBS

0.5% Triton X-100

Stored at room temperature



PBS wash buffer

1xPBS

0.1% Triton X-100

Stored at room temperature

TBS

50mM Tris pH 7.5

150mM NaCl

Stored at room temperature

Protocol materials

- ⊗ L-Glutamine (200 mM) **Gibco - Thermo Fisher Catalog #25030081** Materials, Step 1
- ⊗ DMEM/F-12, HEPES **Thermo Fisher Catalog #11330057** Materials, Step 2
- ⊗ ProLong[®]; Gold Antifade Mountant with DAPI **Thermo Fisher Catalog #P36935** Materials, Step 15
- ⊗ Live Cell Imaging Solution **Thermo Fisher Catalog #A14291DJ** Materials, Step 18
- ⊗ 16% Paraformaldehyde Aqueous Solution EM Grade **Fisher Scientific Catalog #50-980-487** Materials
- ⊗ Dulbecco's Phosphate Buffered Saline **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D1408-500ML** Materials
- ⊗ DMEM, high glucose **Thermo Fisher Catalog #11965118** Materials, Step 1
- ⊗ Penicillin-Streptomycin (10,000 U/mL) **Gibco - Thermo Fisher Catalog #15140122** Materials, Step 1

Safety warnings

! For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).



A. General cell culture

- 1 Culture COS-7 or HeLa (wild-type, parental control, and all knock-in lines) cells at 37 °C and 5 % CO₂ in DMEM (Gibco, #11965118) containing 10 % FBS , 1 millimolar (mM) sodium pyruvate , 100 U/mL penicillin , 100 mg/mL streptomycin (Gibco/LifeTech, #15140-122) and 2 millimolar (mM) L-glutamine (Gibco/LifeTech, #25030-081).
- 2 Culture hTERT-RPE1 (wild-type, parental control, and all knockout lines) cells at 37 °C and 5 % CO₂ in DMEM/F12 (Gibco, #11330057) containing 10 % FBS , 1 millimolar (mM) sodium pyruvate , 100 U/mL penicillin , 100 mg/mL streptomycin and 2 millimolar (mM) L-glutamine .

DMEM, high glucose **Thermo Fisher Catalog #11965118**

Penicillin-Streptomycin (10,000 U/mL) **Gibco - Thermo Fisher Catalog #15140122**

L-Glutamine (200 mM) **Gibco - Thermo Fischer Catalog #25030081**

DMEM/F-12, HEPES **Thermo Fisher Catalog #11330057**

Note

Note: For general maintenance, when cells reach 80-90% confluency, detach from the dish with Trypsin (Gibco, #25300054) and dilute 1:10-20 in a new dish.

Trypsin-EDTA (0.05%), phenol red **Thermo Fisher Scientific Catalog #25300054**

B. Cell transfection

- 3 For live-cell imaging experiments, seed cells on glass-bottom dishes (MatTek; 35mm) at the following concentrations: COS-7, 8.5×10^4 ; HeLa, 7.5×10^4 ; RPE1, 1.0×10^5 .
- 4 Allow cells to adhere for 8-24 hours before being transiently transfected using FuGene HD (Promega) following manufacturer's instructions.

**Note**

Note: Replace Media just before addition of the transfection reagent that contains no antibiotics to lower toxicity.

- 5 Incubate transfection reagents with cells for 08:00:00 - 36:00:00 before imaging.

1d 12h

**C. Immunocytochemistry**

2h 30m

- 6 For fixed-cell imaging experiments, seed cells on to glass coverslips (FisherBrand # 12542A) in a well of a 6-well plate (Corning #3516).

i. Perform seeding at the following concentrations for all fixed-cell experiments where cells are not additionally transfected: HeLa, 1.5×10^5 ; RPE1, 2.0×10^5 .

ii. Perform seeding at the following concentrations for all fixed-cell experiments where cells are additionally transfected: HeLa, 7.5×10^4 ; RPE1, 1.25×10^5 . Perform transfection as described above before fixation.

- 7 Just before fixation, wash cells with 1 mL PBS quickly.



- 8 To each well, add 1 mL fixation buffer + 4% PFA equilibrated at 37 °C . Perform fixation at 37 °C for a total of 00:10:00 .

10m



- 9 Remove Fixation solution, and wash cells with TBS three 10-minute washes to ensure the stoppage and removal of fixation solution.



- 9.1 Remove Fixation solution.

- 9.2 Wash cells with TBS for 00:10:00 .

10m

- 9.3 Wash cells with TBS for 00:10:00 .

10m







- 9.4 Wash cells with TBS for 00:10:00 . 10m
- 10 After the last TBS wash, add 1 mL solubilization buffer to each well for 00:10:00 at Room temperature . 10m
- 11 After solubilization is complete, add 1 mL AbDil buffer for 01:00:00 at Room temperature . 1h
- 12 Dilute antibodies in AbDil buffer. 100 μ L AbDil buffer containing the diluted antibody is used for each coverslip. Pipette the diluted antibody solution as a droplet onto parafilm in a humidity-controlled chamber.
- 13 Invert coverslips onto the antibody droplet (cells facing down into solution) and place at 4 °C for 16:00:00 - 24:00:00 . 1d
- 14 Move coverslips into individual wells of a clean 6-well plate (cells now facing upwards) and wash three times with PBS wash buffer for 00:10:00 each. 10m
- 14.1 Move coverslips into individual wells of a clean 6-well plate (cells now facing upwards).
- 14.2 Wash PBS wash buffer for 00:10:00 . 10m
- 14.3 Wash PBS wash buffer for 00:10:00 . 10m
- 14.4 Wash PBS wash buffer for 00:10:00 . 10m
- 15 After the final wash, partially dry coverslips using an aspirator and invert (cells facing down) onto a drop of mounting medium (Invitrogen #P36935) on a glass slide (Electron Microscopy Sciences # 4951-001T). (Allow slides to cure as directed by the manufacturer of the mounting medium).

ProLong[®] Gold Antifade Mountant with DAPI **Thermo Fisher Catalog #P36935**



- 16 After curing, permanently immobilize coverslips onto the glass slide using nail polish (Sally Hansen Advanced Hair and Nail).

D. Imaging

- 17 Perform all imaging using an Andor Dragonfly system equipped with a plan apochromat objective (60x, 1.4 NA, oil) and a Zyla scientific CMOS camera. 
- 18 Perform all live-cell imaging at  37 °C and  5 % CO₂ in live cell imaging buffer (Invitrogen #A14291DJ). 



Live Cell Imaging Solution **Thermo Fisher Catalog #A14291DJ**