

Jun 12, 2020

Semi-dry western blot and chemiluminescent detection

 [Nature Communications](#)

DOI


dx.doi.org/10.17504/protocols.io.puxdnxn

Douglas Adamoski¹, Sandra Martha Gomes Dias¹

¹Brazilian Biosciences National Laboratory (LNBio), Brazilian Center for Research in Energy and Materials (CNPEM), Zip Code 13083-970, Campinas, Sao Paulo, Brazil.



Douglas Adamoski



OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.puxdnxn>

Protocol Citation: Douglas Adamoski, Sandra Martha Gomes Dias 2020. Semi-dry western blot and chemiluminescent detection. [protocols.io https://dx.doi.org/10.17504/protocols.io.puxdnxn](https://dx.doi.org/10.17504/protocols.io.puxdnxn)

Manuscript citation:

Quintero M, Adamoski D, Reis LM, Ascensão CFR, Oliveira KRS, Gonçalves KA, Dias MM, Carazzolle MF, Dias SMG. Guanylate-binding protein-1 is a potential new therapeutic target for triple-negative breast cancer. **BMC Cancer** (2017) 17:727. DOI: [10.1186/s12885-017-3726-2](https://doi.org/10.1186/s12885-017-3726-2)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: April 30, 2018

Last Modified: June 12, 2020

Protocol Integer ID: 11895

Keywords: western, western blot, protein detection, semi-dry, semi dry blot, chemiluminescent detection, detection



Protocol materials

⊗ SuperSignal™ West Pico PLUS Chemiluminescent Substrate **Thermo Fisher Scientific Catalog #34577**

⊗ Wiper Wypall X60 Quartfold **Kimberly-Clark Catalog #34865**

⊗ Immun-Blot® PVDF Membrane 0,2µm **Bio-Rad Laboratories Catalog #1620177**

⊗ Glycine **Merck MilliporeSigma (Sigma-Aldrich) Catalog #50046**

⊗ Trizma® base **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T1503-1KG**

Troubleshooting

Safety warnings

- ! Methanol: Repeated or prolonged contact with skin may cause dermatitis. The substance may have effects on the central nervous system , resulting in persistent or recurring headaches and impaired vision. Wear safety goggles.



Buffers preparation

1 **Transfer buffer: Alcohol-free Timmons-Dunbar**

Weight reagents as indicated:


- 30.285 g of Tris-Base
- 14.41 g of Glycine


Transfer the reagents to a beaker, add 150 mL of ultrapure water, and stir using a magnetic bar until complete dissolution. Using a volumetric flask, complete the volume to 250 mL with ultrapure water. Store at room temperature.

Please cite Timmons TM, Dunbar BS. Protein blotting and immunodetection. *Methods Enzymol* (1990), 182:679-688. DOI: 10.1016/0076-6879(90)82053-5

Timmons-Dunbar 4X Transfer Buffer			
Reagent	Concentration	Amount	
Tris-Base	1000 mM	30.285 g	
Glycine	768 mM	14.41 g	
Ultrapure water	q.s.	250 mL	

Timmons-Dunbar Alcohol-Free Buffer

 Glycine Merck MilliporeSigma (Sigma-Aldrich) Catalog #50046

 Trizma® base Merck MilliporeSigma (Sigma-Aldrich) Catalog #T1503-1KG

2 **General Washing Buffer: Tris-buffered saline (TBS) 10X**

Weight reagents as indicated:

- 24,2 g of Tris-Base
- 87,6 g of Sodium chloride (NaCl)

Transfer the reagents to a beaker, add 900 mL of ultrapure water, and stir using a magnetic bar until complete dissolution. Using a volumetric flask, complete the volume to 250 mL with ultrapure water. Store at room temperature.



Gradient PAGE Gel preparation

3

Sample preparation and SDS-PAGE

4 Prepare your samples and run SDS-PAGE as usual.

Semi-dry transfer preparation

5 Dilute Timmons-Dunbar Buffer 1X using ultrapure water. Usually, 50 mL for each mini gel (8.4cm x 7cm) to be transferred are enough.


Note


Tip: Perform all steps from this section during SDS-PAGE run

6 Soak 12 WypAll X60 sheets (cutted in gel size) in Timmons-Dunbar Buffer 1X.

 Wiper Wypall X60 Quartfold **Kimberly-Clark Catalog #34865**

7 Pre-wet PVDF membrane (cutted in gel size) for 30 seconds in 100% methanol

 Immun-Blot® PVDF Membrane 0,2µm **Bio-Rad Laboratories Catalog #1620177**

 00:00:30 Membrane activation

8 Discard the methanol and overlay membrane with Timmons-Dunbar Buffer 1X

Note

Membrane will usually float

Semi-dry transfer

9 Remove gel from electrophoresis chamber and disassembly the sandwich


Note

Do **not** wash the gel! Remaining SDS is important for protein transfer.

10 Carefully and removing bubbles with an roller, assembly the following:

1. Bottom electrode
2. 6 x WypAll X60 buffer-soaked sheets
3. Pre-wet PVDF membrane
4. 6 x WypAll X60 buffer-soaked sheets
5. Top electrode

Follow the picture as an guideline:



11 Perform the transfer at 0.325 mA/mm^2 for 7 minutes when using 0.75mm mini gels.

The following table (From Trans-Blot Turbo Blotting System (Cat #170-4155) manual) is quite usefull:

	Protocol Name	MW (kDa)	Time (min)	2 Mini Gels or 1 Midi Gel	1 Mini Format Gel
	0.75 mm thick gels	> 150	10	2.5 A; up to 25 V	1.3 A; up to 25 V
	Mixed	7			
	< 30	5			
	1.5 mm thick gels	Any	10		


- 12 Disassembly the sandwich and proceed with blocking and antibody incubation.

Note


Tip: Membranes can be stained with Ponceau S to check protein transfer and binding. Gels can be checked for protein leftover by Coomassie Blue staining.

Blocking and antibody incubation

- 13 Immediately transfer the membrane to 3% non-fat dry milk diluted in TBS-T (Tris Buffered Saline with 0.05% Tween 20) and incubate with gently rocking for 1 hour.

 01:00:00 Membrane blocking

- 14 Pour off blocking solution and overlay membrane with TBS-T, let it gently rocking for 5 minutes (wash 1/3)

 00:05:00 Wash (1/3)


- 15 Pour off washing solution and overlay membrane with TBS-T, let it gently rocking for 5 minutes (wash 2/3)

 00:05:00 Wash (1/3)

- 16 Pour off washing solution and overlay membrane with TBS-T, let it gently rocking for 5 minutes (wash 3/3)

 00:05:00 Wash (1/3)

- 17 Pour off washing solution and overlay with desired antibody. Incubate from 1 hour (room temperature) up to 16 hours (4°C) with gently rocking.

 4 °C for overnight (16h) incubations

Note

Tip: Follow manufacturer's instructions for primary antibody dilution. As a starting point, 1:1,000 dilution in TBS-T with 1% BSA (Bovine Serum Albumin) and 0.02% sodium azide (in order to prevent microorganism growth and antibody reuse).

- 18 Collect antibody and store at 4°C for further uses. Overlay membrane with TBS-T, let it gently rocking for 5 minutes (wash 1/3).

 00:05:00 Wash (2/3)

- 19 Pour off washing solution and overlay membrane with TBS-T, let it gently rocking for 5 minutes (wash 2/3).

 00:05:00 Wash (2/3)

- 20 Pour off washing solution and overlay membrane with TBS-T, let it gently rocking for 5 minutes (wash 3/3).

 00:05:00 Wash (2/3)

- 21 Pour off washing solution and overlay with horseradish peroxidase conjugated antibody. Incubate from 1-3 hours at room temperature with gently rocking.


Note

Tip: Follow manufacturer's instructions for primary antibody dilution. As a starting point, 1:10,000 dilution in TBS-T.

- 22 Pour off the secondary antibody. Overlay membrane with TBS-T, let it gently rocking for 5 minutes (wash 1/3).

 00:05:00 Wash (2/3)

- 23 Pour off washing solution and overlay membrane with TBS-T, let it gently rocking for 5 minutes (wash 3/3).

 00:05:00 Wash (2/3)

- 24 Pour off washing solution and overlay membrane with TBS-T, let it gently rocking for 5 minutes (wash 3/3).

 00:05:00 Wash (2/3)

Chemiluminescent x-ray film development

- 25 Mix 1 mL of 'Stable Peroxide Solution' with 1 mL of 'Luminol Enhancer Solution' in a clean container.



SuperSignal™ West Pico PLUS Chemiluminescent Substrate **Thermo Fisher Scientific Catalog #34577**

Note

Tip: Any chemiluminescent reagent can be used

- 26 Gently transfer the membrane to this container using tweezers and rock for a few seconds.

Note

Tip: Grab the membrane by the border in order to avoid tweezer marks over the bands

- 27 Transfer the wet membrane to the X-ray cassette. The best approach is 'sandwich' it with two transparent hard plastic foils (as picture bellow).



Note

Tip: X-ray cassettes with white enhancing screens give better results



- 28 Inside the dark room (remember to left your bright-screen cellphone outside), add one X-ray film over the assembled sandwich and expose for 1 minute.

Note

Tip 1: Increase or reduce exposure time as needed.

Tip 2: Cheaper "hospital grade" X-ray films can be used in place of expensive "research grade" films.

- 29 Remove the film from cassette and put in a container with diluted autoradiography developer solution and incubate for 30 seconds.


 00:00:30 Developer

Note

Tip: Replace solution as it becomes dark

- 30 Rinse with tap water

- 31 Put in a container with diluted autoradiography fixer solution and incubate for 1 minute

 00:01:00 Fixer

Note

Tip: Replace solution as it becomes yellowish

- 32 Rinse with tap water and let it dry

 [go to step #28](#) If exposure is underoptimal

- 33 Overlay the dried film over the original cassette sandwich to tick the prestained marker, stick to blank sheet of paper and scan it.