



Aug 12, 2021

Protocol for RNA isolation and RT-PCR confirmation for Neomycin gene expression in transformed *Bodo saltans*

Fatma Gomaa¹, Zhu-Hong Li², Roberto Docampo³, Virginia Edgcomb⁴¹Harvard University; ²The University of Georgia; ³University of Georgia; ⁴Woods Hole Oceanographic Institution

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bxbxpipn

Protist Research to Optimize Tools in Genetics (PROT-G)

 Fatma Gomaa
Harvard University

ABSTRACT

Developing transfection protocol for *Bodo saltans*, using SaCas9/sgRNA ribonucleoprotein (RNP) complex in conjunction with DNA repair template to disrupt the Paraflagellar rod 2 gene (*BsPFR2*) and increase the efficiency of targeted homologous recombination when a repair template DNA is provided. The exogenous repair template is double stranded DNA and it consists of *eGFP* fused with the drug selection gene *nptII/neo* and flanked by 500 bp of the untranslated regions (UTRs) upstream and downstream of the targeted *BsPFR2* as homologous repair arms.

DOI

dx.doi.org/10.17504/protocols.io.bxbxpipn

PROTOCOL CITATION

Fatma Gomaa, Zhu-Hong Li, Roberto Docampo, Virginia Edgcomb 2021. Protocol for RNA isolation and RT-PCR confirmation for Neomycin gene expression in transformed *Bodo saltans*. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bxbxpipn>

LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 12, 2021

LAST MODIFIED

Aug 12, 2021

PROTOCOL INTEGER ID

52311

- 1 -Total RNA was isolated from 4 months old transformed clones and wild type *Bodo saltans* cultures using the Quick-RNA Miniprep Kit (Zymo Research) according to the manufacturer's instructions.
- 2 -RNA underwent two rounds of DNase treatments. First, we performed an on-column DNase digestion, following the protocol of the Quick-RNA Miniprep Kit (Zymo Research). Second, DNase digestion was performed on the purified RNA using gDNA wipeout buffer and the QuantiTect Reverse Transcription Kit (Qiagen, Germany), in the following 14 µl

reaction: 2 µl of gDNA wipeout buffer 7x, 6 µl of RNA, and 6 µl of RNase-free water.

- 3 -cDNA was synthesized using the QuantiTect Reverse Transcription Kit (Qiagen, Germany) according to the manufacturer's instructions, the RT Primer Mix and the 5X Quantiscript RT buffer were mixed in 1:4 ratio in a 20 µl reaction that included 14 µl of RNA template (post-DNA elimination).
- 4 -cDNA was then amplified in a 25 µl PCR reaction, with primers targeting the Neomycin gene; the Neo F primer (5'-CGCTTGGGTGGAGAGGCTATTCGGCTATGAC 3') and Neo R primer (5'-TCCACCATGATATTCGGCAAGCAGGCATC -3').
- 5 -Control PCR reactions was performed using RNA template without the RT to verify the absence of the DNA from the transformed *B. saltans* cell.
- 6 -PCR products were visualized by gel electrophoresis, with purified plasmid as a positive control for the PCR. Amplified PCR products at the expected size of 550-pb were documented (**Figure 1**).

6.1




Figure 1: Agarose gel electrophoresis image showing the RT-PCR results confirming the Neomycin gene expression in *Bodo saltans* transformants at the expected size of 550- bp. (L) L= 100 bp Ladder(New England Biolabs); lanes (1) Neomycin expression profile in *B. saltans* cells transformed with eGFP-Neo-PFRplasmid; (2) Control reaction using the same extracted RNA from transfected *B. saltans* without the RT step to verify the absence of DNA after the DNase treatments(3) Wild-type RT-PCR confirm the absence of the Neomycin expression in wild type samples (4) Control reaction using the same extracted RNA from wild-type *B. saltans* without the RT step to verify the absence of DNA after the DNase treatments (4) PCR positive PCR control using the eGFP-Neo-PFR plasmid DNA; (5) PCR negative control. Forward and reverse PCR primers targeting the Neomycin gene are published in Faktorová et al. (2020). The PCR products were visualized under UV light using 1% agarose gel.

7