



May 11, 2021

Protocol for preparing seagrass-microbiome samples for DNA extraction

Cassie Ettinger¹¹University of California, Riverside**1** Works for me dx.doi.org/10.17504/protocols.io.butqnwmw

EisenLab

Early Career Ocean Scientists

1 more workspace



Cassie Ettinger

University of California, Riverside

ABSTRACT

This protocol details methods for preparing seagrass-associated microbiome samples for DNA extraction.

ATTACHMENTS

Protocol for preparing
samples for DNA
extraction.docx

DOI

dx.doi.org/10.17504/protocols.io.butqnwmw

PROTOCOL CITATION

Cassie Ettinger 2021. Protocol for preparing seagrass-microbiome samples for DNA extraction.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.butqnwmw>



KEYWORDS

Seagrass, Microbiome, DNA extraction

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 07, 2021

LAST MODIFIED

May 11, 2021

OWNERSHIP HISTORY

May 07, 2021  Urmilas

May 11, 2021  Cassie Ettinger University of California, Riverside

PROTOCOL INTEGER ID

49744

Materials for Epiphyte washes:**Redford buffer solution:**

A	B
Tris-HCL	1 M
NaEDTA	0.5 M
Triton-X	1.2%

Mo Bio solution C1

Materials for Endophyte sterilization:

95% EtOH

0.5% NaOCl

70% EtOH

Autoclaved H₂O**Weighing samples:**

1

Before proceeding through wash steps (if necessary), samples should be weighed (this is important if you are considering quantitative PCR or want to standardize sample amounts); prepare balance for samples by cleaning with 70% ethanol solution and obtain sterile weigh boats.

Place weigh boat on balance and tare.

2 Place sample in weigh boat using sterile forceps or scoopulas.

Flame sterilize between samples; procedure: immerse in 1) 10% bleach (0.5% NaOCl) 2) water 3) ethanol 4) flame) – attempt to get approx. constant mass across different sample types – e.g. for sediment ~ **250 mg**.

3 Record sample weight.

4 Place sample in appropriate tube to proceed-sediment directly into DNA extraction tubes, other tissues into tubes to be washed.

Plant tissues (roots, leaves, rhizomes) can now be directly placed in tube for DNA extraction or can be further processed to get epiphyte washes and/or then further surface-sterilized in the following steps.

Epiphyte washes

25m

5m

5 

Using a 1:50 dilution (in autoclaved H₂O) of Redford buffer solution [**1 Molarity (M)** Tris-HCL, **0.5 Molarity (M)** NaEDTA, 1.2% Ctab] from Kembel et al 2014; I used 1.2% Triton-X instead of Ctab as in Kadivar & Stapleton 2003

Add **500 µl** of diluted Redford buffer to samples for **00:05:00** and vortex/agitate.

6 

20m

Pellet in centrifuge @ **4000 x g** for **00:20:00**.

7 Resuspend pellet in MoBio solution C1 (or equivalent) and proceed with DNA extraction.

Endophyte sterilization

5m 5s

8 

Using the surface cleaning method described in Arnold et al 2007 (Mycologica); complete after washing for epiphytes – or – if epiphytes are not wanted, you can proceed here

Rinse tissue with autoclaved H₂O.

9 Immerse tissue in 95% EtOH (**00:00:05**).

5s

10 Immerse tissue in 0.5% NaOCl (~10% bleach) (**00:02:00**).

2m

11 Immerse tissue in 70% EtOH (**00:02:00**).

2m

12 

1m

Rinse by with autoclaved H₂O (**00:01:00**).

13 Repeat step 5 until tissue has been rinsed three times in H₂O.

We can stop after this step if not proceeding directly to DNA extraction and can freeze tissues.

14 Crush tissue (using sterile centrifuge pestle) & proceed with DNA extraction.

For hardy plant tissues, you can sterilize an electric drill and use the centrifuge pestle as a drill bit to crush tissues – just make sure to clean drill in between samples.