



Sep 01, 2022

Proliferation assay

Philippa R Kennedy¹¹University of Minnesota

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.261geojwol47/v1 Philippa R Kennedy
University of Minnesota

ABSTRACT

Assessing the impact of drugs and treatments on natural killer (NK) cell viability and expansion.

DOI

dx.doi.org/10.17504/protocols.io.261geojwol47/v1

PROTOCOL CITATION

Philippa R Kennedy 2022. Proliferation assay. **protocols.io**
<https://protocols.io/view/proliferation-assay-bfcsjiwe>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Apr 20, 2020

LAST MODIFIED

Sep 01, 2022

PROTOCOL INTEGER ID

35954

- 1 PBMCs or enriched NK cells (see *Isolating NK cells from human blood products*) are labeled with a permanent amine-reactive dye (CellTrace Violet Proliferation Kit, Cat. No:C34557, Thermo Fisher) according to the manufacturer's instructions.
- 2 Cells are exposed to various experimental conditions, then harvested after 7 days of culture.

- 3 *Optional:* If required, cells can be stained for apoptotic markers (FITC Annexin V, Cat. No: 556419, BD Biosciences) according to the manufacturer's instructions prior to antibody staining.
- 4 Cells are stained with a dead cell marker (LIVE/DEAD NEAR-IR, Cat. No: L34976, Thermo Fisher), PE-CY7 conjugated anti-CD56 (HCD56, BioLegend), and PE-CF594 conjugated anti-CD3 (UCHT1, BD BioSciences).
- 5 The samples are run on a flow cytometer (LSRII, BD Biosciences) and analyzed using FlowJo software (Tree Star Inc., RRID:SCR_008520).

The proliferation of live (dead cell marker- and annexin V-) NK cells (CD56+/CD3-) and T cells (CD56-/CD3+) is assessed by dilution of the dye.

Relative numbers of live NK cells and T cells for each treatment condition are obtained by running the flow cytometer at a constant flow rate and analyzing the number of each cell type obtained within a specified time limit e.g. 60s.