

2 ▼

Nov 05, 2021

# Recombinant protein expression and purification of codon-optimized Pfu-Sso7d V.2

Maira Rivera<sup>1</sup>, Javiera Reyes<sup>1</sup>, Paula Blazquez-Sanchez<sup>1</sup>,  
Cesar A Ramirez-Sarmiento<sup>1</sup>

<sup>1</sup>Institute for Biological and Medical Engineering, Pontificia Universidad Católica de Chile

1



[dx.doi.org/10.17504/protocols.io.bzusp6we](https://dx.doi.org/10.17504/protocols.io.bzusp6we)

**Reclone.org (The Reagent Collaboration Network)**

Tech. support email: [protocols@recode.org](mailto:protocols@recode.org)

**Click here to message tech. support**



Cesar A Ramirez-Sarmiento  
Pontificia Universidad Catolica de Chile

This protocol has been optimized for the recombinant expression of a codon-optimized Pfu-Sso7d DNA polymerase. This is a fusion protein composed of the Pfu enzyme from *Pyrococcus furiosus* for DNA amplification by PCR fused to a small 7 kDa protein from *Sulfolobus solfataricus* that binds to double-stranded DNA without any preference for specific sequences, thus enhancing polymerization processivity without affecting the catalytic activity or thermal stability of the enzyme.

The goal of this protocol was to eliminate the use of large volumes for dialysis and potential issues with the protein crashing out of the solution due to the use of concentrators for buffer exchange of this enzyme into storage conditions. We also eliminated the use of DTT, which is often found in other similar protocols.

The sequence plasmid encoding the codon-optimized Pfu-Sso7d enzyme used here can be found at <https://benchling.com/s/seq-2TcUPjO2uMbDG5ufTQN4>

DOI

[dx.doi.org/10.17504/protocols.io.bzusp6we](https://dx.doi.org/10.17504/protocols.io.bzusp6we)

Maira Rivera, Javiera Reyes, Paula Blazquez-Sanchez, Cesar A Ramirez-Sarmiento 2021. Recombinant protein expression and purification of codon-optimized Pfu-Sso7d . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bzusp6we>  
Cesar A Ramirez-Sarmiento



## ANID Millennium Science Initiative Program

Grant ID: ICN17\_022

## ANID CONCYTEC

Grant ID: covbio0012

\_\_\_\_\_ protocol ,

Nov 05, 2021

Nov 05, 2021

54898

### MATERIALS

[Sodium phosphate monobasic monohydrate](#) **Sigma**

**Aldrich Catalog #S9638**

[PMSF](#) **Sigma**

**Aldrich Catalog #P7626**

[Sodium phosphate dibasic](#) **Sigma**

**Aldrich Catalog #7558-79-4**

[Imidazole](#) **Sigma Catalog #I5513**

[NaCl](#) **Sigma**

**Aldrich Catalog #53014**

[HiTrap Heparin HP affinity column](#) **Ge Life**

**Sciences Catalog #17040701**

[HisTrap FF Crude Column](#) **Ge**

**Healthcare Catalog #17528601**

[Dextrose](#) **Sigma –**

**Aldrich Catalog #D9434**

[Nonidet P40 Substitute](#) **Sigma**

**Aldrich Catalog # 74385**

[EDTA](#) **Sigma**

**Aldrich Catalog #ED2SS**

### Buffer A, pH 8.0

[M]**50 millimolar (mM)** NaPO<sub>4</sub>, pH 8.0

[M]**50 millimolar (mM)** dextrose

[M]**300 millimolar (mM)** NaCl  
 [M]**1 millimolar (mM)** EDTA  
 [M]**0.1 % volume** Nonidet P-40  
 [M]**40 millimolar (mM)** Imidazole, pH 8.0

#### Buffer B, pH 8.0

[M]**50 millimolar (mM)** NaPO<sub>4</sub>, pH 8.0  
 [M]**50 millimolar (mM)** dextrose  
 [M]**300 millimolar (mM)** NaCl  
 [M]**1 millimolar (mM)** EDTA  
 [M]**0.1 % volume** Nonidet P-40  
 [M]**150 millimolar (mM)** Imidazole, pH 8.0

#### Buffer C , pH 8.0

[M]**50 millimolar (mM)** NaPO<sub>4</sub>, pH 8.0  
 [M]**50 millimolar (mM)** dextrose  
 [M]**300 millimolar (mM)** NaCl  
 [M]**1 millimolar (mM)** EDTA  
 [M]**0.1 % volume** Nonidet P-40  
 [M]**500 millimolar (mM)** Imidazole, pH 8.0

#### Buffer HA, pH 8.0

[M]**50 millimolar (mM)** Tris-HCl, pH 8.0  
 [M]**100 millimolar (mM)** NaCl  
 [M]**0.1 % volume** Nonidet P-40

#### Buffer HB, pH 8.0

[M]**50 millimolar (mM)** Tris-HCl, pH 8.0  
 [M]**2000 millimolar (mM)** NaCl  
 [M]**0.1 % volume** Nonidet P-40

### DAY 1 – Plasmid transformation

1d


- 1 Transform  **100 ng** of plasmid containing codon-optimized into *E. coli* C41 competent cells<sup>2h</sup> using either heat shock or electroporation.
- 2 Spread transformed cells in LB Agar plates supplemented with [M]**0.05 mg/mL** Kan. Grow<sup>12h</sup>

plate overnight at **37 °C**

#### DAY 2 – Preinoculum 1d

- 3 Select a single colony from the LB agar plate to prepare a preinoculum in **10 mL** LB media<sup>1d</sup> supplemented with **0.05 mg/mL** Kan. Grow overnight at **250 rpm, 37°C**

#### DAY 3 – Protein Overexpression 1d

- 4 Use the full volume of the preinoculum to inoculate **1 L** of LB (or TB) media supplemented<sup>4h</sup> with **0.05 mg/mL** Kan (1% inoculation). Grow at **160 rpm, 37°C** until reaching an optical density at 600 nm ( $OD_{600}$ ) = 0.8
- 5 Upon reaching  $OD_{600}$  = 0.8, add IPTG to a final concentration of **0.5 millimolar (mM)**<sup>16h</sup> and incubate overnight at **160 rpm, 18°C**

#### DAY 4A – Protein Purification by IMAC 4h

- 6 Centrifuge the cell culture **4000 x g, 4°C, 00:20:00**. Then, resuspend the cell pellet in<sup>30m</sup> **50 mL** of **Buffer A** freshly supplemented with **0.5 millimolar (mM)** PMSF and **0.2 mg/mL** lysozyme.
- 7 Incubate the resuspended cells at **80 rpm, Room temperature, 00:20:00**<sup>30m</sup>
- 8 Sonicate on ice for **00:04:00** using cycles of **00:00:01** ON and **00:00:04** OFF at<sup>10m</sup> 40% amplitude (Qsonica Q125, 125W).
- 9 On an ultracentrifugation tube, Incubate the unclarified lysate at **70 °C** for **00:30:00**<sup>1h</sup> to precipitate most of *E. coli* proteins, and then place on ice for **00:05:00**. Centrifuge **20000 x g, 4°C, 00:20:00** and collect the supernatant. You might want to collect a small sample for SDS-PAGE afterwards.
- 10 On a **1 mL HisTrap column** preequilibrated with 10 column volumes (c.v.) (here, 10 mL) of **Buffer A**, load the supernatant.<sup>1h</sup>  
Wash with 10-20 c.v. of **Buffer A**.  
Then, elute with 5 c.v. of **Buffer B**, collecting the eluted fractions every **0.5 mL** in 1.5 mL tubes.

- 11 To quickly pool the fractions containing the protein of interest, prepare a 96-well plate or 1.5<sup>5m</sup> mL tubes with **40 µL** of Bradford reagent and **150 µL** of distilled water. Then, add **10 µL** of each protein fraction and compare against a blank reference sample corresponding to **10 µL** of **Buffer B**. You can determine your protein-containing fractions either by absorbance at 595 nm on a plate reader or visually by comparing the blue coloration of each fraction against the blank reference. Pool your fractions and collect a **10 µL** sample for SDS-PAGE


#### DAY 4B – Second purification and buffer exchange by Heparin

2h

- 12 This method was preferred over protein dialysis or Amicon protein concentration to avoid using large buffer volumes and proteins crashing out of the solution. <sup>5m</sup>
- Dilute the pooled fractions 3X in 50 mM Tris-HCl pH 8.0, such that the final concentration of NaCl is 100 mM.
- 13 Next, load the diluted sample onto a **1 ml HiTrap Heparin column** previously equilibrated <sup>30m</sup> with 10 c.v. (here, 10 mL) of **Buffer HA**. Then, elute the protein using a 10 c.v. linear gradient against **Buffer HB**, collecting the eluted fractions every **1 mL** in 1.5 mL tubes. The protein will elute at high concentrations between 300 and 600 mM NaCl
- This linear gradient can be achieved by connecting two containers, one with 5 c.v. **Buffer HA** and the other with 5 c.v. **Buffer HB**, with a syphon or a tube, and withdrawing solution from **Buffer HA** container to the column using a cheap peristaltic pump or by gravity.
- 14 Again, determine your protein-containing fractions using the Bradford assay. Pool your fractions and determine its protein concentration using the same method and collect a **10 µL** sample for SDS-PAGE. <sup>5m</sup>
- 15 For storage, supplement your pooled fraction with **0.2 % volume** Nonidet P-40 and **0.2 millimolar (mM)** EDTA. Then, dilute the sample by adding glycerol up to **50 % volume** to reach final storage conditions: **25 mM Tris-HCl pH 8.0, ~250 mM NaCl, 0.1 mM EDTA, 0.1% Nonidet P-40, 50% glycerol**. <sup>5m</sup>
- 16 Generate **200 µL** aliquots of the enzyme and store it at **-20 °C** until required. Usual final protein concentrations for storage are between **0.2 mg/mL** and **0.6 mg/mL**. <sup>30m</sup>

#### SDS-PAGE Result

17



The unlabelled lane corresponds to the PageRuler protein ladder. His = Pooled fraction from HisTrap IMAC column (20 µl), 1-2-3 correspond to eluted fractions from Heparin column (5 µl).