



AUG 11, 2023

OPEN ACCESS



**Protocol Citation:** Jonathan Breiter 2023. Free-floating Mouse Brain Immunohistochemistry. [protocols.io](https://protocols.io/view/free-floating-mouse-brain-immunohistochemistry-cyeaxtae) <https://protocols.io/view/free-floating-mouse-brain-immunohistochemistry-cyeaxtae>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** In development  
We are still developing and optimizing this protocol

**Created:** Aug 08, 2023

**Last Modified:** Aug 11, 2023

**PROTOCOL integer ID:** 86178

## Free-floating Mouse Brain Immunohistochemistry

Jonathan Breiter<sup>1</sup>

<sup>1</sup>University of Cambridge



Jonathan Breiter  
University of Cambridge


### ABSTRACT

This protocol enables immunohistochemical staining of murine tissue with superior penetration of the tissue by the reagents due to the free-floating approach.

### MATERIALS

10x PBS pH 7.4, Gibco, Cat. No: 70011-036  
Triton X-100, Merck, Cat. No: 648466-50ML  
Bovine Serum Albumin heat shock fraction, Sigma-Aldrich, Cat. No: A9647 - 100G  
Normal Goat Serum, Fisher Scientific, Cat. No: 11819220  
Vectashield Antifade Mounting Medium PLUS, Vector Labs, Cat. No: H-1900  
Mouse mAb to Alpha-synuclein pS129 (81A), Abcam, Cat. No: ab184674, RRID:AB\_2819037  
Donkey anti-Mouse IgG (H+L) Alexa Fluor 568, Thermo Fisher, Cat. No: A10037, RRID: AB\_2534013  
Millex Filter Unit 0.22 um, Merck, Cat. No: SLGP033RS  
12 Well Cell Culture Plate, Corning, Cat. No: 3513  
24 Well Cell Culture Plate, Thermo Scientific, Cat. No: 144530  
Netwell Permeable Supports 15mm Diameter Insert 74 um Polyester Mesh, Costar, Cat. No: 3477  
Micro Slides Single Frosted 75 x 25 mm, Corning, Cat. No: 2948-75X25  
Cover Glass 22 x 50 mm Thickness No. 1, VWR, Cat. No: 631-0137


## Tissue Preparation

- 1 Remove PFA-fixed tissue from storage solution and add to mesh bottom netwell insert inside of 12 well plate that is filled with 0.22  $\mu$ m filtered 1x PBS.
- 2 Place the 12 well plate with the tissue on a horizontal shaker and wash tissue 3x 5 min at approx. 150 rpm, moving the netwell insert to the next well down after each 5 min period to immerse the tissue in new 1x PBS.  Room temperature



## Buffer Preparation


- 3 Per well of tissue make a minimum of 3.5 mL of blocking solution, consider this is needed for blocking, primary and secondary antibody incubation. All reagents should be 0.22  $\mu$ m filtered:

1x PBS  
5% Normal Goat Serum  
2.5% Bovine Serum Albumin  
0.2% Triton-X




Make in excess and keep on ice. Keep at 4 degrees celcius overnight.  4 °C

## Primary Antibody Incubation


- 4 In a new 12 well plate, add 2-3 mL per well of blocking solution and transfer the washed tissue sections inside their netwell inserts into these wells. Incubate for 1h - 2.5 h on horizontal shaker at approx 150 rpm.  Room temperature
- 5 Meanwhile, dilute primary antibodies in blocking solution to appropriate concentrations and keep  On ice .  
  
E.g.  
- Mouse monoclonal anti pS129 (81A) (ab184674) @ 1:750
- 6 Add approx 250 ul of primary antibody solution to the appropriate number of wells of a 24 well plate for the number of brain sections.

- 7 When blocking is finished, move brain tissue sections from the netwell inserts to their appropriate primary antibody well using a fine paintbrush, being careful not to destroy the tissue.
- 8 Incubate on horizontal shaker at approx. 150 rpm overnight.  4 °C

## Secondary Antibody Incubation

- 9 Prepare a new 12 well plate with clean netwell inserts and fill all wells with 0.22 µm filtered 1x PBS.
- 10 Transfer brain tissue sections from 24 well plate into netwell inserts inside 12 well plate and wash 4x 10 minutes on a horizontal shaker at approx. 150 rpm.  Room temperature
- 11 Meanwhile, make appropriate dilutions of secondary antibody in blocking solution. Keep away from light and keep  On ice .  
  
E.g.  
- Donkey anti-Mouse Alexa Fluor 568 (A10037) @ 1:500
- 12 Add approx 250 ul of secondary antibody solution to the appropriate number of wells of a 24 well plate for the number of brain sections.
- 13 When washing step is finished, move brain tissue sections from the netwell inserts to their appropriate secondary antibody well using a fine paintbrush, being careful not to destroy the tissue.  
  
Keep the well plates covered from now to avoid bleaching of fluorophores.
- 14 Incubate 24 well plate on a horizontal shaker at room temperature at approx. 150 rpm for 1 h - 2.5 h.  Room temperature

- 15** Transfer brain tissue sections from 24 well plate into netwell inserts inside 12 well plate and wash 4x 10 minutes on a horizontal shaker at room temperature at approx. 150 rpm.

 Room temperature

## Microscope slide preparation & Imaging

- 16** Plasma-clean microscope cover slips in Argon plasma for 15 minutes.
- 17** Using a fine paintbrush, transfer brain tissue sections from 1x PBS onto a microscope slide, using excess 1x PBS to mount multiple sections next to each other and making sure they are not folded over.
- 18** Leave sections on microscope slides to dry out in the dark (they will turn white).
- 19** Add approx. 150  $\mu$ L of mounting media on top of each tissue section and apply the plasma-cleaned cover slip on top of the tissue, closing the slide.
- 20** Leave the mounting media to dry in the dark.
- 21** Seal the edges of the coverslip with nail varnish and let dry for approx. 30 minutes in the dark.
- 22** Immediately take finished microscope slides to imaging, avoiding unnecessary light exposure.