



MAR 01, 2024

🌐 SOP and tick-sheet for using TDP-43 RNA aptamer (TDP-43APT) to detect pathological TDP-43 in FFPE-preserved human tissue, as described in Spence and Waldron et al., 2024 (Acta Neuropathologica)

Holly Spence¹, Fergal M Waldron¹, Jenna Gregory¹

¹University of Aberdeen

Gregory Lab

OPEN ACCESS



Jenna Gregory

DISCLAIMER

Authors declare no conflicts of interest.



DOI:

dx.doi.org/10.17504/protocols.io.eq2lyjo4mlx9/v1

Protocol Citation: Holly Spence, Fergal M Waldron, Jenna Gregory 2024. SOP and tick-sheet for using TDP-43 RNA aptamer (TDP-43APT) to detect pathological TDP-43 in FFPE-preserved human tissue, as described in Spence and Waldron et al., 2024 (Acta Neuropathologica). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2lyjo4mlx9/v1>

MANUSCRIPT CITATION:

RNA aptamer reveals nuclear TDP-43 pathology is an early aggregation event that coincides with STMN-2 cryptic splicing and precedes clinical manifestation in ALS.

Holly Spence*, Fergal M. Waldron*, Rebecca S. Saleeb, Anna-Leigh Brown, Olivia M. Rifai, Martina Gilodi, Fiona Read, Kristine Roberts, Gillian Milne, Debbie Wilkinson, Judi O'Shaughnessy, Annalisa Pastore, Pietro Fratta, Neil Shneider, Gian Gaetano Tartaglia, Elsa Zacco, Mathew H. Horrocks, Jenna M. Gregory (2024). *Acta Neuropathologica* (in press at the time of release of this SOP on 01/03/2024).

*equal contributions,

‡corresponding author

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: Mar 01, 2024

Last Modified: Mar 01, 2024

PROTOCOL integer ID: 96023

Keywords: ALS, MND, FTD, ALSFTSD, TDP-43, RNA aptamer, Human tissue staining, Immunohistochemistry

Funders Acknowledgement:

Target ALS

Grant ID: BB-2022-C4-L2

ABSTRACT

Here we provide a SOP to outline the correct procedures for performing Immunohistochemistry (IHC) to detect pathological TDP-43 in FFPE-preserved human tissue using the TDP-43 RNA Aptamer (TDP-43^{APT}), as described in Spence and Waldron et al., 2024 published in *Acta Neuropathologica* (in press at the time of release of this SOP on

01/03/2024)

Users with access to Sequenza immunostaining racks and histological facilities (with fume hood) should be able to carry out all steps over two days.

This protocol uses the TDP-43^{APT} sequence published in Zacco *et al.*, 2022. The sequence is: CGGUGUUGCU with a 3' Biotin-TEG modification, purified using HPLC, scale: 1.0 µM synthesis.

Reference for citations of this method

RNA aptamer reveals nuclear TDP-43 pathology is an early aggregation event that coincides with STMN-2 cryptic splicing and precedes clinical manifestation in ALS.

Holly Spence*, Fergal M. Waldron*, Rebecca S. Saleeb, Anna-Leigh Brown, Olivia M. Rifai, Martina Gilodi, Fiona Read, Kristine, Roberts, Gillian Milne, Debbie Wilkinson, Judi O'Shaughnessy, Annalisa Pastore, Pietro Fratta, Neil Shneider, Gian Gaetano Tartaglia, Elsa Zacco, Mathew H. Horrocks, Jenna M. Gregory (2024). *Acta Neuropathologica* (in press at the time of release of this SOP on 01/03/2024).

*equal contributions, ‡corresponding author

ATTACHMENTS

[TDP-43 RNA Aptamer_IHC_SOP_v1.0_01Mar24.pdf](#)

[TDP-43 RNA Aptamer_IHC_TS_v1.0_01Mar24.pdf](#)

IMAGE ATTRIBUTION

Jenna M Gregory

GUIDELINES

This protocol is designed for users with access to Sequenza immunostaining racks and histological facilities (with fume hood) should be able to carry out all steps over two days.

BEFORE START

Please see appendix for materials required to carry out this protocol.

MATERIALS

See SOP appendix for materials.

SAFETY WARNINGS



Safety First

Before starting, please ensure all relevant Health & Safety documentation is in order including the following

- COSHH assessment
- Risk assessment
- Safe System of Work
- SOP read and understood

BEFORE START INSTRUCTIONS

Please see appendix for materials required to carry out this protocol.