



Version 2

Nov 10, 2020

Stable transfection of unicellular relative of animals, *Corallochytrium limacisporum*, using Lonza Nucleofector V.2

Aleksandra Kozyczkowska¹, Sebastian R Najle², Eduard Ocaña-Pallarès², Cristina Aresté³, Iñaki Ruiz-Trillo², Elena Casacuberta⁴

¹Institut de Biologia Evolutiva (UPF-CSIC) Barcelona; ²Institut de Biologia Evolutiva (CSIC-Universitat Pompeu Fabra); ³IIBB;

⁴Institute of Evolutionary Biology

1

Works for me

This protocol is published without a DOI.

Protist Research to Optimize Tools in Genetics (PROT-G)

Multicellgenomelab

1 more workspace



Aleksandra Kozyczkowska

ABSTRACT

This protocol describes the transfection of *Corallochytrium limacisporum*, one of the unicellular relatives of animals. The transfection is performed using Lonza Nucleofector System X and reporter cassettes that contain endogenous regulatory elements of *C. limacisporum*.

PROTOCOL CITATION

Aleksandra Kozyczkowska, Sebastian R Najle, Eduard Ocaña-Pallarès, Cristina Aresté, Iñaki Ruiz-Trillo, Elena Casacuberta 2020. Stable transfection of unicellular relative of animals, *Corallochytrium limacisporum*, using Lonza Nucleofector. **protocols.io**

<https://protocols.io/view/stable-transfection-of-unicellular-relative-of-ani-bpfgmjw>

Version created by Aleksandra Kozyczkowska

KEYWORDS

transfection, electroporation, unicellular organism, stable

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 05, 2020

LAST MODIFIED

Nov 10, 2020

PROTOCOL INTEGER ID

44232

GUIDELINES

Keeps cells on ice during the whole procedure & all reagents must be ice-cold

MATERIALS TEXT

Main materials:

P3 PrimaryCell 4D-Nucleofector Kit S (32 rxn) Kit Catalog# H3V4XP-5032

pUC19 carrier plasmid

highly concentrated reporter construct

1XPBS

Marine Broth

ABSTRACT

This protocol describes the transfection of *Corallochytrium limacisporum*, one of the unicellular relatives of animals. The transfection is performed using Lonza Nucleofector System X and reporter cassettes that contain endogenous regulatory elements of *C. limacisporum*.

Cell Culture preparation

- 1 Inoculate 100 µl of *C. limacisporum* culture in 5ml of Marine Broth in 25cm² flask. Keep at 23°C air incubator for two days.

Cell count

- 2 Scrape gently and count cells to reach 1.5E6 cells / condition using hemocytometer
If there are more conditions, it is recommended to pull the cells together to get more visible pellet

Medium removal

- 3 Spin down the cells at 1500xg for 5'
Discard the supernatant.

Washing step

5m

- 4 Wash the pellet with an ice-cold 1XPBS (100 µl / condition)
Spin down and discard the supernatant.

Prepare DNA mix

- 5 Mix DNA plasmids: 10 µg of an reporter plasmid + 20 µg carrier DNA (pUC19)
It is recommended that plasmid is highly concentrated (2-5 µg/ul) to avoid final buffer dilution with the DNA solution
Our plasmid contains mCherry coding sequence & resistance to puromycin gene (*pac*) under constitutive promoter: actin (Addgene: #104446) or tubulin (Addgene#129560)

Note: 10 µg of reporter plasmids guarantees the highest efficiency, however a successful transfection can be also obtained with 5 and 1 µg of the plasmid

- 5.1 Transfer DNA mix to 20 µl of ice-cold P3 buffer
It is recommended not to keep the cells in P3 buffer for too long

- 6 Add DNA mix to the cell pellet

- 7 Transfer the mix (cells + DNA plasmid) into a strip (from the manufacture kit)
Do it carefully & without creating bubbles
Final volume: 20-25 µl

Transfection of plasmid DNA into *C. limacisporum*

- 8 Insert into Lonza nucleofection system and apply code: EN - 138

Recovery step

- 9 After transfection, immediately add 80 μ l of Marine Broth medium (MB) to a well
Mix up and down

9.1 Transfer to 12-well plate (NUNC) that contains 1 ml of growth medium (MB) and incubate overnight

Check-out for results

- 10 Check for positive cells 24 - 48h post-transfection using fluorescent microscopy

Enrichment & selection of transfected cells

- 11 Replace 500 μ l of medium with fresh MB that contains 300 μ g/ml of puromycin
Increase of mCherry-positive cells can be observed over time.