



Sep 23, 2022

Genomic DNA isolation from fixed cells

Herschel Dhekne¹, Ebsy Jaimon¹, Suzanne R Pfeffer¹¹Department of Biochemistry, Stanford University

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.eq2lynm9qvx9/v1

Suzanne R Pfeffer

ABSTRACT

This protocol details the procedure of genomic DNA isolation from fixed cells.

ATTACHMENTS

[iugzbc32f.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.eq2lynm9qvx9/v1

PROTOCOL CITATION

Herschel Dhekne, Ebsy Jaimon, Suzanne R Pfeffer 2022. Genomic DNA isolation from fixed cells. **protocols.io**
<https://protocols.io/view/genomic-dna-isolation-from-fixed-cells-cfgptjvn>

FUNDERS ACKNOWLEDGEMENT



Aligning Science Across Parkinson's Disease
Grant ID: ASAP-000463

KEYWORDS

DNA isolation, Proteinase K, Econospin, ASAPCRN

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 18, 2022

LAST MODIFIED

Sep 23, 2022

OWNERSHIP HISTORY

Aug 18, 2022  madhavi.d

Aug 23, 2022  Suzanne R Pfeffer

PROTOCOL INTEGER ID

68847

MATERIALS TEXT

Reagents:

- Proteinase K (Zymo #D3001-2-20) (20mg/ml stock in storage buffer).
- RNase A (Sigma #70856) (10mg/ml stock).
- 100% Molecular biology grade Ethanol.
- Silica spin columns (Econospin).
- Qiagen buffer AL.
- Qiagen buffer AW1.
- Qiagen buffer AW2.
- Nuclease free water.

Buffer AL (storage: room temperature (RT))

A	B
Tris-HCl (7.4)	50 mM
Guanidine HCl	5.5 M
EDTA	20 mM
Triton X-100	1.3%

Buffer AW1 (storage: RT)

A	B
Guanidine HCl	1 M
EtOH, pH 5.5	57%

Buffer AW2 / PE (storage: RT)



A	B
Tris-HCl (pH 7.5)	10 mM
Ethanol	80%

Buffer AW2 / PE (storage: RT)

A	B
Tris-HCl (pH 9)	10 mM
EDTA	0.5 mM

Genomic DNA isolation from fixed cells

38m 5s

- 1 Resuspend 1-3 X10⁶ cells in  **250 µL** of PBS and transfer to a 2 mL tube (this would represent one  **3.5 cm** dish of 3T3 cells).

2



Add **200 µg** Proteinase K (from a **20 mg/mL** stock) and **200 µg** RNase A (from a **20 mg/mL** stock).

3



30m

Incubate cells at **37 °C** for **00:30:00** in a water bath.

4



Add **250 µL** Qiagen AL lysis buffer per **250 µL** of the protease and RNAase-containing cell suspension and mix thoroughly.

5



Place tubes in an incubator at **56 °C** with shaking at **800 rpm** **Overnight**, capped.

Each tube is parafilm sealed to ensure safety.

6



5s

Add **250 µL**, 100% molecular biology grade ethanol; mix slowly using a slow vortex for **00:00:05**.

7

With a razor blade, trim the tip of a **1 mL** pipet tip to enlarge the opening. Use this tip to pipet out DNA from the ethanol solution and apply it onto a silica spin DNA binding column (e.g. EconoSpin 1920-250).

8






1m

Spin at **6000 x g** for **00:01:00** in a fixed angle tabletop microfuge; aspirate and discard flow-through.




9





1m

Add  **500 µL** Qiagen **buffer AW1** to the column, spin again at  **6000 x g** for  **00:01:00** , aspirate and discard flow-through.

10   1m

Add  **500 µL** Qiagen **buffer AW2**, spin at  **8000 x g** for  **00:01:00** , aspirate and discard flow-through.


11  1m


Spin once more using microfuge to remove excess ethanol at  **13000 x g** for  **00:01:00** .

12 

Transfer column into a new, 1.5 mL collection tube.

13 Elute with pre-warmed,  **100 µL** nuclease-free **water** or **TE**.

Volume depends on starting number of cells: use  **100 µL** per 1 million cells.


14   2m

Incubate for  **00:01:00** , then spin at  **13000 x g** for  **00:01:00** .

15 


Add another  **50 µL** **nuclease free water** to accomplish a second elution.

16   2m

Incubate at  **Room temperature** for  **00:01:00** , then spin as before at  **13000 x g** for  **00:01:00** .

17 **The two flow-through fractions contain the genomic DNA.**

18 Perform Nanodrop and Qubit HS DNA estimation to calculate yield.

Theoretically, 1×10^6 cells should yield  **6 µg** DNA.