



SOLUTION- 05 - Wash solution (PBS/FBS) for PBMC

Marco Cosentino¹, Elisa Storelli¹, Alessandra Luini¹, Massimiliano LM Legnaro¹, Emanuela Rasini¹, Marco Ferrari¹, Franca Marino¹

¹Center for Research in Medical Pharmacology, University of Insubria (Varese, Italy)

1 Works for me dx.doi.org/10.17504/protocols.io.biankade



Elisa Storelli
Center for Research in Medical Pharmacology, University of I...

DOI

dx.doi.org/10.17504/protocols.io.biankade

DOCUMENT CITATION

Marco Cosentino, Elisa Storelli, Alessandra Luini, Massimiliano LM Legnaro, Emanuela Rasini, Marco Ferrari, Franca Marino 2020. SOLUTION- 05 - Wash solution (PBS/FBS) for PBMC. **protocols.io**
dx.doi.org/10.17504/protocols.io.biankade

LICENSE

————— This is an open access document distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 07, 2020

LAST MODIFIED

Jul 23, 2020

DOCUMENT INTEGER ID

38958

Sterile PBS/FBS 2% (v/v)

Prepare with sterile PBS and FBS (Fetal Bovine Serum). Both can be found in Fridge 2 - (Room PS03).

Prepare 50 ml of PBS/FBS 2% by adding 1 ml of FBS to 49 ml of RPMI, work under the laminar flow hood (Room PS03).

Storage:  **4 °C** Fridge 2 - (Room PS03).

FBS: catalog number ECS0180L, Euroclone, Italy.