

Dec 19, 2024

# Intracardiac perfusion with buffer for anatomical studies

 Forked from [Intracardiac perfusion with fixative for anatomical studies](#)

DOI

**[dx.doi.org/10.17504/protocols.io.3byl4kpp8vo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl4kpp8vo5/v1)**

Janet R Keast<sup>1</sup>, Peregrine B Osborne<sup>1</sup>, Nicole Wiedmann<sup>1</sup>

<sup>1</sup>University of Melbourne

SPARC

Tech. support email: [info@neuinfo.org](mailto:info@neuinfo.org)



Nicole Wiedmann

University of Melbourne

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.3byl4kpp8vo5/v1](https://dx.doi.org/10.17504/protocols.io.3byl4kpp8vo5/v1)**

**Protocol Citation:** Janet R Keast, Peregrine B Osborne, Nicole Wiedmann 2024. Intracardiac perfusion with buffer for anatomical studies. [protocols.io](https://dx.doi.org/10.17504/protocols.io.3byl4kpp8vo5/v1) **<https://dx.doi.org/10.17504/protocols.io.3byl4kpp8vo5/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** November 01, 2021

**Last Modified:** December 19, 2024

**Protocol Integer ID:** 54689

**Keywords:** tissue preservation, tissue fixation

**Funders Acknowledgements:**

NIH SPARC

Grant ID: 3OT2OD023872



## Abstract

This protocol is suitable for preserving tissues for anatomical studies of organs and major pelvic ganglia in adult rats. The protocol is performed under anesthesia and should incorporate all local requirements for standards of animal experimentation.

## Materials

### MATERIALS

⊗ Ketamine **Lyppard Catalog #KETAI1**

⊗ Xylazine **Lyppard Catalog #L10605**

⊗ Heparin sodium **Ellar Laboratories**

⊗ Paraformaldehyde fixative: 4% paraformaldehyde in phosphate buffered saline (PBS) **Contributed by users**

## Preparation for perfusion

- 1 Make up the following solutions:

**Perfusion prewash solution:** To 300 ml 0.9% sodium chloride (w/v) add 3.75 ml 1% sodium nitrite (w/v) and 0.11 ml heparin (5000 IU/ml). This is made up immediately prior to use.

**Fixative:** 4% paraformaldehyde in 0.1M phosphate buffer, pH 7.4. This is made up no longer than 48h prior to use, stored at 4C and brought to room temperature on the day of perfusion.

## Perfusion

- 2 Induce anesthesia by an intraperitoneal injection of ketamine (100 mg/kg) and xylazine (10 mg/kg).
- 3 After opening the chest cavity, inject the left ventricle with mixture of 0.25 ml heparin (5000 IU/ml) and 0.5 ml 1% sodium nitrite.
- 4 A needle connected to tubing (T-connector to prewash and perfusion solutions), and a peristaltic pump (setting: 50 ml/min) is then inserted into the left ventricle and clamped into place using hemostats. Make a small incision in the right atrium to drain blood and perfusate during the procedure.
- 5 Perfuse with pre-wash solution until the fluid flowing from the right atrium is clear, the liver and extremities are pale (typically 2-3 min).
- 6 Dissect tissues required for analysis and place in fixative overnight at 4 degrees.

## Post-perfusion

- 7 Wash tissues with phosphate-buffered saline (PBS; 0.1 M, pH7.2), 3 x 30 min washes.
- 8 Store in PBS containing 0.1% sodium azide at 4C until used for immunohistochemistry or other microscopic analysis.