

Sep 26, 2024

## Bacteria Transformation

DOI

[dx.doi.org/10.17504/protocols.io.n2bvjn3kpgk5/v1](https://dx.doi.org/10.17504/protocols.io.n2bvjn3kpgk5/v1)

Carolina Lopez<sup>1</sup>

<sup>1</sup>Washington University



Cecilia Escudero

wustl

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.n2bvjn3kpgk5/v1](https://dx.doi.org/10.17504/protocols.io.n2bvjn3kpgk5/v1)

**Protocol Citation:** Carolina Lopez 2024. Bacteria Transformation. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n2bvjn3kpgk5/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 28, 2024

**Last Modified:** September 26, 2024

**Protocol Integer ID:** 100763

## Abstract

General protocol for bacteria transformation



## Materials

### **Reagents:**

- Competent E. coli cells (NEB, C3040H)
- LB plates with necessary antibiotic
- 50% Glycerol



## Bacteria Transformation:

1h 24m

- 1
  1. Pipette  $2\ \mu\text{L}$  of product to  $50\ \mu\text{L}$  of E. coli competent cells and pipette slowly up and down 5 times to mix. For plasmids, use  $2\ \mu\text{L}$  ul of  $0.005\ \mu\text{g}/\mu\text{L}$  plasmid solution. For ligation reactions, DNA concentrations will vary so use  $2\ \mu\text{L}$  ul of reaction.
  2. Incubate for 00:20:00 on ice.
  3. Heat shock bacteria in  $42\ ^\circ\text{C}$  waterbath for 00:01:00 .
  4. Incubate on Ice for 00:03:00 .
  5. Add  $100\ \mu\text{L}$  of SOC media and shake in warm room for 01:00:00 .
  6. Pipette  $100\ \mu\text{L}$  of bacteria onto LB-Antibiotic Plate and spread cells out with plate spreader to get individual colonies.
  7. Incubate overnight in  $37\ ^\circ\text{C}$  warm room.
  8. Pick colonies for miniprep growth and sequencing.

## Glycerol Bacteria Stocks:

- 2
  1. In screw top  $2\ \text{mL}$  tube. Mix  $500\ \mu\text{L}$  of Bacteria Solution (Either from Miniprep or Maxiprep before centrifugation) and  $500\ \mu\text{L}$  of 50% (v/v) Glycerol Solution.
  2. Store at  $-80\ ^\circ\text{C}$  .
  3. To pull out bacteria, use bacteria loop to scrap some of still frozen glycerol stock and streak bacteria onto LB-Antibiotic plate. Then pick colonies for growth.