

May 21, 2024

Multiplex Labeling with Tyramide Fluorophores for Detecting CIAP-Resistant PSER129 and Proteinase K-Resistant aSyn in situ (Killinger 2024)

DOI

dx.doi.org/10.17504/protocols.io.x54v92dq4l3e/v1

Solji Choi¹, Bryan Killinger²

¹Rush University Medical Center; ²Rush University

ASAP Collaborative Rese...

Killinger



Bryan Killinger

Rush University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.x54v92dq4l3e/v1

Protocol Citation: Solji Choi, Bryan Killinger 2024. Multiplex Labeling with Tyramide Fluorophores for Detecting CIAP-Resistant PSER129 and Proteinase K-Resistant aSyn in situ (Killinger 2024). protocols.io

<https://dx.doi.org/10.17504/protocols.io.x54v92dq4l3e/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 15, 2024

Last Modified: May 21, 2024

Protocol Integer ID: 100194



Funders Acknowledgement:

Michael J. Fox Foundation

Grant ID: ASAP-021030

NIH-NINDS

Grant ID: 1R01NS128467


Abstract

This protocol aims to examine the association of calf-intestinal alkaline phosphatase (CIAP)-resistant alpha-synuclein phosphorylated at serine 129 (PSER129) and proteinase K (PK)-resistant alpha-synuclein (aSyn) in the mouse brain, particularly in M83 transgenic mice treated with preformed fibrils. M83 lines exhibit a notably higher abundance of endogenous PSER129 compared to wild-type mice.

Materials

▪ Dilution media:

A	B
Tris-HCl, pH 7.4	50 mM
NaCl	150 mM
Triton- X100	0.5%

 Alkaline Phosphatase, Calf Intestinal (20 u/μl) **Promega Catalog #M2825**

▪ CIAP buffer:

A	B
NaCl	100 mM
Tris-HCl	50 mM
MgCl ₂ , pH 7.9	10 mM
Autoclave and store RT	

▪ Blocking buffer:

A	B
Dilution media	100 mL
Normal serum	3 mL
BSA	2 g
Triton X100	0.4 mL
Mix well so the Triton is completely dissolved	

▪ Borate buffer:

A	B
Borate buffer, pH 8.5	0.05 M
DI H ₂ O	300 mL
Sodium tetraborate decahydrate	5.72 g
Mix well to dissolve completely. Adjust to pH 8.5	

▪ Components:

A	B
Borate buffer	10 mL




A	B
H2O2	1 uL
TF	5 uL

▪ **Sodium Citrate Buffer, pH 6.0 (1L):**

A	B
Sodium citrate-Trisodium salt (Dihydrate) in 1000 mL DI water	2.94 g
Tween-20	0.5 mL
Mix well	

▪ **PBS:**

A	B
Tris-HCl, pH 7.2	50 mM
NaCl	158 mM

 Proteinase K **Thermo Fisher Scientific Catalog #EO0491**

Protocol materials

 Proteinase K **Thermo Fisher Scientific Catalog #EO0491** Materials, Step 25

 Alkaline Phosphatase, Calf Intestinal (20 u/μl) **Promega Catalog #M2825** Materials, Step 5



Day 1

1d 1h 10m

- 1 Wash free-floating tissue (3 x 10 minutes) in dilution media (DM).

**Dilution media:**

A	B
Tris-HCl, pH 7.4	50 mM
NaCl	150 mM
Triton- X100	0.5%

- 1.1 Wash free-floating tissue for 00:10:00 in dilution media (DM). (1/3)

10m



- 1.2 Wash free-floating tissue for 00:10:00 in dilution media (DM). (2/3)

10m



- 1.3 Wash free-floating tissue for 00:10:00 in dilution media (DM). (3/3)

10m



- 2 Incubate the samples with 1% Triton X-100 in DM for 00:10:00 .

10m



- 3 Wash in DM for 00:10:00 .

10m



- 4 Wash the tissues in CIAP buffer (2x10 minutes).

**CIAP buffer:**

A	B
NaCl	100 mM
Tris-HCl	50 mM
MgCl ₂ , pH 7.9	10 mM
Autoclave and store RT	

4.1 Wash the tissues in CIAP buffer for 00:10:00 . (1/2)

10m



4.2 Wash the tissues in CIAP buffer for 00:10:00 . (2/2)

10m



5 Incubate the tissues with CIAP at a dilution of 1:333 for 24:00:00 at 37 °C on a shaker.

1d



- Alkaline Phosphatase, Calf Intestinal (20 u/μl) **Promega Catalog #M2825** .
- In 500uL CIAP buffer, add 1.5 μl CIAP (30 units).

Day 2

8h 30m

6 Wash in DM (3 x 10 minutes).



6.1 Wash in DM for 00:10:00 . (1/3)

10m



6.2 Wash in DM for 00:10:00 . (2/3)

10m



6.3 Wash in DM for 00:10:00 . (3/3)

10m



7 Endogenous peroxidase inhibition and serum blocking step (1-hour incubation):

- 0.3% H₂O₂+0.1% Sodium Azide in 50 mL blocking buffer.
- **Blocking buffer:**

A	B
Dilution media	100 mL
Normal serum	3 mL
BSA	2 g
Triton X100	0.4 mL



A	B
Mix well so the Triton is completely dissolved	

8 Dilute primary antibody in blocking buffer. Incubate Overnight at 4 °C .

8h



- Recombinant Anti-Alpha-synuclein (phospho S129) antibody (EP1536Y, ab51253), dilution factor: 1:50K
- In 30 mL blocking buffer, add 0.6 µL PSER129 antibody.

Day 3

17h 20m

9 Wash in DM (3 x 10 minutes).



9.1 Wash in DM for 00:10:00 . (1/3)

10m



9.2 Wash in DM for 00:10:00 . (2/3)

10m



9.3 Wash in DM for 00:10:00 . (3/3)

10m



10 HRP-Secondary antibody incubation 1:1000 dilution (01:00:00).

1h

- Solvent is 100 mL DM + 1 mL normal serum + 1 g BSA

11 Wash in DM (2 x 10 minutes).



11.1 Wash in DM for 00:10:00 . (1/2)

10m

11.2 Wash in DM for 00:10:00 . (2/2)

10m



12 Wash in borate buffer for 00:10:00 .

10m



▪ **Borate buffer:**

A	B
Borate buffer, pH 8.5	0.05 M
DI H2O	300 mL
Sodium tetraborate decahydrate	5.72 g
Mix well to dissolve completely. Adjust to pH 8.5	

13 Incubate with tyramide fluorophore (TF) for 00:30:00 while blocking light. After this step, always protect the tissues from light.

30m



▪ **Components:**

A	B
Borate buffer	10 mL
H2O2	1 μ L
TF	5 μ L

14 Wash in DM (2 x 10 minutes).



14.1 Wash in DM 00:10:00 . (1/2)

10m



14.2 Wash in DM 00:10:00 . (2/2)

10m



15 View under the microscope to confirm successful staining.



16 Heat water bath to 80 °C - 85 °C for 01:30:00 before the primary antibody elution step.

1h 30m

17 Place the dish containing sodium citrate buffer in the water bath and heat it for 00:10:00 .

10m

**▪ Sodium Citrate Buffer, pH 6.0 (1L):**

A	B
Sodium citrate-Trisodium salt (Dihydrate) in 1000 mL DI water	2.94 g
Tween-20	0.5 mL
Mix well	

18 Wash the tissues in sodium citrate buffer for 00:10:00

10m



19 Incubate the tissues in the heated sodium citrate buffer for 00:30:00 .

30m



20 Cool the dish containing tissues to Room temperature (at least 00:20:00).

20m

21 Wash in DM for 10 min x 2 times.



21.1 Wash in DM for 00:10:00 . (1/2)

10m



21.2 Wash in DM for 00:10:00 . (2/2)

10m



22 Mount the tissues on Superfrost Plus Microscope Slides (Fisherbrand), and completely dry it for at least 02:00:00 .

2h

23 Heat water bath to 37 °C for 00:30:00 before the PK digestion step.

30m

24 Place the dish containing PBS in the water bath and heat it for 00:10:00 .

10m

▪ PBS:

A	B
Tris-HCl, pH 7.2	50 mM



A	B
NaCl	158 mM

25 Add the PK to the PBS at a dilution of 1:666 and mix well.



- Proteinase K **Thermo Fisher Scientific Catalog #EO0491**
- 30 mL PBS, add 45 µL PK.

26 Incubate the mounted tissues in the PK containing PBS for 00:30:00 .

30m



27 Wash the slide in PBS (2 x 5 minutes).



27.1 Wash the slide in PBS for 00:05:00 . (1/2)

5m



27.2 Wash the slide in PBS for 00:05:00 . (2/2)

5m



28 Incubate the slide in 4% PFA for 00:30:00 at Room temperature on a shaker.

30m



29 Wash in DM (2 x 5 minutes).



29.1 Wash in DM for 00:05:00 . (1/2)

5m



29.2 Wash in DM for 00:05:00 . (2/2)

5m



30 Block the tissues on slide using Bloxall endogenous blocking solution (Vector Laboratories) for 00:10:00 at Room temperature on a shaker.

10m



31 Dilute primary antibody in blocking buffer. Incubate Overnight at 4 °C .

8h



- Recombinant Anti-Alpha-synuclein antibody (EPR20535, ab212184), dilution factor: 1:20K
- In 30 mL blocking buffer, add 1.5 µL aSyn antibody.

Day 4

3h 50m

32 Wash in DM (3 x 10 minutes).



32.1 Wash in DM for 00:10:00 . (1/3)

10m



32.2 Wash in DM for 00:10:00 . (2/3)

10m



32.3 Wash in DM for 00:10:00 . (3/3)

10m



33 HRP-Secondary antibody incubation 1:1000 dilution (01:00:00).

1h

- Solvent is 100 mL DM + 1 mL normal serum + 1 g BSA

34 Wash in DM (2 x 10 minutes).



34.1 Wash in DM for 00:10:00 . (1/2)

10m



34.2 Wash in DM for 00:10:00 . (2/2)

10m




35 Wash in borate buffer for 00:10:00 .

10m

- **Borate Buffer:**



A	B
Borate buffer, pH 8.5	0.05 M
DI H ₂ O	300 mL
Sodium tetraborate decahydrate	5.72 g
Mix well to dissolve completely. Adjust to pH 8.5	

35.1 Incubate with tyramide fluorophore (TF) for  00:30:00 while blocking light.

30m



35.2 Wash in DM (2 x 10 minutes).

20m



▪ Wash in DM  00:10:00 . (1/2)

▪ Wash in DM  00:10:00 . (2/2)

35.3 **Components:**


A	B
Borate buffer	10 mL
H ₂ O ₂	1 uL
TF	5 uL

36 View under the microscope to confirm successful staining.




37 Wash in PBS (2 x 10 minutes).



37.1 Wash in PBS for  00:10:00 . (1/2)

10m



37.2 Wash in PBS for  00:10:00 . (2/2)

10m






38 Counterstain with DAPI for  00:20:00 at  Room temperature .

20m

- 1:2000 dilution in ddH₂O or PBS.


39 Wash the tissues in PBS (2 x 10 minutes).



39.1 Wash the tissues in PBS for  00:10:00 . (1/2)


10m



39.2 Wash the tissues in PBS for  00:10:00 . (2/2)

10m



40 Cover the slide with Fluoroshield, and coverslip. Seal with nail polish on all sides of the coverslip. Always protect the slides from light. Slides can be stored at  4 °C .