



SEP 05, 2023

OPEN ACCESS



**DOI:**  
[dx.doi.org/10.17504/protocols.io.x54v9p2b1g3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v9p2b1g3e/v1)

**Protocol Citation:** Daniel Manrique-Castano 2023. Staining of CD31 and CD13 in PDGFR-B/Td-tomato brain sections. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v9p2b1g3e/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


**Protocol status:** Working  
 We use this protocol and it's working

**Created:** Sep 05, 2023

## 🌐 Staining of CD31 and CD13 in PDGFR-B/Td-tomato brain sections

Daniel Manrique-Castano<sup>1</sup>

<sup>1</sup>Université Laval



Daniel Manrique-Castano  
 Université Laval

### ABSTRACT

This protocol is suitable for the staining of CD31 and CD13 in PDGFR-B/Td-Tomato fixed mouse brain sections.

### GUIDELINES

Read the entire protocol before starting the procedure.

Note that this protocol uses free-floating sections to enhance CD31 staining.

This protocol uses sequential staining to avoid cross-reactivity between CD31 and CD13.

Do not submit the tissue to antigen retrieval or other acidic solutions to preserve the Td-tomato signal.

Works with the listed antibodies and dilutions. For other references, preliminary tests are highly recommended.

### MATERIALS

Last Modified: Sep 05, 2023

PROTOCOL integer ID: 87394

Keywords: Staining, GFAP, IBA1, Immunofluorescence

Equipment	
ImmEdge® Hydrophobic Barrier PAP Pen	NAME
Vector	BRAND
H-4000	SKU
<a href="https://vectorlabs.com/products/histology/immedge-hydrophobic-barrier-pen">https://vectorlabs.com/products/histology/immedge-hydrophobic-barrier-pen</a>	LINK
Immunohistochemistry / Immunocytochemistry, Immunofluorescence, In situ hybridization	SPECIFICATIONS

Permeabilization solution =

[M] 0.3 % (v/v)

Triton

Triton X-100 Sigma

Aldrich Catalog #T8787-50ML

/

=

[M] 0.1 % (v/v)

TWEEN 20 Sigma

Aldrich Catalog #P7949

/

[M] 0.3 Molarity (m)

Glycine

Contributed by users

in 1x PBS

Antibody buffer =

[M] 0.1 % (v/v)

TWEEN 20 Sigma

Aldrich Catalog #P7949

/

[M] 1 % (v/v)

Normal Donkey Serum Merck MilliporeSigma (Sigma-Aldrich) Catalog #D9663-10ml

[M] 1 Mass / % volume

Bovine Serum Albumin (BSA) Sigma

Aldrich Catalog #A7906

Primary antibodies:

Mouse Aminopeptidase N/CD13 Antibody R&D Systems Catalog #AF2335

Purified Rat Anti-Mouse CD31 BD Biosciences Catalog #550274

Secondary antibodies:

Donkey anti-Goat IgG (H+L) Thermo Fisher Scientific Catalog #A-11055

protocols.io | <https://dx.doi.org/10.17504/protocols.io.x54v9p2b1g3e/v1>

Oct 5 2023

2



Alexa Fluor® 647 AffiniPure Donkey Anti-Rat IgG (H+L) Jackson ImmunoResearch Laboratories, Inc. Catalog #712-605-153




DAPI (46-Diamidino-2-Phenylindole Dilactate) Invitrogen - ThermoFisher Catalog #D3571

#### Mounting media:




Fluoromount-G Electron Microscopy Sciences Catalog #17984-25


## PROTOCOL MATERIALS

 Donkey anti-Goat IgG (H L) Thermo Fisher Scientific Catalog #A-11055


Materials

 Purified Rat Anti-Mouse CD31 BD Biosciences Catalog #550274


Materials

 Alexa Fluor® 647 AffiniPure Donkey Anti-Rat IgG (H L) Jackson ImmunoResearch Laboratories, Inc. Catalog #712-605-153


Materials

 Mouse Aminopeptidase N/CD13 Antibody R&D Systems Catalog #AF2335


Materials

 Fluoromount-G Electron Microscopy Sciences Catalog #17984-25

Materials, Step 12


 Triton X-100 Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML

Materials


 TWEEN 20 Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7949

In Materials, Materials and [7 steps](#)

 Glycine **Contributed by users** Materials

 Normal Donkey Serum Merck MilliporeSigma (Sigma-Aldrich) Catalog #D9663-10ml

Materials

 Bovine Serum Albumin (BSA) Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7906

Materials

## SAFETY WARNINGS





DAPI is highly toxic. Handle it with care.

## Tissue preparation and blocking

20m



- 1 Extract the sections from the anti-freeze media and rinse them in PBS using a 24-well plate.

- 2 Aspirate the PBS and incubate the sections in the **Blocking solution** for  01:00:00 at  Room temperature .

1h

## Antibody incubation

12h

- 3 When blocking is finished, aspirate the buffer (no washing is required) and incubate CD31 (BD Bioscience, 550274, 1:200) in **antibody buffer** for  Overnight at  4 °C .

3h



### Note


Since the sections are already permeabilized, please note that antibody buffers do not contain Triton X-100

Please note that the proposed antibody incubation includes only **CD31**. **CD13** will be incubated after.


- 4 When primary antibody incubation is finished, aspirate the media and wash the sections

5m




 00:05:00 5x

with  0.05 % (v/v)


 TWEEN 20 Sigma  
Aldrich Catalog #P7949

in PBS.

- 5 Incubate **Donkey α Goat 488 (Invitrogen, A-11055, 1:500) secondary antibody** in  0.1 % (v/v)

1h




 TWEEN 20 Sigma  
Aldrich Catalog #P7949


for  01:00:00 at  Room temperature

- 6 When secondary antibody incubation is finished, aspirate the media and wash the sections

5m

 00:05:00 5x

with  0.05 % (v/v)

 TWEEN 20 Sigma  
Aldrich Catalog #P7949

in PBS.

## Tissue blocking

- 7 Perform a second blocking step using **5% Donkey serum** in [M] 0.1 % (v/v)



TWEEN 20 Sigma  
Aldrich Catalog #P7949

in PBS.

## Antibody incubation

12h

- 8 When blocking is finished, aspirate the buffer (no washing is required) and incubate CD13 (R&D, AF2335, 1:300) in **antibody buffer** for [🕒] 03:00:00 at [🌡️] Room temperature .

3h



### Note

Please note that this antibody incubation is performed differently from the first one (3 hours at room temperature).

- 9 When primary antibody incubation is finished, aspirate the media and wash the sections

5m



00:05:00 5x

with

[M] 0.05 % (v/v)



TWEEN 20 Sigma  
Aldrich Catalog #P7949

in PBS.

- 10 Incubate **Donkey α Rat 647 (Jackson, 712-605-153, 1:300) secondary antibody** in

1h



[M] 0.1 % (v/v)



TWEEN 20 Sigma  
Aldrich Catalog #P7949

for



[🕒] 01:00:00

at



[🌡️] Room temperature

. Add DAPI (**Invitrogen, D3571, 1:5000**)

## Wash and mounting

12h

- 11 When secondary antibody incubation is finished, wash the sections [🕒] 00:05:00 3x with

10m



[M] 0.05 % (v/v)



TWEEN 20 Sigma  
Aldrich Catalog #P7949

in PBS. Follow this with



[🕒] 00:05:00 2x

washes with

**PBS** to remove all detergent traces.

- 12 Clean the remaining buffer on the slides using absorbent tissue and mount the sections with a drop of



Fluoromount-G Electron Microscopy  
Sciences Catalog #17984-25