



Dec 15, 2020

CuPCR SARS-Cov-2

B Ronacher¹, Christoph Reschreiter²¹cube dx GmbH; ²Cube Dx GmbH

1

Works for me

dx.doi.org/10.17504/protocols.io.bpqwmmxe

Viral Testing Inc.

Christoph Reschreiter
Cube Dx GmbH

ABSTRACT

RNA based diagnostic test of different human samples (nasopharyngeal swabs, saliva) to detect SARS-Cov-2. The test requires RNA/DNA extraction and a follow-up RT-PCR (reverse transcriptase polymerase chain reaction). The test covers the N-gene and RdRp-gene and covers a human control as well. The human control safeguards the validity of sample taking and test procedure.

EXTERNAL LINK

<http://www.cubedx.com>

DOI

dx.doi.org/10.17504/protocols.io.bpqwmmxe

EXTERNAL LINK

<http://www.cubedx.com>

PROTOCOL CITATION

B Ronacher, Christoph Reschreiter 2020. CuPCR SARS-Cov-2 . **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bpqwmmxe>

EXTERNAL LINK

<http://www.cubedx.com>

KEYWORDS

RT-PCR SARS-Cov-2, N-gene, RdRp-gene, covid19, corona test

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 15, 2020

LAST MODIFIED

Dec 15, 2020

PROTOCOL INTEGER ID

44534

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

RNA based diagnostic test of different human samples (nasopharyngeal swabs, saliva) to detect SARS-Cov-2. The test requires RNA/DNA extraction and a follow-up RT-PCR (reverse transcriptase polymerase chain reaction). The test covers the N-gene and RdRp-gene and covers a human control as well. The human control safeguards the validity of sample taking and test procedure.

RT-PCR test

30m

30m

1 ⌚ 00:30:00 RNA / DNA extraction

Kingfisher Flex 96
RNA / DNA extraktor
Thermo Fisher Scientific 5400630
RNA/DNA extraktor processing up to 96
samples in a single run.



100µl eluation buffer (e.g. water) containing RNA / DNA from the sample

2 ⌚ 01:20:00 RT-PCR

1h 20m

Prepare mastermix:
5µl of primer mix + 5µl of enzyme mix (per sample)

Put mastermix and sample into PCR tube:
10µl of mastermix + 10µl of sample (eluate)

Select (or program) RT-PCR protocol:

- ☐ 2 minutes 50 °C
- ☐ 2 minutes 95 °C
- ☐ 45 cycles:
 - 5 seconds 95 °C
 - 1 minute 60 °C
 - Scan
- ☐ 25°C for hold

Load samples and start PCR

Analyse amplification curves:

Human control: ct < 34 --> valid
N-gene SARS-Cov-2: ct > 0 --> positive
RdRp-gene SARS-Cov-2: ct > 0 --> positive

qTower3

Realtime PCR

Analytic Jena 844-00563-2



Amplification curves for: human control, N-gene of SARS-Cov-2, RdRp-gene of SARS-Cov-2