



Version 2

Mar 24, 2021

FCMPASS - Cataloguing light scatter reference materials

V.2

Joshua A Welsh¹, Sean M Cook¹, Jennifer Jones¹¹Translational Nanobiology Section, Laboratory of Pathology, Center for Cancer Research, National Cancer Institute, National Institutes of Health1 Works for me dx.doi.org/10.17504/protocols.io.bniamcae

Translational Nanobiology Section



Sean Cook

SUBMIT TO PLOS ONE

ABSTRACT

This protocol outlines the steps required to catalogue light scatter reference materials using the FCMPASS software. This is one of a number of protocols in the pipeline for performing small particle calibration using the fcmpass software package.

DOI

dx.doi.org/10.17504/protocols.io.bniamcae

PROTOCOL CITATION

Joshua A Welsh, Sean M Cook, Jennifer Jones 2021. FCMPASS - Cataloguing light scatter reference materials. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bniamcae>

Version created by [Joshua Welsh](#)

WHAT'S NEW

Upd

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 16, 2020

LAST MODIFIED

Mar 24, 2021

PROTOCOL INTEGER ID

43298

MATERIALS TEXT

FCMPASS software can be accessed at <https://nanopass.ccr.cancer.gov>.

DISCLAIMER:

This protocol summarizes key steps for a specific type of assay, which is one of a collection of assays used for EV analysis in the NCI Translational Nanobiology Section at the time of submission of this protocol. Appropriate use of

this protocol requires careful, cohesive integration with other methods for EV production, isolation, and characterization. By using the FCMPASS software you agree to the following terms and conditions.

Terms & Conditions of use for FCMPASS software.

Definitions: The term "SOFTWARE" throughout this agreement means the machine readable, binary, object code form, and the related documentation for FCMPASS, a software package that is designed to allow flow cytometer calibration for small particles. The term "RECIPIENT" means the party that downloads the software. The term "PROVIDER" means the National Cancer Institute (NCI), a participating institute of the National Institutes of Health (NIH), and an agency of the United States Government. By downloading or otherwise receiving the SOFTWARE, RECIPIENT may use the SOFTWARE subject to RECIPIENT's agreement to the following terms:

1. THE SOFTWARE SHALL NOT BE USED IN THE TREATMENT OR DIAGNOSIS OF HUMAN SUBJECTS. RECIPIENT is responsible for compliance with all laws and regulations applicable to the use of the SOFTWARE.
2. RECIPIENT shall not distribute the SOFTWARE, in whole or in part without express advance written approval of PROVIDER.
3. The SOFTWARE may be used for research, academic, and educational purposes only. The SOFTWARE may not be used for commercial purposes. RECIPIENT will not license or sell or use the SOFTWARE for commercial purposes or applications.
4. The SOFTWARE that is distributed pursuant to this Agreement has been created by United States Government employees. In accordance with Title 17 of the United States Code, section 105, the SOFTWARE is not subject to copyright protection in the United States. Other than copyright, all rights, title and interest in the SOFTWARE shall remain with the PROVIDER.
5. RECIPIENT shall not modify, extend, decompile, make derivatives of, merge, publish, reverse engineer or distribute the SOFTWARE without written permission from PROVIDER.
6. RECIPIENT may publish or otherwise publicly disclose the results of using the SOFTWARE. RECIPIENT agrees to acknowledge PROVIDER's contribution of the SOFTWARE in all written publications containing any data or information regarding or resulting from use of the SOFTWARE.
7. THE SOFTWARE IS PROVIDED "AS IS" AND ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE AND NONINFRINGEMENT ARE DISCLAIMED. IN NO EVENT SHALL THE PROVIDER OR THE INDIVIDUAL DEVELOPERS BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING, BUT NOT LIMITED TO, PROCUREMENT OF SUBSTITUTE GOODS OR SERVICES; LOSS OF USE, DATA, OR PROFITS; OR BUSINESS INTERRUPTION) HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OF THIS SOFTWARE, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGE. PROVIDER makes no representations that the use of SOFTWARE will not infringe any patent or proprietary rights of third parties.
8. RECIPIENT may, on an as-needed basis, send to PROVIDER reports regarding the application of the SOFTWARE and the effectiveness and problems encountered in using the SOFTWARE, without disclosing RECIPIENT's confidential information. Information from general reports may be used by the PROVIDER to enhance the capabilities of the SOFTWARE. Reports can be forwarded to the PROVIDER at one of the following addresses: joshua.welsh@nih.gov or jennifer.jones2@nih.gov
9. No indemnification for any loss, claim, damage, or liability is intended or provided by either Party under this Agreement. Each Party shall be liable for any loss, claim, damage, or liability that said Party incurs as a result of said Party's activities under this Agreement, except that Provider, as an agency of the United States, assumes liability only to the extent as provided under the Federal Tort Claims Act (28 U.S.C. Chapter 171 Sections 2671-2680).
10. RECIPIENT agrees not to claim, infer, or imply endorsement by the United States Government, or any of its organizational units, contractors or employees. RECIPIENT agrees not to use any trademarks, service marks, trade names, logos or product names of NCI or NIH to endorse or promote products derived from the SOFTWARE without specific, prior and written permission.
11. Title in the SOFTWARE shall remain with the PROVIDER. It is understood that nothing herein will be deemed to constitute, by implication or otherwise, the grant to either party by the other of any license or other rights under any patent, patent application or other intellectual property right or interest. PROVIDER reserves the right to distribute the SOFTWARE to others and to use it for PROVIDER's own purposes. The United States Government explicitly retains all rights to use the SOFTWARE for any purpose, to have it used on the Government's behalf or to allow others to use it.

- 1 Open FCM_PASS.

FCMPASS 3 [↗](#)

- 2 Click 'Catalogue' in the top menu bar

Bead Catalogue Basic Protocol

- 3 Under the 'Light Scatter' tab entry fields exist for each of the pertinent metadata for reporting with light scatter calibration.

3.1 Diameter CV should be the percent coefficient of variation of the mean diameter provided on the certificate of analysis

3.2 Refractive Index should be the provided refractive index of the bead population on certificate of analysis

If a refractive index is not available an approximate guide for polystyrene refractive index is 1.59 at 589 nm. Silica tends to vary more in refractive index than polystyrene but tends to be ~1.45 at 589 nm.

3.3 'RI Measurement Wavelength' is the wavelength at which the refractive index was measured and should be provided on the certificate of analysis. This tends to be 589 nm.

3.4 Composition can be selected as polystyrene, silica, or other. If polystyrene or silica are selected, changes in detection wavelength e.g. 488 nm to 405 nm are accounted for using the appropriate Sellmeier equations. If 'Other' is selected then the refractive index change is made propositionally to the sheath refractive index.

3.5 Manufacturer, Catalogue Number, and Lot Number should all be completed appropriately.

- 4 Once the entry fields are filled out as desired, press the 'Add Bead' Button to add the bead information to the table. Each button has an in-depth explanation at the end of the protocol.

- 5 Continue to add beads until all the desired beads have been added to the table. Next, select the desired beads that should be in the bead set. This is done from the first column of the table which is labeled 'Selection' and is marked by a red box below.

	Selection	Diameter (nm)	Diameter CV (%)	Refractive Index (RI)	RI Wavelength (nm)	Catalogue Number	Manufacturer	Lot Number	Set
1	<input type="checkbox"/>	81	11.7000	1.5900	589	3080A	ThermoFisher Scientific	228748	Example ...
2	<input type="checkbox"/>	100	7.8000	1.5900	589	3100A	ThermoFisher Scientific	204935	Example ...
3	<input type="checkbox"/>	152	3.3000	1.5900	589	3150A	ThermoFisher Scientific	202026	Example ...
4	<input type="checkbox"/>	203	2.6000	1.5900	589	3200A	ThermoFisher Scientific	205131	Example ... jhv
5	<input type="checkbox"/>	269	1.6000	1.5900	589	3269	ThermoFisher Scientific	202729	Example ... jhv
6	<input type="checkbox"/>	303	1.6000	1.5900	589	3300A	ThermoFisher Scientific	204665	Example ... jhv
7	<input type="checkbox"/>	345	1.9000	1.5900	589	3350A	ThermoFisher Scientific	199283	Example ...
8	<input type="checkbox"/>	401	1.3000	1.5900	589	3400A	ThermoFisher Scientific	203859	Example ...
9	<input type="checkbox"/>	453	1.7000	1.5900	589	3450A	ThermoFisher Scientific	204047	Example ...
10	<input type="checkbox"/>	480	4.2000	1.4500	589	8050	ThermoFisher Scientific	203277	Example ...
11	<input type="checkbox"/>	508	1.7000	1.5900	589	3500A	ThermoFisher Scientific	204667	Example ...
12	<input type="checkbox"/>	600	1.7000	1.5900	589	3600A	ThermoFisher Scientific	205833	Example ...
13	<input type="checkbox"/>	730	4.1000	1.4500	589	8070	ThermoFisher Scientific	207434	Example ...

The 'Selection' column allows for the selection of beads, this functionality is used by the Bead Functions and Set Functions which will be discussed in further detail. When a bead is selected, its information is written in the entry fields above.

- Press the 'Create Set' button when all the desired beads have been selected. A pop-up dialog box will ask for a unique name. Upon completion, the new bead set will now be accessible from the Display Set Dropdown.

If you have created your desired bead set, you can now continue using the FCMPASS software.

However, the following steps in this protocol outline the complete functionality of the Bead Catalogue.

Bead Functions

- The first bead function button is the '**Add Bead**' button marked by the blue box. When pressed, the 'Add Bead' button adds the entry fields information into a row in the table below.

Note: The new bead is added to the set listed in the Display Set Dropdown that is marked by the red box. Thus, beads can be added to any existing bead sets by changing the Display Set Dropdown. If a bead is added to an existing bead set, then it is also added to the master 'All' bead set in the background.

Reference Bead Catalogue

View Tutorial

Light Scatter | Fluorescence

0 Diameter (nm) Manufacturer Polystyrene Composition

0 Diameter CV (%) Lot Number Catalogue Number

0 Refractive Index (RI) 0 RI Measurement Wavelength (nm)

Bead Functions

Add Bead (highlighted)

Delete Bead

Append Bead

Set Functions

Create Set

Delete Set

Append Set

Example Set ▼

Display Set

All ▼

	Selection	Diameter (nm)	Diameter CV (%)	Refractive Index (RI)	RI Wavelength (nm)	Catalogue Number	Manufacturer	Lot Number	Set
1	<input type="checkbox"/>	81	11.7000	1.5900	589	3080A	ThermoFisher Scientific	228748	Example ...
2	<input type="checkbox"/>	100	7.8000	1.5900	589	3100A	ThermoFisher Scientific	204935	Example ...
3	<input type="checkbox"/>	152	3.3000	1.5900	589	3150A	ThermoFisher Scientific	202026	Example ...
4	<input type="checkbox"/>	203	2.6000	1.5900	589	3200A	ThermoFisher Scientific	205131	Example ... jhw
5	<input type="checkbox"/>	269	1.6000	1.5900	589	3269	ThermoFisher Scientific	202729	Example ... jhw
6	<input type="checkbox"/>	303	1.6000	1.5900	589	3300A	ThermoFisher Scientific	204665	Example ... jhw
7	<input type="checkbox"/>	345	1.9000	1.5900	589	3350A	ThermoFisher Scientific	199283	Example ...
8	<input type="checkbox"/>	401	1.3000	1.5900	589	3400A	ThermoFisher Scientific	203859	Example ...
9	<input type="checkbox"/>	453	1.7000	1.5900	589	3450A	ThermoFisher Scientific	204047	Example ...
10	<input type="checkbox"/>	480	4.2000	1.4500	589	8050	ThermoFisher Scientific	203277	Example ...
11	<input type="checkbox"/>	508	1.7000	1.5900	589	3500A	ThermoFisher Scientific	204667	Example ...
12	<input type="checkbox"/>	600	1.7000	1.5900	589	3600A	ThermoFisher Scientific	205833	Example ...
13	<input type="checkbox"/>	730	4.1000	1.4500	589	8070	ThermoFisher Scientific	207434	Example ...

- 8 The second bead function button is the '**Delete Bead**' button which is immediately below the 'Add Bead' button. To delete beads, they must be selected from the first column labeled selection. You may select as many beads as desired to delete. The beads will only be deleted from the currently selected bead set from the Display Set Dropdown.


If beads are deleted from the master 'All' bead set, the beads are deleted from all bead sets. This is because the 'All' bead set serves as a master list of all stored bead information.

- 9 The last bead function button is the '**Append Bead**' button which is immediately below the 'Delete Bead' button. When appending a bead, a single bead must be selected from the selection column. Upon pressing the 'Append Bead' button, the bead information is overwritten by the information in the entry fields above. You can append any bead in any bead set, but you cannot append more than one bead at a time.

As discussed earlier in the note for step 5, by selecting a bead, the bead information is written into the entry fields above. This allows for easy appending of beads that may only need a single entry field changed.

Set Functions

- 10 The first set function button is the '**Create Set**' button marked by the red box below. First select the desired beads that should be included in the new bead set. Then, press the 'Create Set' button and give the bead set a unique name. The new bead set will now be accessible from the Display Set Dropdown.




A bead set can be created from any other bead set.

- 11 The second set function button is the **'Delete Set'** button which is immediately below the 'Create Set' button. The 'Delete Set' button deletes the currently displayed bead set. This is the bead set selected from the Display Set Dropdown. The master 'All' bead set cannot be deleted.

The beads will not be deleted, only the bead set. The beads will remain in their other bead sets including the master 'All' bead set.

- 12 The last set function is the **'Append Set'** button which is immediately below the 'Delete Set' button. The 'Append Set' button allows for the addition of beads to existing bead sets. **The 'Append Set' button uses the Dropdown directly below, and the selected beads are appended to the selected bead set from that Dropdown (marked by a blue box below).**



A bead cannot be appended to a bead set that it is already in.