



May 28, 2020

# Standard Volvox Medium

Raymond Goldstein<sup>1</sup>, Kyriacos Leptos<sup>1</sup><sup>1</sup>University of Cambridge

1 Works for me

This protocol may be deleted by the owner

Kyriacos Leptos

## ABSTRACT

This is a protocol for the execution of a Standard Volvox Medium recipe. Starr's (1969) modification of Provasoli & Pintner's Medium (1959), as further modified by Starr, Kochert and Kirk at various times.

## GUIDELINES

### Description

Class of Chemicals	Stock Solution	Stock Concentration	Final Concentration
<b>Macro minerals, N-source, and S-source</b>	Ca(NO <sub>3</sub> ) <sub>2</sub>	500 mM	0.5 mM
	MgSO <sub>4</sub>	160 mM	0.16 mM
	KCl	670 mM	0.67 mM
	Na <sub>2</sub> CO <sub>3</sub>	190 mM	0.19 mM
<b>N-source</b>	Urea	500 mM	0.5 mM
<b>P-source</b>	Na <sub>2</sub> Glycerophosphate	163 mM	0.16 mM
<b>pH Regulator</b>	HEPES	250 mM	2.5 mM
<b>Trace minerals</b>	Na <sub>2</sub> EDTA	4.47 mM	13.4 μM
	FeCl <sub>3</sub>	717 μM	2.15 μM
	MnCl <sub>2</sub>	414 μM	1.24 μM
	ZnCl <sub>2</sub>	73 μM	0.22 μM
	CoCl <sub>2</sub>	17 μM	0.05 μM
	Na <sub>2</sub> MoO <sub>4</sub>	38.8 μM	0.10 μM
<b>Vitamins</b>	Thiamine (Vitamin B1)	1 mg/mL	1 μg/mL
	Biotin (Vitamin B7)	2.5 μg/mL	0.25 ng/mL
	Cyanocobalamin (Vitamin B12)	1.5 μg/mL	0.15 ng/mL

Description of the SVM Medium compounds and their final concentration.

### Molar Composition

- [M] **2.2 mM** Na<sup>+</sup>
- [M] **670 μM** K<sup>+</sup>
- [M] **160 μM** Mg<sup>2+</sup>
- [M] **500 μM** Ca<sup>2+</sup>
- [M] **679 μM** Cl<sup>-</sup>
- [M] **160 μM** SO<sub>4</sub><sup>2-</sup>
- [M] **2.5 mM** [HEPES]<sup>-</sup>
- [M] **1 mM** NO<sub>3</sub><sup>-</sup>
- [M] **0.5 mM** Urea
- [M] **0.16 mM** [Glycerophosphate]<sup>2-</sup>
- [M] **13 μM** EDTA<sup>4-</sup>
- [M] **2.2 μM** Fe<sup>2+</sup>
- [M] **0.22 μM** Zn<sup>2+</sup>

- [M] **1.2  $\mu\text{M}$**   $\text{Mn}^{2+}$
- [M] **0.1  $\mu\text{M}$**   $\text{MoO}_4^{2-}$
- [M] **0.05  $\mu\text{M}$**   $\text{Co}^{2+}$

#### MATERIALS TEXT

##### Required Materials:

- Stock solutions (see below)
- NaOH
- 1 x 2 L Beaker
- 1 x 1 L Measuring cylinder
- 2 x 1 L Bottles

##### Stock solutions:

###### *Macro minerals, N-source, and S-source*

- [M] **500 mM**  $\text{Ca}(\text{NO}_3)_2$
- [M] **160 mM**  $\text{MgSO}_4$
- [M] **670 mM** KCl
- [M] **190 mM**  $\text{Na}_2\text{CO}_3$

###### *N-source*

- [M] **500 mM** Urea

###### *P-source*

- [M] **163 mM**  $\text{Na}_2\text{Glycerophosphate}$

###### *pH Regulator*

- [M] **250 mM** HEPES

###### *Trace Minerals ("P IV Metal solution" 2X)*

- [M] **4.47 mM**  $\text{Na}_2\text{EDTA}$
- [M] **717  $\mu\text{M}$**   $\text{FeCl}_3$
- [M] **414  $\mu\text{M}$**   $\text{MnCl}_2$
- [M] **73  $\mu\text{M}$**   $\text{ZnCl}_2$
- [M] **38.8  $\mu\text{M}$**   $\text{Na}_2\text{MoO}_4$

###### *Vitamins*

- [M] **1 mg/mL** Thiamine (Vitamin B<sub>1</sub>)
- [M] **2.5  $\mu\text{g/mL}$**  Biotin (Vitamin B<sub>7</sub>)
- [M] **1.5  $\mu\text{g/mL}$**  Cyanocobalamin (Vitamin B<sub>12</sub>)



#### Storage and Shelf Live of Stock Solutions

⌄ Room temperature

- $\text{Na}_2\text{CO}_3$
- P IV Metal Solution

#### 4 °C

- $\text{Ca}(\text{NO}_3)_2$
- $\text{MgSO}_4$
- KCl urea
- Phosphoglycerate HEPES
- Thiamine (wrapped in foil)

#### -20 °C

- Biotin
- Cyanocobalamin (in dark e.g. a cryobox)

#### Additional Notes:

1. All stocks except **biotin** and **cyanocobalamin** should be made fresh at least every 2 months.
2. The **trace metal** stocks should be made fresh every month, and sooner if there is any problem with the cultures that cannot be explained.

#### SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

#### BEFORE STARTING

The amounts in this protocol are shown for **1 L SVM** but the protocol can be scaled up or down as desired if more or less SVM is needed.

#### Adding Solvent

- 1 Add **950 ml Initial MilliQ water** to a 2L beaker.

#### Adding Macro minerals, N-source, and S-source

- 2 Add **1 ml  $\text{Ca}(\text{NO}_3)_2$**  .
- 3 Add **1 ml  $\text{MgSO}_4$**  .
- 4 Add **1 ml KCl** .
- 5 Add **1 ml  $\text{Na}_2\text{CO}_3$**  .

#### Adding N-source

6 Add  **1 ml Urea** .

Adding P-source

7 Add  **1 ml Na<sub>2</sub>Glycerophosphate** .

Adding pH Regulator

8 Add  **10 ml HEPES** .

Adding Trace Minerals


9 Add  **3 ml "P IV Metal solution" 2X** .

Adjusting pH and final volume

10 Adjust the pH to  **7.8** with NaOH.



Make sure this is done under mild mixing.

11 Adjust the final volume to  **1000 ml** .

Autoclaving

12 Dispense the mixture as 500 ml batches into desired 1L bottles.


13 Autoclave the medium.

14 Keep the medium at  **4 °C** until you are ready to use it.

Adding Vitamins (Just before use)

15






Always add vitamins fresh and keep the medium at  **4 °C** .



Be sure to work under a laminar flow hood of this step.

Using pre-plugged **sterile** tips mix the following as a cocktail in a 1.5 ml **sterile** microcentrifuge:

-  **1000 µl Thiamine (Vitamin B1)**
-  **100 µl Biotin (Vitamin B7)**
-  **100 µl Cyanocobalamin**

 **15.1** Add  **1000 µl Thiamine (Vitamin B1)** .

 **15.2** Add  **100 µl Biotin (Vitamin B7)** .

 **15.3** Add  **100 µl Cyanocobalamin** .

 **15.4** Mix.

 **16** Add the vitamins to the SVM-vit.