



## Oct 18, 2021

## Preparation of BG-11 liquid media

## Ashwinuday 1

<sup>1</sup>Indian Institute of Science Education and Research, Pune



dx.doi.org/10.17504/protocols.io.by6qpzdw

iGEM IISER Pune India 2021

Ashwinuday

This is for making liquid media of BG-11.

DOI

dx.doi.org/10.17504/protocols.io.by6qpzdw

Ashwinuday 2021. Preparation of BG-11 liquid media . **protocols.io** https://dx.doi.org/10.17504/protocols.io.by6qpzdw

•

\_\_\_\_\_ protocol,

Oct 18, 2021

Oct 18, 2021

54192

- 1 Take the medium in powdered form and mix it with doubly distilled water to just 1/10 th of the total volume.
- 2 Add a common stock of micronutrients made in ddH2O and mix well.
- 3 Make up the volume to the total and mixed well.



4	Adjust the pH on a calibrated pH meter using 1N HCl and /or 1N NaOH. The media was mixed using a magnetic stirrer continuously to homogenize the pH and to prevent precipitation of certain constituents.
5	This mix was autoclaved.