



Sep 28, 2021

# Western blot protocol for detecting ATP10B in mouse/rat brain

María Sanchiz Calvo<sup>1</sup>, eduard.bentea<sup>1</sup>, Veerle Baekelandt<sup>1</sup><sup>1</sup>Katholieke Universiteit Leuven

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.byhfpt3n](https://dx.doi.org/10.17504/protocols.io.byhfpt3n)

eduard.bentea

## ABSTRACT

Protocol for detection of ATP10B in rat and mouse brain tissue by Western blotting

## DOI

[dx.doi.org/10.17504/protocols.io.byhfpt3n](https://dx.doi.org/10.17504/protocols.io.byhfpt3n)

## PROTOCOL CITATION

María Sanchiz Calvo, eduard.bentea , Veerle Baekelandt 2021. Western blot protocol for detecting ATP10B in mouse/rat brain. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.byhfpt3n>

## KEYWORDS

ATP10B, western blot, mouse, rat, brain

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 24, 2021

## LAST MODIFIED

Sep 28, 2021

## PROTOCOL INTEGER ID

53511

## MATERIALS TEXT

Recipe for 100 mL 1x RIPA buffer :

- o 5 mL Tris pH 7.4
- o 1 mL Triton x100
- o 1 g Deoxycholate
- o 3 mL 5M NaCl
- o 200 µL 0.5M EDTA (pH 8.0)
- o 1 mL 10% SDS

Add water until 100 mL (MilliQ water)

Add protease and phosphatase inhibitors (Pi + PPi) right before use

Reagents used during Western blotting :

- o 6X SDS buffer
- o PageRuler™ Plus Prestained Protein Ladder
- o 3-8% NuPAGE™ Tris-Acetate gel
- o NuPAGE™ Tris-Acetate SDS Running Buffer
- o Trans-Blot® Turbo™ PVDF Membrane
- o 5% milk powder dissolved in PBS-T 0.1% (PBS + 0.1% Triton-X100)
- o ATP10B HPA034574 (Sigma) primary antibody
- o Goat anti-rabbit HRP secondary antibody (Dako)
- o ECL Prime chemiluminescence kit (GE Healthcare)

Recipe for 50 mL 6x SDS buffer (160 mM TrisHCl pH 6.8, 2% SDS, 200 mM DTT, 40% glycerol, bromophenol blue):

- o 5 mL 1.6M TrisHCl pH 6.8
  - o 1 g SDS
  - o 1.54 g DTT
- Add water until 25 mL (MilliQ water) and dissolve
- o 20 mL glycerol
- Add water until 50 mL (MilliQ water)
- Add spatula tip of bromophenol blue









## BEFORE STARTING

Perform protein extraction from snap-frozen brain tissue:
















- weigh tissue
- add RIPA buffer (see Materials) : 10 X of the weight (40 mg = 400 µL)
- homogenize samples using sample homogenizer
- sonicate samples at 4 degrees C, 3 times 15 seconds (keep the samples on ice between each sonication)
- centrifuge samples at 6000 g for 10 minutes at 4 degrees C
- collect supernatant and measure protein concentration
- aliquot protein extracts and store at -20 degrees C

Day 1 4h 5m

- 1 Prepare samples for Western blotting : 30m
  - 30 µg proteins in 12 µL of volume (adjust with milliQ water)
  - 2.4 µL of 6X SDS buffer + 10% of β-mercaptoethanol
- 2 Vortex samples and spin down
- 3 Boil the samples for 00:10:00 at 98 °C 10m

- 4 Load samples, together with  **7 µl** of mass marker (PageRuler™ Plus Prestained Protein Ladder) on a 3-8% NuPAGE™ Tris-Acetate gel. Use NuPAGE™ Tris-Acetate SDS Running Buffer for migration. 30m
- 5 Run the gel at 80V for  **00:10:00** 10m
- 6 Run the gel at 150V for another  **00:45:00** 45m
- 7 Transfer the proteins on PVDF membrane (Trans-Blot® Turbo™ PVDF Membrane) using Trans-Blot Turbo Transfer System (Bio-Rad), using the pre-programmed protocol STANDARD SD:  **00:30:00** , up to 1.0 A, 25 V. 30m
- 8 Block membranes for  **01:00:00** using 5% milk dissolved in PBS-T 0,1% at  **Room temperature** 1h
- 9 Incubate membranes  **Overnight** in ATP10B HPA034574 (Sigma) diluted 1/500 in 5% milk PBS-T 0,1% at  **4 °C** 30m

Day 2 3h 10m

- 10 Wash membranes with PBS-T 0,1% for 10 minutes, 5 times at  **Room temperature** : 50m
  - PBS-T  **00:10:00**
  - PBS-T  **00:10:00**
  - PBS-T  **00:10:00**
  - PBS-T  **00:10:00**
  - PBS-T  **00:10:00**
- 11 Incubate membranes for  **02:00:00** in secondary antibody solution Goat anti-rabbit HRP (Dako) diluted 1/10000 in 5% milk PBS-T 0,1% at  **Room temperature** 2h
- 12 Wash membranes with PBS-T 0,1% for 10 minutes, 5 times at  **Room temperature** : 50m
  - PBS-T  **00:10:00**
  - PBS-T  **00:10:00**
  - PBS-T  **00:10:00**
  - PBS-T  **00:10:00**
  - PBS-T  **00:10:00**
- 13 Develop membranes using ECL Prime (GE Healthcare) for at least  **00:10:00** 10m

