


Jun 17, 2021

Expression and purification protocol of FIP200 full length or FIP200 Δ 641-779

Xiaoshan Shi¹¹University of California, Berkeley

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.bvjjn4kn

ULK1 complex purification



Xiaoshan Shi

University of California, Berkeley, Shanghai Institute of Bi...

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

This protocol describes the expression and purification of FIP200 full length or FIP200 Δ 641-779.

ATTACHMENTS

Expression and purification
protocol of FIP200 full
length or FIP200 Δ 641-
779.pdf

DOI

dx.doi.org/10.17504/protocols.io.bvjjn4kn

PROTOCOL CITATION

Xiaoshan Shi 2021. Expression and purification protocol of FIP200 full length or FIP200 Δ 641-779.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.bvjjn4kn>

KEYWORDS

FIP200 full length or FIP200 Δ 641-779, Expression and purification protocol of FIP200 full length or FIP200 Δ 641-779, FIP200 full length, expression, purification, FIP200

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 06, 2021

LAST MODIFIED

Jun 17, 2021

OWNERSHIP HISTORY

Jun 06, 2021  zahara93.zs

Jun 16, 2021  Xiaoshan Shi University of California, Berkeley

PROTOCOL INTEGER ID

50507

MATERIALS TEXT

Lysis Buffer: 50 mM Hepes pH=7.4; 1% Triton X-100, 200 mM NaCl, 1 mM MgCl₂, 10% Glycerol, 1 mM TCEP, EDTA-free Protease Inhibitors (Roche).

Wash Buffer 1: 50 mM Hepes pH=7.4, 500 mM NaCl, 1 mM MgCl₂, 1 mM TCEP, 1% Triton X-100, 10% Glycerol;

Wash Buffer 2: 20 mM Hepes pH=8, 200 mM NaCl, 1 mM MgCl₂, 1 mM TCEP;

Elution Buffer 1: 50 mM Hepes pH=7.4, 500 mM NaCl, 1 mM MgCl₂, 1 mM TCEP, 1% Triton X-100, 10% Glycerol, 50 mM Glutathione;

Elution Buffer 2: 20 mM Hepes pH=8, 200 mM NaCl, 1 mM MgCl₂, 1 mM TCEP, 50 mM Maltose;

Resin:

- Glutathione Sepharose 4B (GE Healthcare)
- Amylose resin (New England Biolabs).

SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

BEFORE STARTING

General information: expression system: human embryonic kidney (HEK) 293 GnTI suspension cells, medium: freestyle media + Anti-Anti; plasmid origin: Addgene 171410 or 17141; backbone: pCAG; resistance: Amp; insert: *Homo sapiens* FIP200 (synthetic gene); tags & cleavage sites: GST-TEVcs-FIP200-MBP, Ext coeff: 183470 M⁻¹cm⁻¹, MW 252 kDa

1- Protein expression: 2d 11h 10m

- 1 Transfect DNA in cells at a concentration of 2.5–3×10⁶/mL using polyethylenimine (Polysciences) and harvest after ^{2d}
⌚ 48:00:00 expression, harvest and lyse cells with lysis buffer.

- 2  10m

Pellet the harvested cells at ⚙ 500 x g for ⌚ 00:10:00 at 🌡 4 °C, wash with PBS once, and then store at

8 -80 °C .

- 3 Resuspend the cell pellet in Lysis Buffer, gently shake at 4 °C for 00:30:00 and clear it at 16.000 rpm at 4 °C for 00:30:00 . 1h

2- Protein Purification:

10h


4 

10h

Glutathione Sepharose 4B (GE Healthcare) followed by MBP-tag affinity purifications.

Incubate supernatant (added 6M NaCl to a final concentration of 0.5 M) with Glutathione Sepharose 4B resin at 4 °C with gentle shaking for 10:00:00 , apply to a gravity flow column, and wash extensively with Wash Buffer 1.

- 5 Elute protein of interest with Elution Buffer 1 and then further apply onto Amylose resin for a second step of affinity purification. Wash Amylose resin with Wash Buffer 2, and elute Protein of interest with Elution Buffer 2.
- 6 Pool peak fractions containing pure protein, concentrate, snap-freeze in liquid nitrogen and store at -80 °C .



Yield: about 0.2 mg per litter culture