



Feb 15, 2021

Methanol-based HIV-Flow

Marion Pardons¹, Basiel Cole¹¹HIV Cure Research Center, Department of Internal Medicine and Pediatrics, Ghent University Hospital, Ghent University, 9000 Ghent, Belgium

1

Works for me

This protocol is published without a DOI.

Marion Pardons

SUBMIT TO PLOS ONE

PROTOCOL CITATION

Marion Pardons, Basiel Cole 2021. Methanol-based HIV-Flow. **protocols.io**
<https://protocols.io/view/methanol-based-hiv-flow-bpedmja6>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 04, 2020

LAST MODIFIED

Feb 15, 2021

PROTOCOL INTEGER ID

44197

MATERIALS TEXT

- PMA (Sigma, P8139)
- Ionomycin (Sigma, I9657)
- ARV: 3TC/Raltegravir
- Fixable Viability Stain 510 (ThermoFisher Scientific, L34957)
- p24 KC57-FITC (Beckman Coulter, 6604665)
- p24 28B7-APC (MediMabs, MM-0289-APC)

Day 1 - CD4 enrichment and stimulation

- 1 Thaw cryopreserved PBMCs (maximum 2 vials per 50mL tube) from the HIV-infected individuals and uninfected control in FCS (1mL FCS per 50M cells).
Centrifuge (1,500rpm, 5min), and discard supernatant.
Resuspend the cells in 25mL of cRPMI and count the cells.
Centrifuge (1500rpm, 5min), and discard supernatant.
- 2 Proceed to the negative selection of CD4 T cells
-According to cell count:
 - a. If more than 50.10⁶cells: Resuspend the pellet with 800µL of MACS Buffer and complete up to a concentration of 50x10⁶ cells/mL.
 - b. If 10x10⁶-50x10⁶cells: Resuspend the pellet with 800µL of MACS Buffer-Transfer the cell suspension in a 5mL FACS tube (max of 150M of cells per tube = 3mL).
-Add 25µL/mL of cells of Ab cocktail.
-Incubate 5 min at RT.
-Vortex beads 30s and add 25µL of rapidospheres per mL of cells.
-Add MACS buffer up to 3mL.
-Place the tube (without the cap) in the BD magnet for 5min.

-Pipet cautiously the cells and transfer them to a 15mL tube containing cRPMI up to 10mL.
-Centrifuge the cells (5min, 1,500rpm), and discard the supernatant

- 3 Resuspend the cells in 800µL of cRPMI (+ ARV) and count the cells
- 4 Resuspend cells at 2×10^6 cells/mL in cRPMI (+ ARV) and incubate at 37°C 5% CO₂.
- 5 Rest the cells for at least 1 hour before proceeding to PMA/ionomycin stimulation
-If ART-treated individuals: PMA 162nM, Ionomycin 1µg/mL, 24h
-If untreated individuals: PMA 25nM, Ionomycin 1µg/mL, 18h
- 6 Incubate cells for 18h or 24h (for samples from untreated VS treated individuals, respectively) at 37°C, 5% CO₂

Day 2 - Methanol-based HIV-Flow

- 7 Collect cells in 5mL tubes and wash the wells with cold PBS. Do not exceed 5×10^6 cells per tube.
- 8 Centrifuge at 1,800rpm for 5min.
- 9 Discard the supernatant and resuspend the cells in 200µL of PBS + Live/Dead stain (1µL L/D stain per mL of PBS). Incubate the cells for 25min at RT.
- 10 Wash the cells with 800µL of PBS/FCS2%. Centrifuge (1,800 rpm, 5min). Discard the supernatant.
- 11 Perform the extracellular staining in PBS/FCS2% (final volume = 100µL/tube) (25min, 4°C).
- 12 Wash with 1.8mL of PBS/FCS2%. Centrifuge (1,800rpm, 6min, **4°C**). Discard the supernatant. As soon as the supernatant is discarded, transfer the tubes on ice.
- 13 Vortex the pellet to avoid clumps (keep the tubes on ice!).
- 14 Gently add 1mL of cold methanol to the cells.
As soon as methanol is added to each tube, keep the tube on ice.
Incubate 15min, on ice.
- 15 Centrifuge at 2100rpm for 5min, 4°C.

- 16 Discard supernatant and wash the cells with 1mL of PBS/FCS2%.
- 17 Centrifuge at 2100rpm for 5min, 4°C.
- 18 Discard supernatant and proceed to intracellular staining in PBS/FCS2% (final volume = 100µL/tube) (45min, **RT**).
Antibodies P24 KC57 FITC 1/500
P24 28B7 APC (Stock: 25µg in 125µL H2O) 1/500

Dilution A: 1.5uL KC57 + 1.5uL p24 28B7 + 12uL PBS/FCS2%
Dilution B: 10uL Dilution A + 490uL PBS/FCS2%

Do not vortex the antibodies and the mixes (homogenize with 1mL pipette)
- 19 Wash with 1.8mL of PBS/FCS2%. Centrifuge (2,100 rpm, 5min). Discard the supernatant.
- 20 Resuspend cells in the adequate volume of PBS (NO FCS!).