



JAN 29, 2024

🌐 Labeling lysosomes in the adult *Drosophila* brain using LysoTracker

Mel Feany^{1,2}

¹Department of Pathology, Brigham and Women's Hospital, Harvard Medical School, Boston, Massachusetts 02115;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815

ASAP Collaborative Research Network

OPEN  ACCESS



Beatrice Weykopf

ABSTRACT

This protocol describes the labeling and quantification of Lysosomes in the adult *Drosophila* brain using the LysoTracker dye

DOI:

dx.doi.org/10.17504/protocols.io.j8nlkojjwv5r/v1

Protocol Citation: Mel Feany 2024. Labeling lysosomes in the adult *Drosophila* brain using LysoTracker. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkojjwv5r/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


Protocol status: Working
We use this protocol and it's working

Created: Jan 29, 2024

Last Modified: Jan 29, 2024

PROTOCOL integer ID: 94349

Keywords: ASAPCRN

- 1 Prepare the working solution, 0.05 μ M LysoTracker Green DND26, using 1X PBS
- 2 Brains from 10-day-old animals are dissected in ice-cold 1X PBS
- 3 The dissected brain is transferred to a glass slide, and the excess 1X PBS removed
- 4 The brains are re-incubated with 0.05 μ M LysoTracker for  00:05:00 at  Room temperature 
- 5 After 5 minutes, the excess solution is removed, a drop of DAPI Fluoromount added, and the brain mounted with a coverslip
- 6 Brains re imaged immediately using a Zeiss confocal microscope under a 63X objective lens
- 7 The number of puncta per 1000 μ m² is counted in Zen software

