



SEP 27, 2023

KAPP-Sen TMC: Dissociation of Pancreatic Islets (recovered)

Jessica

Juliana Alcoforado Diniz¹, Garofalo¹,

Dylan Baker¹, Paul Robson^{1,2,3}

¹The Jackson Laboratory for Genomic Medicine, Farmington, CT, USA;

²Department of Genetics and Genome Sciences, University of Connecticut School of Medicine, Farmington, CT, USA;

³Institute for Systems Genomics, University of Connecticut, Farmington, CT, USA

Cellular Senescence Network (SenNet) Method Development Community

KAPP-Sen TM

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.x54v9pyq4g3e/v1

Protocol Citation: Juliana Alcoforado Diniz, Jessica Garofalo, Dylan Baker, Paul Robson 2023. KAPP-Sen TMC: Dissociation of Pancreatic Islets (recovered). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v9pyq4g3e/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Aug 01, 2023

Last Modified: Sep 27, 2023

PROTOCOL integer ID:
 85794



Ashley M Raynock

UConn Health, UConn Center on Aging

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

The dispersed samples were shipped cold from **PRODOLABS**. Prior to scRNA-seq dispersed samples from brain dead donor's pancreatic islets were recovered and dissociated as follows.

Cell Dissociation with Accutase

1

Note

Before beginning cell dissociation coat all the materials (pipettes, tubes, etc.) with PIM-S001GMP media to prevent sticking.

- 1.1 Transfer cell suspension of pure islets to a new 50ml tube. Use additional media to rinse original container.
- 1.2 Centrifuge at room temperature 130g for 5 mins.
- 1.3 Aspirate the supernatant and add media to the appropriate concentration of 1,000 islets/1 ml.
--> Ex: If receiving 7,000 islets, use 7 ml media.
- 1.4 To recover cells, add half of the islets to a coated flask and incubate at 37°C overnight.
- 1.5 Coat pipette and add islets from flask to a 50ml tube. Wash flask with media to make sure all islets have been collected.
- 1.6 Centrifuge at 130 g, room temperature for 2 mins.
- 1.7 Aspirate media and resuspend in 4 ml accutase. Incubate at 37C for 8 mins, mixing with pipette every 2 mins.
Check at 6 mins. -->1 mL accutase/1,000 islets
islets

- 1.8 Add CMRL 1066 (Cat. 11530037) media to 25 ml then centrifuge at 230 g.
- 1.9 Aspirate supernatant and resuspend in 1.5 ml of CMRL
- 1.10 Filter through a 40 µm Flowmi.
- 1.11 Count cells using AO/PI (acridine orange/propidium iodide) Cell Viability Kit for Luna-FL automated cell counter.
- 1.12 Proceed to cell fixation.

Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling

- 2 Cells were fixated prior to scRNAseq according to <https://dx.doi.org/10.17504/protocols.io.x54v9py5zg3e/v1>