



OCT 16, 2023

Microscopy-based GSH bead protein-protein interaction assay

Minghao Chen¹, Xuefeng Ren¹

¹Department of Molecular and Cell Biology, University of California, Berkeley, Berkeley, CA 94720, USA.

ASAP Collaborative Research Network



Minghao Chen

OPEN  ACCESS



ABSTRACT

Microscopy-based GSH bead protein-protein interaction assay

DOI:

dx.doi.org/10.17504/protocols.io.4r3l27xdxg1y/v1

Protocol Citation: Minghao Chen, Xuefeng Ren 2023. Microscopy-based GSH bead protein-protein interaction assay. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l27xdxg1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jun 05, 2023

Last Modified: Oct 16, 2023



PROTOCOL integer ID:
82919

Keywords: ASAPCRN

Funders

Acknowledgement:

Aligning Science Across
Parkinson's (ASAP) initiative
Grant ID: ASAP-000350

- 1 Wash pre-blocked glutathione sepharose beads (GE Healthcare) with the reaction buffer (25 mM HEPES at pH 7.5, 150 mM NaCl, 1 mM MgCl₂ and 1 mM TCEP) three times.
- 2 Make a mixture of 1 μ M purified GST tagged protein and 500 nM purified fluorescent protein in total 70 μ L
- 3 Incubate at  Room temperature for  00:30:00, samples were mixed with additional 100 μ L reaction buffer. 30m
- 4 Transferred to the observation chamber for imaging.
- 5 Acquire images on a Nikon A1 confocal microscope with a Nikon Plan APO VC 20x/0.75 NA UV Microscope Objective.