

Oct 03, 2024

## Zymo OneStep PCR Inhibitor Removal Kit

DOI

**dx.doi.org/10.17504/protocols.io.5qpvokmrbl4o/v1**

Shannon Fitz<sup>1</sup>, Alex Shaw<sup>1</sup>, michael Owusu<sup>2</sup>, Dilip Abraham<sup>3</sup>

<sup>1</sup>Imperial College London; <sup>2</sup>Kwame Nkrumah University of Science and Technology; <sup>3</sup>CMC Vellore

GEMS - Genomic Environ...



**Shannon Fitz**

Imperial College London

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.5qpvokmrbl4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvokmrbl4o/v1)

**Protocol Citation:** Shannon Fitz, Alex Shaw, michael Owusu, Dilip Abraham 2024. Zymo OneStep PCR Inhibitor Removal Kit. protocols.io <https://dx.doi.org/10.17504/protocols.io.5qpvokmrbl4o/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 03, 2024

**Last Modified:** October 03, 2024

**Protocol Integer ID:** 108988

**Keywords:** PCR, PCR inhibitor removal, Purification

**Funders Acknowledgement:**

**The Gates Foundation**

**Grant ID:** INV-049092 PA4273



## Abstract

### Features:


- For high quality DNA or RNA that is free of enzymatic inhibitors including polyphenolics, humic/fulvic acids, tannins, melanin, etc.
- Fast, one-step procedure for “cleaning” impure samples prior to PCR, sequencing, RT, etc.

### Description:

The OneStep™ PCR Inhibitor Removal Kit contains all the components needed to efficiently remove contaminants from DNA/RNA preparations that can inhibit downstream enzymatic reactions such as PCR and RT. The column matrix has been specifically designed for efficient removal of polyphenolic compounds, humic/fulvic acids, tannins, melanin, etc. from most impure DNA and RNA preparations. Sample cleanup is as simple as: applying, spinning and recovering a sample from the column.

The steps and abstract for this protocol consist entirely of the manufacturer's instructions with no deviations, as per [https://files.zymoresearch.com/protocols/\\_d6030\\_onestep\\_pcr\\_inhibitor\\_removal\\_kit.pdf](https://files.zymoresearch.com/protocols/_d6030_onestep_pcr_inhibitor_removal_kit.pdf).

## Materials

 Zymo OneStep PCR Inhibitor Removal Kit **Zymo Research Catalog #D6030**

### Total cost per sample:

£0.5

### Approximate time needed per sample:

<10 minutes



## Column Preparation

- 1 Insert column into a Collection Tube.

*\*Please note that matrix in the column may appear dehydrated, or powdery. This is normal.*

- 2 Open the cap, add 600 µl of Prep-Solution and centrifuge at 8,000 x g for 3 minutes.

- 3 Inhibitor Removal:

Transfer the prepared column to a clean 1.5 ml microcentrifuge tube. Add 50 - 200 µl DNA or RNA (in water, TE, or similar) to the Zymo-Spin™ III-HRC Column and centrifuge at 16,000 x g for 3 minutes.

The filtered DNA (or RNA) is suitable for PCR, (RT), and other downstream applications.

## Protocol references

Zymo Research Version 2.0.2-last update, *OneStep™ PCR Inhibitor Removal Kit Cat. No. D6030 (50 spin columns/purifications) Protocol*.

Available: [https://files.zymoresearch.com/protocols/\\_d6030\\_onestep\\_pcr\\_inhibitor\\_removal\\_kit.pdf](https://files.zymoresearch.com/protocols/_d6030_onestep_pcr_inhibitor_removal_kit.pdf) [2024, October 3].