



FEB 01, 2024

## 🌐 Live-cell imaging; cell death assay

Minee-Liane Choi<sup>1</sup>

<sup>1</sup>UCL Institute of Neurology



Minee-Liane Choi

UCL Institute of Neurology, The Francis Crick Institute

### ABSTRACT

This protocol contains instructions for examining cell death in live-cell imaging.

OPEN  ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.rm7vzy1d5lx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzy1d5lx1/v1)

**Protocol Citation:** Minee-Liane Choi 2024. Live-cell imaging; cell death assay. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.rm7vzy1d5lx1/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** May 26, 2022

Last Modified: Feb 01, 2024

PROTOCOL integer ID: 63250

Keywords: ASAPCRN

**1** Cell death is detected using Propidium iodide (PI, Thermo Fisher Scientific) or SYTOX™ Green (SYTOX, Thermo Fisher Scientific) which is excluded from viable cells but exhibits red (PI) and green (SYTOX) fluorescence following a loss of membrane integrity and Hoechst 33342 (Hoechst, Thermo Fisher Scientific).

**2** 20 µM PI or 500 nM SYTOX and 10 uM Hoechst are directly added to the dishes.

**3** Cells are incubated for 15 min.

**4** Live-cell imaging is performed using a confocal microscope.

**Note**

Live-cell imaging is performed using a confocal microscope (Zeiss LSM 710 or 880 with an integrated META detection system). For confocal microscopes, illumination intensity is limited to 0.1-0.2% of laser output to prevent phototoxicity, and the pinhole is set to allow optical slice at approximately 1-2 µm.

**5** Hoechst and PI are excited by 405 nm and 543 nm laser lines with the emission between 405 nm to 470 nm and 570 nm to 640nm respectively.

**6** SYTOX is excited by a 488 nm laser with emissions between 488 nm and 516 nm.

