



Mar 04, 2021

# Real-Time PCR protocol to screen for SARS-COV-2 variants of concern (B.1.1.7, P.1 and B.1.1.35)

Camila M Romano<sup>1</sup>, Jaqueline De Jesus<sup>1</sup>, Alvina Clara Felix<sup>1</sup>, Anderson V de Paula<sup>1</sup>, Pâmela S Andrade<sup>1</sup>, Franciane M de Oliveira<sup>1</sup>, Darlan Cândido<sup>2</sup>, Nuno Faria<sup>2</sup>, William M de Souza<sup>2</sup>, Ester C Sabino<sup>1</sup>

<sup>1</sup>Institute of Tropical Medicine, University of Sao Paulo; <sup>2</sup>Department of Zoology, University of Oxford, Oxford, United Kingdom

*In Development* This protocol is published without a DOI.

Camila Romano  
Institute of Tropical Medicine, University of Sao Paulo, Cli...

SUBMIT TO PLOS ONE

## ABSTRACT

The recently emerged variants of concern (VOC) B.1.1.7 in UK, B.1.1.35 in South Africa and P.1 in Brazil have increased transmissibility and are rapidly spreading worldwide. These variants have several substitutions that allow them to bypass the immune responses and are up to 10 times more contagious than the previous circulating SARS-COV-2 lineages. While in some regions P.1 predominates (as several states in Brazil), in other countries B.1.1.7 or B.1.1.35 are most prevalent. Therefore, real-time surveillance of VOC is the key for rapid public health responses. In this matter, RT-PCR protocols that are capable of detecting VOCs easily and rapidly are needed. We developed and validated a protocol that detect simultaneously the 3 most widespread VOC. Although this protocol does not discriminate among the VOCs, it can reliably detect mixed infection. Importantly, our protocol successfully identified a co-infection by two SARS-COV-2 lineages (P.1 and B.1), confirmed through complete genome sequencing.

## PROTOCOL CITATION

Camila M Romano, Jaqueline De Jesus, Alvina Clara Felix, Anderson V de Paula, Pâmela S Andrade, Franciane M de Oliveira, Darlan Cândido, Nuno Faria, William M de Souza, Ester C Sabino 2021. Real-Time PCR protocol to screen for SARS-COV-2 variants of concern (B.1.1.7, P.1 and B.1.1.35) . **protocols.io** <https://protocols.io/view/real-time-pcr-protocol-to-screen-for-sars-cov-2-va-bszbnf2n>

## KEYWORDS

SARS-COV-2, variants of concern, RT-PCR

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Mar 03, 2021

## LAST MODIFIED

Mar 04, 2021

## PROTOCOL INTEGER ID

47875

## GUIDELINES

This protocol was not set for diagnoses purpose.

## MATERIALS TEXT

Here are the sequences of primers and probes.

A	B	C
Name	Definition (primer/probe)	sequence
NSP6_Forw	primer F	GTGATGCGTATTATGACATGGT
NSP6_Rev	primer R	GATTAGTAACACTACAGCTGATGC
NSP6_Wild	probe for non-VOC	FAM-TGTCTGGTTTAAAGCTAAAAGACTGTG-BHQ
NSP6_Delete	probe that detects VOC	HEX-GGATATGGTTGATACTAGTTTGAAGCT-BHQ

We recommend the [4x TaqMan Fast Virus 1-Step Master Mix \(Applied Biosystems\)](#).

## DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

- 1 Reconstitute the oligos/probes to 100uM with TE or nuclease-free water. Vortex briefly and allow to sit for about 2 minutes.
- 2 Prepare 10uM working solutions (primers and probes). Suggestion: if the original concentration is 100uM, add 10ul of the reagents in 90ul of ultrapure nuclease-free water.
- 3 The reaction was not adjusted for multiplex use, so prepare the reactions separately. We recommend running the probes NSP6 Wild and NSP6 Delete in the same plate/experiment for better results.
- 4 Prepare your Master mix as the following (use primers and probes at 10uM)

A	B	C
Reagents	Volume 1 reaction (ul)	x100
4x TaqMan Fast Virus 1-Step Master Mix	5	500
NSP6_Forw (300nM/reaction)	0.5	50
NSP6_Rev (300nM/reaction)	0.5	50
Probe (300nM/reaction)	0.5	50
RNA	8.5	-
Total	15	-

- 5 Add 6.5ul of the mix to each well and 8.5ul of the RNA template. Don't forget using positive and negative controls.

6 In case of bubbles, spin the plate before cycling.

7 Run the following cycling at the Real time machine:

A	B	C
50°C	2 min	
95°C	10 min	
95°C	15 sec	40x
60°C	1 min	

*Set the thermocycler for the proper FAM and HEX fluorophores.*

8 Interpreting results:

A	B	C
NSP6_Wild	NSP6_Delete	Result
CT ≤ 38	CT > 39/undetected	non- VOC
CT > 39 /undetected	CT ≤ 38	VOC (B.1.17; P.1 or B.1.1.35)
CT ≤ 38	CT ≤ 38	mixed infection
CT > 39/undetected	CT > 39/undetected	invalid/negative