



Mar 01, 2022

DigiGait

Haley Geertsma¹

¹University of Ottawa

1



dx.doi.org/10.17504/protocols.io.b5stq6en

Haley Geertsma
University of Ottawa

This protocol is used to test mice in the DigiGait behavioural test.

DOI

dx.doi.org/10.17504/protocols.io.b5stq6en

Haley Geertsma 2022. DigiGait. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b5stq6en>



protocol ,

Mar 01, 2022

Mar 01, 2022

58931

Mouse Specifics Inc. DigiGait treadmill
DigiGait Imager software
DigiGait Analysis software

- 1

Habituate mice in the testing room for 60 minutes.

1h
- 2

Place mice on the DigiGait treadmill and start the treadmill at 18-22cm/s. Record their movement for 3 seconds then place them back into their home cage.

1m
- 3

Analyze 3-second clips with DigiGait Analysis Software.

1m

