

Aug 30, 2024

## JAX-Sen: Mouse heart dissociation for single-cell RNA sequencing

DOI

[dx.doi.org/10.17504/protocols.io.rm7vzj55xlx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzj55xlx1/v1)

Ramalakshmi Ramasamy<sup>1</sup>, Juliana Alcoforado Diniz<sup>1</sup>, Patrick Fleming<sup>1</sup>, Jessica Garofalo<sup>1</sup>, Paul Robson<sup>1,2</sup>

<sup>1</sup>The Jackson Laboratory for Genomic Medicine, Farmington, CT, USA;

<sup>2</sup>Department of Genetics and Genome Sciences, University of Connecticut School of Medicine, Farmington, CT, USA

Cellular Senescence Net...



Harshpreet Chandok

The Jackson Laboratory

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.rm7vzj55xlx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzj55xlx1/v1)

**Protocol Citation:** Ramalakshmi Ramasamy, Juliana Alcoforado Diniz, Patrick Fleming, Jessica Garofalo, Paul Robson 2024. JAX-Sen: Mouse heart dissociation for single-cell RNA sequencing. [protocols.io](https://dx.doi.org/10.17504/protocols.io.rm7vzj55xlx1/v1) <https://dx.doi.org/10.17504/protocols.io.rm7vzj55xlx1/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 02, 2024

**Last Modified:** August 30, 2024

**Protocol Integer ID:** 102774

**Keywords:** JAX-Sen



**Funders Acknowledgement:**

**National Institute on Aging  
(NIA) JAX-Sen Senescence  
Tissue Mapping Center Grant  
Grant ID: ID: U54 AG079753**

## Abstract

These samples are part of the JAX-Sen project in the SenNet Consortium. We aim to study and characterize senescence in the C57Bl/6 mouse heart. We compare young (6 months old) and old (24 months old) mouse hearts using scRNAsequencing. This protocol describes the single-cell dissociation of the hearts before library preparation and sequencing.

## Reagents and Materials

- 1
  - 5 ml Protein lobind tubes
  - 1.5 or 2 ml Protein lobind tubes
  - Micro scissors/ scalpel
  - Ice-cold PBS
  - Wide-bore pipette tips
  - 0.25% Trypsin/EDTA (Gibco, 25200056). - 2 ml/sample
  - 20 mg/mL collagenase A and B (Sigma, 10103578001, 11088807001) – 2 ml/sample
  - DMEM + 10% FBS
  - 100 µm cell strainer (Corning, 431752)
  - HBSS (Ca<sup>2+</sup>/Mg<sup>2+</sup> free) (Gibco, 14170120) – 2 ml/sample
  - ACK lysis buffer – 2 ml/sample
  - SS\_buffer – 5 ml/sample
  - 40 um cell strainer
  - 70 uM strainer

### SS\_buffer

A	B	C	D	E	F	G	H	I
	Stock co nc.	Working conc.	50 ml	100 ml	250 ml	500 ml	20 ml	
EDTA	0.5 M	2 mM	0.2	0.4	1	2	0.08	ml
BSA	10%	2%	10	20	50	100	4	ml
FBS	100%	15%	7.5	15	37.5	75	3	ml
DPBS	1X	make up	32.2	64.4	161	322	12.88	ml

## Procedure

- 2 Collect shipment from the warehouse.
- 2.1 Fill up ice buckets and prep them with reagents:  
Bucket 1 (small): PBS, HBSS, DMEM+FBS, SS\_buffer  
Bucket 2 (large): For mincing tissue
- 2.2 Confirm if the sample is cold. Take pictures.



- 2.3 Place tubes on ice and mince the tissue (<1mm) with sterile micro scissors.
- 2.4 After mincing, transfer the sample into labeled 5 ml tubes.
- 2.5 Centrifuge at **500 x g, 4°C, 5 min** to get rid of the supernatant.
- 2.6 Resuspend the pellet in **ice-cold PBS (Ca<sup>+</sup> and Mg<sup>+</sup> free)** and pipette up and down gently with wide-bore pipette tip.
- 2.7 Centrifuge at **500 x g, 4°C, 5 min** to get rid of the supernatant.
- 2.8 Resuspend the pellet in **ice-cold PBS (Ca<sup>+</sup> and Mg<sup>+</sup> free)** and pipette up and down gently with wide-bore pipette tip.
- 2.9 Centrifuge at **500 x g, 4°C, 5 min** to get rid of the supernatant.
- 2.10 Add 2 ml of **0.25% Trypsin/EDTA** to the pellet and enzymatically digest at **37°C for 10 min** with slow rocking. ***Check at 5 mins for dissociation, under the microscope.***
- 2.11 Stop the digestion by adding an equal volume (1ml) of **DMEM/10% FBS**.
- 2.12 Centrifuge at **500 x g, 4°C, 5 min** to get rid of the supernatant.
- 2.13 Resuspend the pellet 2 ml of **20 mg/mL collagenase A and B** mixture and incubate the samples at 37 °C until the tissues were into single cells.
- 2.14 Vigorously pipette samples with 1 ml wide-bore pipette tips or if needed, regular tips.
- 2.15 Filter cells through **MACS SmartStrainers (70um)** and wash 1 ml **HBSS (Ca<sup>2+</sup>/Mg<sup>2+</sup> free)** through the strainer.
- 2.16 Centrifuge at **500 x g, 4°C, 5 min** to get rid of the supernatant.



- 2.17 Resuspend the cell pellet in 2 mL **ACK Lysing buffer**, and incubate on **ice for 3 min** to exclude RBCs.
- 2.18 Add 5 mL **SS\_buffer** to stop ACK Lysing buffer reaction, then collect cells by centrifugation **500 x g, 4°C, 5 min**.
- 2.19 Resuspend the pellet in **1ml** of **HBSS (Ca<sup>2+</sup>/Mg<sup>2+</sup> free)** (Gibco, 14170120).
- 2.20 Filter the cell suspension through **Flowmi 40um cell strainers** tube to remove debris and non-dissociated tissue fragments.
- 2.21 Collect cells by centrifugation **500 x g, 4°C, 5 min** and proceed to Flex protocol.